

Power NON-Articulating Legrest Kits POWELR18/POWELR20 (9/15/01 thru 4/30/02)

NOTE: Power NON-articulating legrest use sector style front rigging hangers and do not come included with kits. See page titled A.S.B.A. Seat Frame and Hardware for hanger assemblies.

NOTE: Torque SP with RII electronics must be upgraded to MKIV A electronics.

NOTE: For 3G bases, a new controller shroud cover may be required. See page titled Controller Shroud Covers.

POWELR18/ POWELR20 with 4-way Switch and Non Quad Link Mounting Plate Assembly.

Part Number	Description	Quantity
1109317-NLA See Note(s) A, B	Power Non-Articulating Legrest Kit	1
1109313-NLA See Note(s) A	Legrest, Right (13"-20" Adj. Length)	1
1109314-NLA See Note(s) A	Legrest, Left (13"-20" Adj. Length)	1
1101520-NLA	Power Leg Harnesses and Mounting Hardware	1
1104248-NLA	T.A.C. Controller with Hardware	1
1080116-NLA	4-Way Switch	1
1077519	4-Way Switch Mounting Plate with Hardware	1
2001682	Cable Tie (8" Long)	10
1092659-NLA	Controller Bracket with Fastener Strips	1
1104911	Instruction Sheet, Power Articulating Legrest Kit	1
NOTE: A - Power NON-articulating legrest and KITS before 9/15/01 are NLA. See power articulating legrest or KITS dated after 7/26/04. Legrest must be ordered in pairs to match. B - Kit includes the 9 items listed below the kit part number		

Power NON-Articulating Legrest Kits POWELR18/POWELR20 (9/15/01 thru 4/30/02)

NOTE: Power NON-articulating legrest use sector style front rigging hangers and do not come included with kits. See page titled A.S.B.A. Seat Frame and Hardware for hanger assemblies.

NOTE: Torque SP with RII electronics must be upgraded to MKIV A electronics.

NOTE: For 3G bases, a new controller shroud cover may be required. See page titled Controller Shroud Covers.

POWELR18/ POWELR20 with 4-way Switch and Quad Link Mounting Plate Assembly

Part Number	Description	Quantity
1109318-NLA See Note(s) A, B	Power Non-Articulating Legrest Kit	1
1109313-NLA See Note(s) A	Legrest, Right (13"-20" Adj. Length)	1
1109314-NLA See Note(s) A	Legrest, Left (13"-20" Adj. Length)	1
1101520-NLA	Power Leg Harnesses and Mounting Hardware	1
1104248-NLA	T.A.C. Controller with Hardware	1
1080116-NLA	4-Way Switch	1
1077520	4-Way Switch Mounting Plate with Hardware	1
2001682	Cable Tie (8" Long)	10
1092659-NLA	Controller Bracket with Fastener Strips	1
1104911	Instruction Sheet, Power Articulating Legrest Kit	1

NOTE: A - Power NON-articulating legrest and KITS before 9/15/01 are NLA. See power articulating legrest or KITS dated after 7/26/04. Legrest must be ordered in pairs to match.

B - Kit includes the 9 items listed below the kit part number