

## Power Articulating Legrest Kits ELRPW (5/1/02-7/26/04)

**NOTE: Power articulating legrest use sector style front rigging hangers and do not come included with kits. See page titled A.S.B.A. Seat Frame and Hardware for hanger assemblies.**

**NOTE: Torque SP with RII electronics must be upgraded to MKIV A electronics.**

**NOTE: For 3G bases, a new controller shroud cover may be required. See page titled Controller Shroud Covers.**

### ELRPW with 4-Way Switch and Non-Quad Link Mounting Plate Assembly

Part Number	Description	Quantity
1111906-NLA  See Note(s) A, B1, B	Kit, ELRPW Power Articulating Legrest Kit	1
1110743-NLA See Note(s) B1, B	Legrest, Right (13"-20" Adj. Length)	1
1110744-NLA See Note(s) B1, B	Legrest, Left (13"-20" Adj. Length)	1
1101520-NLA	Power Leg Harnesses and Mounting Hardware	1
1104248-NLA	T.A.C. Controller with Hardware	1
1080116-NLA	4-Way Switch	1
1077519	4-Way Switch Mounting Plate with Hardware	1
2001682	Cable Tie (8" Long)	10
1092659-NLA	Controller Bracket with Fastener Strips	1
1104911	Instruction Sheet, Power Articulating Legrest Kit	1
<p><b>NOTE:</b> A - Kit includes the 9 items listed below the kit part number            B - Due to a product enhancement that provides a higher pivot point at the knee, a customer must order a complete right and left side legrest assembly. A kit is required when ordering Power Articulating Legrest for the first time.            B1 - See page titled "Power Articulating Legrest Kits After 7/26/04".</p>		

## Power Articulating Legrest Kits ELRPW (5/1/02-7/26/04)

**NOTE: Power articulating legrest use sector style front rigging hangers and do not come included with kits. See page titled A.S.B.A. Seat Frame and Hardware for hanger assemblies.**

**NOTE: Torque SP with RII electronics must be upgraded to MKIV A electronics.**

**NOTE: For 3G bases, a new controller shroud cover may be required. See page titled Controller Shroud Covers.**

### ELRPW with 4-Switch and Quad Link Mounting Plate Assembly

Part Number	Description	Quantity
1111724-NLA See Note(s) B1, B	Kit, ELRPW Power Articulating Legrest Kit	1
1110743-NLA See Note(s) B1, B	Legrest, Right (13"-20" Adj. Length)	1
1110744-NLA See Note(s) B, B1	Legrest, Left (13"-20" Adj. Length)	1
1101520-NLA	Power Leg Harnesses and Mounting Hardware	1
1104248-NLA	T.A.C. Controller with Hardware	1
1080116-NLA	4-Way Switch	1
1077520	4-Way Switch Mounting Plate with Hardware	1
2001682	Cable Tie (8" Long)	10
1092659-NLA	Controller Bracket with Fastener Strips	1
1104911	Instruction Sheet, Power Articulating Legrest Kit	1
<p><b>NOTE:</b> B - Due to a product enhancement that provides a higher pivot point at the knee, a customer must order a complete right and left side legrest assembly. A kit is required when ordering Power Articulating Legrest for the first time. B1 - See page titled "Power Articulating Legrest Kits After 7/26/04".</p>		