


MK5 SAC-E and (MK6i) Elevate Integrated Actuator Module (EIAM)



MK5 SAC-E and (MK6i) Elevate Integrated Actuator Module (EIAM)

Item Number	Note	Part Number	Description	Quantity Per	
				Package	Assembly
1		1126033-NLA	NLA - Control, MK5 SAC-E	1	1
1		1140036	Module, Elevate Integrated Actuator (EIAM) - MK6i (54" L)	1	1
2	A	1125685-NLA	NLA - Kit, MK5 SAC-E/ MK6i EIAM Mounting Hardware	1	1
3			Screw, Phillips Pan Head Tap (#10-32 x 1/2")	1	2
4			NLA - Bracket, Mounting, Black	1	1
5			NLA - Stud Plate, Quad Link	1	2
6	B	1119418	Cable, Auxiliary Power	1	1
NOTE: A - Includes items 3-5					
B - Auxiliary Power cable not used with MK6i electronics					