


MK5 TRECM, Mercury Free Switches and Mounting Hardware

Conventional Tilt and Tilt w/ BPO Option (Before 1/17/05)

As of 7/1/2009, switches are mercury free.



MK5 TRECM, Mercury Free Switches and Mounting Hardware

Conventional Tilt and Tilt w/ BPO Option (Before 1/17/05)

As of 7/1/2009, switches are mercury free.

Item Number	Note	Part Number	Description	Quantity Per	
				Package	Assembly
1		1123263-NLA	NLA - Controller, MK5 TRECM	1	1
2		1080634-NLA	NLA - Harness, Drive Lockout Mercury Free Ball Switch	1	1
3		1123268-NLA	NLA - Harness, Tilt Mercury Free Ball Switch	1	1
4		40157X019-NLA	NLA - Spacer, Coved (1/4 x 7/8 x 7/16")	1	2
5		1074642-NLA	NLA - Screw, Hex Head (1/4-20 x 1-3/8")	1	2
6		1125792-NLA	NLA - Bracket, Mercury Free Switch Mounting, Black	1	1
7	A,C	1125730-NLA	NLA - Kit, Drive Lockout/Recline Switch Mounting Hardware	2	1
8			Screw, Hex Head w/ Patch (1/4-20 x 1/2")	1	4
9			NLA - Bracket, Switch Mounting, Black	1	2
10			Washer (1/4 x 1/2 x 1/16")	1	4
11			Locknut (1/4-20)	1	4
12	B	1128348-NLA	NLA - Kit, MK5 TRECM Mounting Hardware For Conventional Tilt Only	1	1
13			NLA - Bracket, Seat Pan Wiring, Black	1	2
14			Screw, Phillips Pan Head Tap (#10-32 x 1/2")	1	4
15			NLA - Bracket, Controller Mounting, Black	1	1
16			NLA - Screw, Phillips Pan Head w/ Patch (#10-32 x 3/8")	1	2
17			NLA - Screw, Hex Head (1/4-20 x 1-1/8")	1	2
18			Half Clamp, Non-Threaded Holes	1	1
19			NLA - Half Clamp, Threaded Holes	1	1
20		1119418	Cable, Auxiliary Power	1	1

NOTE: A - Includes items 8-11
 B - Includes items 13-19
 C - Kit can be used with the tilt switch as well