

Battery Boxes & Trays R2 jr. (14" Wide)

Item Number	Note	Part Number	Description	Quantity Per	
				Package	Assembly
1		1079203	Decal, Warning Pinch Point	1	1
2		1026032	Washer, Nylon (11/32 x 11/16 x 1/16")	1	2
3		1086166-NLA	NLA - Screw, Hex Head (5/16-18 x 3/4")	1	2
4		1074133	Washer, Steel (5/16 x 3/4 x 1/16")	1	2
5		1028769	Locknut (5/16-18)	1	2
6		1090275-NLA	NLA - Tray, Battery with Buckle and Warning Decal	1	1
7		1086160-NLA	NLA - Bracket, Pivot, Black - Right	1	1
8	A	1086167-NLA	NLA - Battery Tray Assembly	1	1
9	B,C	1057941-NLA	NLA - Battery Box Top Assembly, Dual w/ Single Extended Connector, UI	1	2
10		1009268-NLA	NLA - Battery Box Bottom Assembly, U1	1	2
11	D	1057942-NLA	NLA - Battery Box Assembly, Dual w/ Single Extended Connector, UI	1	2
12		1058520-NLA	NLA - Fastener, Panel Snap-In	1	5
13		1079247	Washer, Flat (21/64 x 53/64 x 7/64")	1	1
14		1090276-NLA	Bracket, Pivot with Warning Decal, Black - Left	1	1
15		1050161-NLA	NLA - Screw, Hex Head w/ Patch (5/16-18 x 1")	1	3
16		1026158	Locknut (5/16-18)	1	2
	C,F	1089385	Fuse, MIDI (Non-G/B & G/B) (80 AMP)	1	2
	E,F	1060318-NLA	NLA - Hardware, Battery Terminal, UI Batteries (Pair)	2	1

NOTE: A - Includes Items 2-7.
 B - It is recommended that new terminal hardware be ordered when a new battery box top is ordered.
 C - As of 11/05/99 the battery harness fuse changed from a Bussman 60AMP to a MIDI 80AMP fuse. This battery box top is compatible for all chairs manufactured before 11/05/99 and must be ordered whenever a replacement Bussman 60AMP fuse is requested.
 D - Includes Items 9 and 10.
 E - Hardware includes small battery terminal caps (gray and orange), 1/4-20 x 7/8" hex flange screws, 1/4-20 hex flange nuts, cable ties, locking tab washers and instruction sheet for two batteries.
 F - Not Shown.