



# Front Rigging Hanger Assembly


Before 5/22/06



Plug Button used when  
Front Rigging Hanger  
is not present



**Hanger Brackets for  
all Swingaway Front Riggings**



**Hanger Brackets for  
all Lift Off Front Riggings**

## Front Rigging Hanger Assembly Before 5/22/06

Item Number	Note	Part Number	Description	Quantity Per	
				Package	Assembly
1		1074640	Screw, Allen Head (1/4-28 x 1-1/4")	1	2
2		1085111	Nut, Hex (1/4-28)	1	2
3	A	1054866	Kit, Sector Block Hardware (Pair)	1	1
4			Screw, Hex Head w/ Patch (1/4-20 x 2")	1	2
5			Washer (1/4 x 1/2 x 1/16")	1	2
6			NLA - Block, Sector - Right	1	1
6			NLA - Block, Sector - Left	1	1
7		1086389	Tube, Support, Front Rigging Hanger, Black - Right	1	1
7		1086390	Tube, Support, Front Rigging Hanger, Black - Left	1	1
8		1127676	Package, Plug Button, Black (1")	10	1
9	B	1142975	Kit, Sector Block Hardware, Lift Off (Pair)	1	1
10			Screw, Hex Head (1/4-20 x 1-5/8")	1	2
11			Washer (1/4 x 1/2 x 1/16")	1	2
12			Block, Sector - Right	1	1
12			Block, Sector - Left	1	1

NOTE: A - Used for all Swingaway Front Riggings. Includes items 4-6.  
B - Used for all Lift Off Front Riggings. Includes items 10-12