


Base Assembly

RPS440 Power Stand-Up Lift



Base Assembly

RPS440 Power Stand-Up Lift

Item Number	Note	Part Number	Description	Quantity Per	
				Package	Assembly
1		1143814-NLA	NLA - Footplate, Removable	1	1
2	A	1143812-NLA	NLA - Kit, Front Caster with Hardware	1	2
3	B	1143813-NLA	NLA - Kit, Rear Caster with Hardware	1	2
NOTE: A - Includes (1) Caster, (1) Screw, (1) Washer, (1) Spacer and (1) Nut. B - Includes (1) Caster with Brake, (1) Bolt, (1) Bushing, (2) Washers and (1) Nut.					