


# Rollator (Four Wheeled) - Model 65450



Wheel assembly shown on the rollator is no longer available as a service part. See wheel assembly item #12 shown below as the replacement.

This model wheel assembly is No Longer Available.

**Wheel Assembly Replacement for Model 65450**

## Rollator (Four Wheeled) - Model 65450

Item Number	Note	Part Number	Description	Quantity Per	
				Package	Assembly
1		1142615	Hand Grip (Pair)	1	1
2	A	1142616-NLA	NLA - Brake Assembly	1	1
3	B	1142618	Lever, Hand (Pair)	2	1
4		1142622-NLA	NLA - Seat w/ Hardware	1	1
5		1142626-NLA	NLA - Connection Linkage Assembly	1	2
6		1142625-NLA	NLA - Fork, Front	1	2
7		1142627-NLA	NLA - Basket	1	1
8		1142619-NLA	NLA - Folding Lock with Bracket	1	1
9		1142620-NLA	NLA - Bracket, Black	1	1
10	D	1142624-NLA	NLA - Caps, End, Black	1	1
11	E	1142617-NLA	NLA - Pad, Brake	1	2
12	C	1142621	Wheel Assembly w/ Hardware (7-1/2")	1	4
	F	1142623-NLA	NLA - Tray	1	1

NOTE: A - Includes Brake Handle with Cable, Brake Pad is not included.  
 B - Includes (2) Levers and Hardware.  
 C - Wheel Assembly that originally came with Rollator model 66450 is no longer available as a service part. Must verify which wheel assembly type the customer has. A quantity of (4) ea. of part #1142621 may be required to match.  
 D - Includes (2) Fork Caps and (4) Wheel Caps.  
 E - Includes (1) Brake Pad and (1) Screw.  
 F - Not Shown.