

Rollator (Four Wheeled) - Model 65850



Rollator (Four Wheeled) - Model 65850

Item Number	Note	Part Number	Description	Quantity Per	
				Package	Assembly
1	A	1142628	Grip, Hand (Pair)	1	1
2		1142640	Brake Handle with Cable	1	2
3		1142643-NLA	NLA - Back Support Assembly	1	1
4		1142644-NLA	NLA - Package, Spring Button	10	2
5		1142633-NLA	NLA - Seat w/ Hardware	1	1
6		1142638-NLA	NLA - Basket	1	1
7		1142639-NLA	NLA - Connection Linkage Assembly	1	2
8		1142641-NLA	NLA - Brake Pad Assembly - Right	1	1
9		1142642-NLA	NLA - Brake Pad Assembly - Left	1	1
10		1142646-NLA	NLA - Kit, Wheel Assembly - Front	1	2
11	D	1142634-NLA	NLA - Fork Assembly, Front	1	2
12		1142637-NLA	NLA - Package, Strap	10	1
13	C	1142631-NLA	NLA - Package, Knob	1	1
14		1141759-NLA	NLA - Kit, Wheel w/Hardware - Rear	1	2

NOTE: A - Includes (2) Hand Grips.
C - Includes (2) Adjustable Knobs with Hardware, (2) Anti-Rattles and (2) Basket Guides.
D - Replacement wheel has different number of spokes that does not impact performance. It is recommended to replace Fork Assembly in pairs to maintain appearance.