


Frame Hardware

Lynx & Zoom-3


Note: Frame is not sold as a service part.

Note: All Front and Rear Frames shown for reference only.

Lynx SX3P/SX3
Lynx LX3/LX3^{PLUS}



Zoom-3



Frame Hardware

Lynx & Zoom-3

Note: Frame is not sold as a service part.

Item Number	Note	Part Number	Description	Quantity Per	
				Package	Assembly
1		1088660-NLA	NLA - Frame Hardware, Frame, Front - Includes: 1 each Safety Pin Assembly 2 each Cable Tie 1 each Chain 1 each Pin, Safety 1 each Bolt, Hex Head (M6 x 1.0 x40) 1 each Nut, Self-Lock (6 x 10)	1	1
2		1088792-NLA	NLA - Hardware, Frame - Rear - Lynx SX3P/SX3/LX3/LX3Plus and Panther LX4 Includes: 1 each Rubber	1	1
3		1088767-NLA	NLA - Knob, Seat Post - Lynx SX3P/SX3/LX3/LX3Plus and Panther LX4	1	1
4		1105914-NLA	NLA - Hardware, Frame - Zoom-3 Includes: 1 ea. Knob, Seat Post 2 ea. Rubber Cover 2 ea. Wheel, Anti-Tip 2 ea. Screw, Hex Cap (M6 x 1.0 x 40) 2 ea. Nut, Self-Lock (6 x 10) 1 ea. Pin, Connecting	1	1