


Back Upholstery, Nylon Sling (CGE) (16"-22" Wide)



Back Upholstery Hardware

Item Number	Note	Part Number	Description	Quantity Per	
				Package	Assembly
	A	See Chart	Single Liners	1	1
	A	See Chart	Dual Liners - Nylon	1	1
	A	1026298	Package, Screw with Washer (#10 x 7/8") (Pair)	2	1

NOTE: A - Not Shown

Single Liners

Width	16"-17" High	18"-19" High	20" High	21"-22" High
16" W	1045175	1045207	1068947	1045159
17" W	1045176	1045208	1068948	1045160
18" W	1045177	1045209	1068949	1045161

**Back Upholstery, Nylon Sling
(CGE) (16"-22" Wide)
Dual Liners**

Width	16"-17" High	18" - 19" High	20" High	21" - 22" High
16" W	1045180	1045212	1068952	1045164
17" W	1045181	1045213	1068953	1045165
18" W	1045182	1045214	1068954	1045166
19" W	1045183	1045215	1068955	1045167
20" W	1045184	1045216	1068956	1045168
21" W	1045185	1045217	1068957	1045169
22" W	1045186	1045218	1068958	1045170

NOTE: Dual liner back upholstery is used with chairs that have a seat deeper than 17" or back wider than 18".