


Spreader Bar Assembly

Contoura (CGR) and (CGTR/CGTRE) (16"-24" Wide)



Spreader Bar Assembly

Contoura (CGR) and (CGTR/CGTRE) (16"-24" Wide)

Item Number	Note	Part Number	Description	Quantity Per	
				Package	Assembly
1	A	1142957	Spreader Bar Assembly, Contoura Recline and Tilt/ Recline	1	1
2		1125765	Kit, Spreader Bar Mounting Hardware For Extruded Back Canes	1	1
3			Screw, Hex Head w/ Patch (1/4-20 x 1/2")	1	8
4			Washer (1/4 x 1/2 x 1/16")	1	8
5			T-Nut (1/4-20)	1	8
NOTE: A - Includes items 3-5					