


Front Shrouds


Lynx & Zoom-3

Note: Shroud Assemblies include Carpet with Hook and Loop and Hardware.


Lynx SX3P/SX3



Lynx LX3/LX3^{Plus}



Zoom-3



Front Shrouds Lynx & Zoom-3

Item Number	Note	Part Number	Description	Quantity Per	
				Package	Assembly
1		1088687-NLA	NLA - Carpet with Hook and Loop	1	1
1		1088688-NLA	NLA - Carpet with Hook and Loop	1	1
1		1105916-NLA	NLA - Carpet with Hook and Loop, Zoom-3	1	1
2		1088710-NLA	NLA - Boot	1	1
3		1088711-NLA	NLA - Hardware, Shroud - Includes: 2 each Washer (6 x 18 x 1.6) 2 each Washer, Spring (6 x 10) 2 each Screw (M6 x 1.0 x 20)	1	1
4		1088678-NLA	NLA - Shroud with Carpet and Hardware, Candy Red	1	1
4		1088679-NLA	NLA - Shroud with Carpet and Hardware - Metallic Green	1	1
4		1088680-NLA	NLA - Shroud with Carpet and Hardware, Deep Blue	1	1
4		1104157-NLA	NLA - Shroud with Carpet and Hardware, Magic Black - Lynx LX3Plus	1	1
5		1088793-NLA	NLA - Hardware, Shroud - Lynx SX3P/SX3 Includes: 2 each Washer (6 x 18 x 1.6) 2 each Washer, Spring (6 x 10) 2 each Screw (M6 x 1.0 x 16)	1	1
6		1088674-NLA	NLA - Shroud with Carpet and Hardware, Gray (Before S/N 00F623035)	1	1
6		1102246-NLA	NLA - Shroud with Carpet and Hardware - Lynx SX3P (After S/N 00F623034) - Ivory	1	1
6		1088675-NLA	NLA - Shroud with Carpet and Hardware, Candy Red	1	1
6		1088677-NLA	NLA - Shroud with Carpet and Hardware, Deep Blue	1	1
7		1105917-NLA	NLA - Hardware, Shroud - Zoom-3 4 ea. Washer (6 x 18 x 1.6) 4 ea. Plate Nut, U-Type 4 ea. Screw (M6 x 1.0 x 16)	1	1
8		1105915-NLA	NLA - Shroud with Carpet and Hardware - Zoom-3 - Candy Red	1	1
		1076178-NLA	NLA - Seat Hardware	--	--