

Rollator Bariatric (Four Wheeled) - Model 66500



Rollator Bariatric (Four Wheeled) - Model 66500

| Item Number | Note | Part Number | Description | Quantity Per | |
|-------------|------|-------------|----------------------------------|--------------|----------|
| | | | | Package | Assembly |
| 1 | A | 1141748-NLA | NLA - Kit, Hand Grip | 1 | 1 |
| 2 | B | 1141749-NLA | NLA - Brake Assembly (Pair) | 2 | 1 |
| 3 | | 1141750-NLA | NLA - Wheel with Hardware - Rear | 1 | 1 |
| 4 | C | 1141751-NLA | NLA - Front Caster Assembly | 1 | 1 |
| 5 | | 1141752-NLA | NLA - Back Rest Assembly | 1 | 1 |
| 6 | | 1141753-NLA | NLA - Basket | 1 | 1 |
| 7 | | 1141754-NLA | NLA - Seat w/ Hardware | 1 | 1 |
| | D,E | 1141755-NLA | NLA - Kit, Frame Hardware | 1 | 1 |

NOTE: A - Includes (2) Hand Grips.
 B - Includes (2) Brake Handle Assemblies, (2) Cables, (2) Break Pads and Bottom Caps.
 C - Includes (1) Fork and (1) Wheel with Hardware.
 D - Includes (2) Anti-Rattles, (2) Push Buttons, (2) Adjustable Knobs and (2) Brake Cord Holders.
 E - Not Shown.