


Shrouds and Attaching Hardware



Shrouds and Attaching Hardware

Item Number	Note	Part Number	Description	Quantity Per	
				Package	Assembly
1		1026114	Screw, Phillips Pan Head (#10-32 x 1/2")	1	2
2		1145420	Shroud, Controller, Black	1	1
3		1143701-NLA	NLA - Decal, Model - M41	1	1
4		1028484	Screw, Phillips Pan Head Tap (#10-32 x 1/2")	1	4
5		1035900	Decal, Transportation	1	1
6		1119281-NLA	NLA - Plate, Name, Patent Numbers	1	1
7		1114276	Screw, Phillips Pan Head (#8-16 x 5/8")	1	6
8		1028516	Washer, Steel (3/16 x 7/16 x 1/16")	1	6
9	A	1145423	Top Shroud Assembly, Red	1	1
9	A	1145424	Top Shroud Assembly, Blue	1	1
10		1145419	Shroud, Black - Front	1	1
11		1025735	Locknut (#10-32)	1	2
12		1095489	Decal, Invacare 3-D Logo (2-1/2" L)	1	1

NOTE: A - Shroud includes decals, items 3, 5-6 & 12.