


Recline Interface Support Frame Hardware

For 3G Storm & TDX SP/SP-GT



Recline Interface Support Frame Hardware

For 3G Storm & TDX SP/SP-GT

Item Number	Note	Part Number	Description	Quantity Per	
				Package	Assembly
1	A	1144589-NLA	NLA - Kit, Recline Interface Support Frame Hardware	1	1
2			Washer, Nylon (3/8 x 7/8 x 1/16")	1	4
3			NLA - Washer, Nylon (13/32 x 1 x 1/4")	1	2
4			Screw, Socket Head (3/8-16 x 2-1/4")	1	2
5			Washer, Nylon (3/8 x 13/16 x 1/16")	1	4
6			NLA - Locknut, High Profile (3/8-16)	1	4
7			Screw, Socket Head w/ Patch (3/8-16 x 1-1/2")	1	2
NOTE: A - Includes items 2-7					