


## Seat Frame Mounting Hardware Elevate (16"-22" Wide) (CGE) For 3G Storm



Item Number	Note	Part Number	Description	Quantity Per	
				Package	Assembly
1	A	1144590-NLA	NLA - Kit, Seat Frame Mtg. Hdwr. 16"-22" Wide (3G Storm Elevate Only)	1	1
2			NLA - Screw, Socket Flat Head (5/16-18 x 1-1/2")	1	4
3			Washer, Steel (11/32 x 11/16 x 1/16")	1	4
4			Locknut (5/16-18)	1	4

NOTE: A - Includes items 2-4