

**AVIVA® FX
Models:
IFX-20R
IFX-20SP
IFX-20MP
Power Wheelchair**

**For Motion Concepts Ultra Low Maxx Rehab
and Power Positioning System**



Yes, you can.

Parts Catalog Usage Guide

The information contained in this document is subject to change without notice. This catalog contains service parts for Form No. 19_118. If you have any questions regarding this catalog, please contact one of our Customer Service Representatives.

For Customer Service Support Call:

1-800-333-6900

For Faxed In Orders Call:

1-800-678-4682 or

1-440-366-0709

For Technical Service Call:

1-800-832-4707

Literature Can Be Found at:

www.invacare.com or

1-800-828-6272

SYMBOL DEFINITION

Throughout this catalog you will notice letters to draw your attention to specific footnotes which are located at the bottom of the page or below specified charts.

Examples:

A,B,C etc. = Special Comments

NOTE: = General information pertaining to that specific page

NLA = No Longer Available

HARDWARE MEASUREMENTS

The measurements for the hardware will be listed as follows:

Screws, Bolts Diameter X Thread Count/Pitch X Length

Washers Inside Diameter X Outside Diameter X Thickness

Spacers Inside Diameter X Outside Diameter X Length

For your convenience, we have supplied an area for you to list your Invacare Representatives.

Name	Phone Numbers

© 2022 Invacare Corporation. Elyria, Ohio. Invacare, Yes, you can., the Blue Medallion Design, SureStep, Dynamic is a registered trademark of Dynamic Controls Limited. Super Lube is a registered trademark of Synco Chemical Corporation.

TABLE of CONTENTS

Description	Page Number
Accessories - Positioning Straps	211
Accessories - TRRO Pelvic Belt and Hardware	213
Accessories - Vent Tray Assembly w/ Shroud	215
Anti-Tip Assembly	23
Battery Box Assembly, Front and Rear Plates	15
Battery Charger & Battery Harnesses GP24	13
Battery Harness, Base, AVIVA FX	17
Caster & Forks - 8" Rear - Single Sided Fork, Nylon	43
Caster & Forks - 8" Rear - Single Sided Fork, Polycarbonate	45
Circuit Breaker 63A	19
Electronics - 4-Way Toggle Switch (FWTM) and Mounting Hardware	201
Electronics - Accessories - Remote (Joystick) Tops	209
Electronics - Accessories - USB Charger and LiNX Access Key (LAK)	207
Electronics - ASL 611 Mount For The Egg Switch	195
Electronics - Backrest Bracket Mounting Hardware, Recline for (DLX-ACT400, DLX-IN500, DLX-OUT500, GLM-CONX4)	176
Electronics - Backrest Bracket Mounting Hardware, Recline for (DLX-ACT400, GLM-CONX4)	178
Electronics - Backrest Bracket Mounting Hardware, Sling Back Only for (DLX-IN200, DLX-IN500, DLX-OUT500, GLM-CONX4)	182
Electronics - Backrest Bracket Mounting Hardware, Tilt or ASBA-2 for DLX-IN200, DLX-ACT400, DLX-IN500, DLX-OUT500, GLM-CONX4	180
Electronics - C629 Junction Box, Shroud Cover and Hardware	186
Electronics - Cables - Bus (PWH), Adapter, Sensor, Actuator	188
Electronics - Dual (2PB)/ Quad (QPB) Push Button Switch w/ Multi-Mount Hardware	199

TABLE of CONTENTS

Description	Page Number
Electronics - Dual Toggle Switch (2SW), Monoport "Y" Phono Cable (MPY)	197
Electronics - G-Trac-LX - LiNX Gyro Module with Mounting Hardware	174
Electronics - Joystick Tubes, Offset and Short	162
Electronics - Kit, Tilt Angle Sensor Assembly w/ Hardware	193
Electronics - LiNX CONX4 4-way Expansion Block, Shroud Cover and Hardware	184
Electronics - LiNX DLX-ACU200 Compact Attendant Remote	203
Electronics - LiNX DLX-CR400 Compact Remote (CREM) and Mounting Hardware	166
Electronics - LiNX DLX-CR400LF Compact Remote Low Force (CREM-LF) and Mounting Hardware	168
Electronics - LiNX DLX-SW100 10-Way Switch and Mounting Hardware	205
Electronics - LiNX Power Module 120AL and Mounting Hardware	172
Electronics - Quad Link Assembly, Remote	164
Electronics - Remotes and Displays	158
Electronics - Swingaway Mounting Hardware, DLX-REM 500	170
Electronics - Toggle Switch Kit (SW400TGL) For LiNX DLX-REM400 Remote	160
Electronics - Wire Guide Bracket and Mounting Hardware, ULMT/ULMTR with AELM	191
Fender w/ Hardware, Drive Wheel, Nylon	33
Fender w/ Hardware, Drive Wheel, Polycarbonate	35
Front Riggings - 70 Degree HD Swing Away Front Rigging (M70HD)	130
Front Riggings - Accessories, Air Mesh Calf Panel CPMD/ CPLG for M70HD	156
Front Riggings - Adjustable Angle Footplate, Pivot and Multi-Axis, Independent Front Riggings	154
Front Riggings - Basic Fixed Center Mount Platform 70/90/97 Degrees (FCMT70/ 90/ 97)	132

TABLE of CONTENTS

Description	Page Number
Front Riggings - Center Mount Legrest Power (LNX)	143
Front Riggings - Center Mount Rigging (CTMT-FJ/CTMT-F)	134
Front Riggings - Center Mount Rigging Mounting Hardware (CTMT-FJ/CTMT-F)	139
Front Riggings - Center Mount Rigging Support Hanger Brackets and Mounting Hardware (CTMT-FJ/F, FCMT70/ 90/ 97, LNX)	141
Front Riggings - Independent Front Rigging Cam Lock Receiver, Straight (M70HD/ MELR/ MPELR)	148
Front Riggings - Maxx Style Pivot Plus Elevating Legrest, Manual (MELR) and Power (MPELR)	150
Front Riggings - Maxx Style Pivot Plus Power Elevating Legrest (MPELR) - Phono Jack Adapter Box Kit	152
Front Riggings - Independent Front Rigging Cross Brace Assembly (M70HD/ MELR/ MPELR), (16" W - 22" W)	146
Light System	37
Motors - 4-Pole - SSD Gearbox and Mounting Hardware	39
Packaging - Literature (English) & Decals	9
Packaging - Literature (French-Canadian) & Decals	11
Shroud and Rim Insert Colors	7
Shrouds, Nylon	25
Shrouds, Polycarbonate	27
Shrouds, Rear Top and Lower Rear, Nylon	29
Shrouds, Rear Top and Lower Rear, Polycarbonate	31
Suspension Assembly	21
Ultra Low Maxx - Lift-Off (LHSFL/ LHSFR) and Fixed Lateral Hip Support Hardware - Right/Left	124
Ultra Low Maxx - Anterior/Posterior Tilt (5FAA/ 5SA)	59
Ultra Low Maxx - Arm Assembly, Dual Post, Quick Height Adjust, Reduced Extension (ADPAH-M)	92

TABLE of CONTENTS

Description	Page Number
Ultra Low Maxx - Arm Assembly, Dual Post, Quick Height Adjust, Reduced Extension (ADPAH-S)	90
Ultra Low Maxx - Arm Assembly, Dual Post, Quick Height Adjust, Reduced Extension (ADPAH-T)	94
Ultra Low Maxx - Arm Assembly, Recline Adjustable Height, Angle Cantilever, Reduced Extension (MXFB-C)	102
Ultra Low Maxx - Arm Assembly, Recline Adjustable Height, Angle Cantilever, Reduced Extension (MXFB-S)	100
Ultra Low Maxx - Arm Assembly, Recline, Flip Back, Adjustable Height, Reduced Extension (REFB-L/ T)	96
Ultra Low Maxx - Arm Assembly, Seat Mount Flip Back Cantilever, Reduced Extension (ULMFB)	98
Ultra Low Maxx - Arm Pad Assembly, Modular (MISK-LD/ RD/ LF/ RF; MSTEX-LD/ RD/ LF/ RF)	104
Ultra Low Maxx - Arm Pad Assembly, Standard (STAP-LF/ LD/RF/RD)	108
Ultra Low Maxx - Arm Pad Assembly, Waterfall , Desk (WFALL-RD/LD), Full Length (WFALL-RF/LF)	110
Ultra Low Maxx - Arm Spacer Kits (1/2" and 1") For Tilt or Recline Armrest (OAS1, OAS12) .	112
Ultra Low Maxx - Arm Trough Assembly, Ergonomic (ERGOL/ R; FHJL-LT/ RT)	106
Ultra Low Maxx - Backrest Recliner Actuator, Shroud (L-Cover) and Mounting Hardware (ULMR/ ULMTR)	78
Ultra Low Maxx - Backrest Recliner Actuator, Shroud and Mounting Hardware (ULMR/ ULMTR)	80
Ultra Low Maxx - Canes, Back and Back Mount Rails and Hardwae (ASBA-2/ ULMT)	63
Ultra Low Maxx - Elbow Block Bar and Clamp Assembly (MASL/ MASR)	114
Ultra Low Maxx - Elbow Block/ Trunk/ Lateral Pad Assembly (EBLT/ EBRT)	116
Ultra Low Maxx - Elbow Block/ Trunk/ Lateral Pad Assembly Fixed/Swing Away Hardware (FLATL/ FLATR/ LMSA/ RMSA)	118
Ultra Low Maxx - Headrest Assembly and Headrest Pads (MEHSP-6,10,14), (MEHOP-9, 12), (MEH-4P), (CAH)	84
Ultra Low Maxx - Headrest Extension Kit, Elite, Low Profile Recline	86
Ultra Low Maxx - Headrest Extension Kit, Reclining Rehab Back	88
Ultra Low Maxx - High Back Backrest Assembly (HB16/ 18/ 20" H)	68

TABLE of CONTENTS

Description	Page Number
Ultra Low Maxx - Lateral Hip Support Pads (RHS-XS/ S/ M/ L) (LHS-XS/ S/ M/ L)	126
Ultra Low Maxx - Matrix Fixed/ Swing Away Lateral Hardware, Standard and Offset, Elite TR Back (PBLSL/R) (PBOL-L/R) (SAEHLB/R) (OSAH/L/R)	122
Ultra Low Maxx - MAXX Style Quick Release Removable Hardware For Hip Supports (MQRHL/ MQRHR)	128
Ultra Low Maxx - Mini-Set Clamp Only Hardware Kit (ULMT/ ASBA-2/ SLINGBK/ 16-22" Wide)	65
Ultra Low Maxx - Recessed Planar Back Assembly	76
Ultra Low Maxx - Recline Sharps Covers Kit (ULMR/ ULMTR)	82
Ultra Low Maxx - Rehab Back Assembly (RHBK)	71
Ultra Low Maxx - Seat Plates Positioning Hardware (also mounts to Tilt Assembly) (ULMT/ AELM)	61
Ultra Low Maxx - Upholstery - Sling Back - Nylon (ASBA-2 w/ SLNGBK)	66
Ultra Low Maxx -Matrix Swing Away Lateral Hardware, Standard, Elite/Elite TR Back (SAEHLB/ SAEHBR)	120
Ultra Low Maxx-"ELOW" - Main System Mounting and Interface Hardware (ULMR/ASBA-2) .	57
Ultra Low Maxx-"ELOW" - Seat Tilt/Seat Tilt w/ Elevate, Main System Mounting and Interface Hardware (ULMT/ ULMTR/AELM)	55
Ultra Low Maxx-"STANDARD" - Main System Mounting and Interface Hardware (ULMR/ASBA-2)	53
Ultra Low Maxx-"STANDARD" - Seat Tilt/Seat Tilt w/ Elevate, Main System Mounting and Interface Hardware (ULMT/ ULMTR/AELM)	51
Wheels - Drive Wheels & Mounting Hardware (14" x 3"), Nylon	47
Wheels - Drive Wheels & Mounting Hardware (14" x 3"), Polycarbonate	49

Shroud and Rim Insert Colors

Shroud Colors

Shroud Color Number	Shroud Color Description	Shroud Material
164P	Standard Black - NLA	Nylon
165P	Anthracite Black - NLA	Nylon
166P	Mercury Blue - NLA	Nylon
167P	Bengal Red - NLA	Nylon
168P	Silver Star - NLA	Nylon
169P	Midnight Purple - NLA	Nylon
170P	Obsidian Black	Polycarbonate
171P	Ruby Red	Polycarbonate
172P	Peacock Blue	Polycarbonate
173P	Pink Flamingo	Polycarbonate
174P	Purple Passion	Polycarbonate
175P	Denim Navy	Polycarbonate
176P	Lime Twist	Polycarbonate
177P	Artic Blue	Polycarbonate


Touch Up paints are not available.

Shroud and Rim Insert Colors

Rim Insert Colors

Rim Insert Color Number	Rim Insert Color Description	Rim Insert Material
164PR	Standard Black - NLA	Nylon
165PR	Anthracite Black - NLA	Nylon
166PR	Mercury Blue - NLA	Nylon
167PR	Bengal Red - NLA	Nylon
168PR	Silver Star - NLA	Nylon
169PR	Midnight Purple - NLA	Nylon
170PR	Obsidian Black	Polycarbonate
171PR	Ruby Red	Polycarbonate
172PR	Peacock Blue	Polycarbonate
173PR	Pink Flamingo	Polycarbonate
174PR	Purple Passion	Polycarbonate
175PR	Denim Navy	Polycarbonate
176PR	Lime Twist	Polycarbonate
177PR	Artic Blue	Polycarbonate
Touch Up paints are not available.		

Packaging, Literature (English) and Decals




Packaging, Literature (English) and Decals

Item Number	Note	Part Number	Description	Quantity Per	
				Package	Assembly
1		1036027	Carton, Tray	1	2
2		1036028	Carton, Cover	1	1
3		1129590	Carton, Corner Reinforcement	1	4
4		40131X161	Insert, Plain Pad	1	3
	A	60124850	Decal, AVIVA, White	1	1
	A	60124851	Decal, AVIVA, Black	1	1
	A	60124852	Decal, AVIVA FX20, White	1	1

NOTE: A - Not Shown

When adjusting or replacing parts, please refer to your Owner's Manual for assistance. The following Owner's Manuals can be found at www.invacare.com/homecare under Literature Search for the following: 60101916 "AVIVA FX Base", 60128796 "Ultra Low Maxx by Motion Concepts Modular Positioning System", 60126082 "LiNX Control System REM500 Remote", 60126081 "LiNX Control System REM400 Touch Screen Remote", 60126083 "LiNX Control System Non-Expandable/Expandable Remotes". Service Manuals are: 60108795 "LiNX For AVIVA Control Systems", 60107431 "AVIVA FX Base", 60108797 "Ultra Low Maxx by Motion Concepts Modular Positioning System".

Packaging, Literature (French-Canadian) and Decals




Packaging, Literature (French-Canadian) and Decals

Item Number	Note	Part Number	Description	Quantity Per	
				Package	Assembly
1		1036027	Carton, Tray	1	2
2		1036028	Carton, Cover	1	1
3		1129590	Carton, Corner Reinforcement	1	4
4		40131X161	Insert, Plain Pad	1	3
	A	60124850	Decal, AVIVA, White	1	1
	A	60124851	Decal, AVIVA, Black	1	1
	A	60124852	Decal, AVIVA FX20, White	1	1

NOTE: A - Not Shown



When adjusting or replacing parts, please refer to your Owner's Manual for assistance. The following Owner's Manuals (French-Canadian) can be found at www.invacare.com/homecare under Literature Search for the following: 60108792 "AVIVA FX Base", 60126882 "Ultra Low Maxx by Motion Concepts Modular Positioning System", 60126821 "LiNX Control System REM500 Remote", 60126820 "LiNX Control System REM400 Touch Screen Remote", 60126822 "LiNX Control System Non-Expandable/Expandable Remotes".

Battery Charger & Battery Harnesses GP24



AVIVA FX Base Assembly
shown for reference only

Batteries shown for
reference only




CTE Battery Charger

SHEP_337A

Battery Charger & Battery Harnesses GP24

Item Number	Note	Part Number	Description	Quantity Per	
				Package	Assembly
1	A	60124371	Kit, Battery Harnesses w/ Hardware, GP24	1	1
2		60102598	Kit, Charger, Battery, 8 Amp. 115vac (Off Board) Single Mode Round Plug w/ AC Power Line Cord	1	1
3		60102582	Kit, Cord, AC Power, 115VAC, 3 Prong On-Board (6' 6-1/2" L)	1	1
NOTE: A - CTE charger is standard on-chair. Used on Gel Cell or AGM (Absorbed Glass Mat) batteries.					


Battery Box Assembly, Front and Rear Plates



Battery Box Assembly, Front and Rear Plates

Item Number	Note	Part Number	Description	Quantity Per	
				Package	Assembly
1	A	60124325	Kit, Plate, Battery Box Rear w/ Hardware	1	1
2			Plate, Battery Box Rear	1	1
3			Screw, Flat Head Cap Serrated (M6 x 1 x 16mm)	1	2
			Sheet, Insert Parts Service	1	1
4	B	60124326	Kit, Plate, Battery Box Front GP24 w/ Hardware	1	1
3			Screw, Flat Head Cap Serrated (M6 x 1 x 16mm)	1	2
5			Plate, Battery Box Front	1	1
6			Label, Warning Battery GP24 Terminal	1	1
7	C	60123251	Label, Warning Battery GP24 Wiring	1	1
			Sheet, Insert Parts Service	1	1
4			Kit, Plate, Battery Box Front GP24 w/ Hardware FC	1	1
3			Screw, Flat Head Cap Serrated (M6 x 1 x 16mm)	1	2
5			Plate, Battery Box Front	1	1
6			Label, Warning Battery GP24 Terminal	1	1
7			Label, Warning Battery GP24 Wiring FC	1	1
			Sheet, Insert Parts Service	1	1
NOTE: A - Includes items 2 and 3. B - Includes items 3, 5-7. Labeling in English. C - Includes items 3, 5-7. Labeling in French-Canadian.					



Battery Harness, Base, AVIVA FX



Battery Harness, Base, AVIVA FX

Item Number	Note	Part Number	Description	Quantity Per	
				Package	Assembly
1	A	60124375	Kit, Battery Base Harness w/ Hardware	1	1
2			Harness, Base Battery	1	1
3			Nut, Hex (M3.5 x 0.6)	1	6
4			Washer, Lock Int-Tooth (M3.5)	1	6
5			Screw, Phillips Pan Head (M3.5 x 0.6 x 20mm)	1	6
6			Clip, Wire Locking (5-9mm)	1	2
7			Grommet, Rubber	1	2
8			Bushing, Open/Close	1	2
			Sheet, Insert Parts Service	1	1
NOTE: A - Includes items 2-8					

Circuit Breaker 63A




Assembled view of Circuit Breaker Assembly
shown for reference

SHEP_336A

Circuit Breaker 63A

Item Number	Note	Part Number	Description	Quantity Per	
				Package	Assembly
1	A	60124377	Kit, Circuit Breaker, 63A	1	1
2			Circuit Breaker, 63A	1	1
3			Foam, PE	1	1
NOTE: A - Includes items 2 and 3					

Suspension Assembly





SHEP_320A

Suspension Assembly

Item Number	Note	Part Number	Description	Quantity Per	
				Package	Assembly
1	A	60124329	Kit, Shock, Damper w/ Hardware (Pair)	1	1
2			Shock, Damper	1	2
3			Screw, Socket Head Cap (M10 x 1.5 x 45mm)	1	2
4			Screw, Socket Head Cap (M10 x 1.5 x 35mm)	1	2
5			Washer, Flat (10.5 x 20 x 2mm)	1	4
			Sheet, Insert Parts Service	1	1
6	B	60124328	Kit, Bumper, Rubber w/ Hardware (Pair)	1	1
7			Bumper, Rubber (M6 x 1mm)	1	2
8			Locknut (M6)	1	2
			Sheet, Insert Parts Service	1	1
9	C	60124350	Kit, Spring, Suspension, Black (Pair)	1	1
			Spring, Suspension, Black	1	2
			Sheet, Insert Parts Service	1	1
NOTE: A - Includes items 2-5 B - Includes items 7 and 8 C - It is recommended to replace Suspension Springs in Pairs					


Anti-Tip Assembly



Anti-Tip Assembly

Item Number	Note	Part Number	Description	Quantity Per	
				Package	Assembly
1	A	60125828	Kit, Anti-Tip Arms with Hardware (Pair)	1	1
2			Anti-Tip Arm - Right	1	1
3			Anti-Tip Arm - Left	1	1
4			Screw, Socket Head Cap w/ Patch (M8 x 1.25 x 16mm)	1	4
5			Plate, Anti-Tip	1	2
6			Bushing, Flanged (22 x 18)	1	4
7			Insert, Shock Absorber	1	2
8			Stop, Anti-Tip Solid	1	2
9	B	60124330	Kit, Wheel, Anti-Tip w/ Hardware (5 x 1") (Pair)	1	1
10			Wheel Assembly, Composite, Urethane (5" x 1")	1	2
11			NLA - Washer, Steel (8.4mm x 16mm x 1.6mm)	1	4
12			Screw, Socket Button Head Cap (M8 x 1.25 x 40mm)	1	2
			Sheet, Insert Parts Service	1	1
NOTE: A - Includes items 2-8 B - Includes items 10-12					

Shrouds, Nylon



Shrouds, Nylon

Item Number	Note	Part Number	Description	Quantity Per	
				Package	Assembly
1	A	60124331-NLA	NLA - Kit, Shroud, Top w/ Labels and Hardware	1	1
4		See Chart	Item #4 Side Rail Shroud Kit	1	1
5		60124335	Kit, Shroud, Cap Front - Left	1	1
6		60124336	Kit, Shroud, Cap Front - Right	1	1
7	B	60124332-NLA	NLA - Kit, Shroud, Panel Side w/ Labels and Hardware - Left	1	1
10	C	60124333-NLA	NLA - Kit, Shroud, Panel Side w/ Labels and Hardware - Right	1	1
12		60124334-NLA	NLA - Kit, Shroud, Front w/ Decal	1	1
13		See Chart	Item #13 Rim Cap Kit For Caster and Drive Wheel	1	1
NOTE: A - Includes items 2 and 3 B - Includes items 8 and 9. Also includes the Switch Plate already attached for the Circuit Breaker (not shown). C - Includes items 9 and 11					


Item #4 Side Rail Shroud Kit

Description	Standard Black	Anthracite Black	Mercury Blue	Bengal Red	Silver Star	Midnight Purple
Kit, Shroud, Rail Side, Right or Left	60124337-NLA See Note(s) A	60124338-NLA See Note(s) A	60124339-NLA See Note(s) A	60124340-NLA See Note(s) A	60124341-NLA See Note(s) A	60124342-NLA See Note(s) A
NOTE: A - Side Rail Shroud can accommodate the right or left side. Nylon Colors are No Longer Available. Select from Polycarbonate Colors.						

Item #13 Rim Cap Kit For Caster and Drive Wheel

Description	Standard Black	Anthracite Black	Mercury Blue	Bengal Red	Silver Star	Midnight Purple
Rim Cap, Caster and Drive Wheel Kit	60124853-NLA See Note(s) A	60124854-NLA See Note(s) A	60124855-NLA See Note(s) A	60124856-NLA See Note(s) A	60124857-NLA See Note(s) A	60124858-NLA See Note(s) A
NOTE: A - Front Drive Wheel Center Rim Cap remains Black in color. Individual Rims and Caps are not sold separately. Nylon Colors are No Longer Available. Select from Polycarbonate Colors.						

Shrouds, Polycarbonate



Shrouds, Polycarbonate

Item Number	Note	Part Number	Description	Quantity Per	
				Package	Assembly
1	A	60130170	Kit, Shroud, Top w/ Labels and Hardware	1	1
4		See Chart	Item #4 Side Rail Shroud Kit	1	1
4		See Chart	Item #4 Side Rail Shroud Kit (continued)	1	1
5		60124335	Kit, Shroud, Cap Front - Left	1	1
6	B	60124336	Kit, Shroud, Cap Front - Right	1	1
7		60130279	Kit, Shroud, Panel Side w/ Labels and Hardware - Left	1	1
10		60130278	Kit, Shroud, Panel Side w/ Labels and Hardware - Right	1	1
12		60130234	Kit, Shroud, Front w/ Decal	1	1
NOTE: A - Includes items 2 and 3 B - Includes items 8 and 9. Also includes the Switch Plate already attached for the Circuit Breaker (not shown). C - Includes items 9 and 11					


Item #4 Side Rail Shroud Kit

Description	Obsidian Black	Ruby Red	Denim Navy	Peacock Blue	Purple Passion	Lime Twist
Kit, Shroud, Rail Side, Right or Left	60130280	60130281	60130282	60130283	60130341	60130342
	See Note(s) A	See Note(s) A	See Note(s) A	See Note(s) A	See Note(s) A	See Note(s) A
NOTE: A - Side Rail Shroud can accommodate the right or left side.						


Item #4 Side Rail Shroud Kit (continued)

Description	Pink Flamingo	Arctic Blue
Kit, Shroud, Rail Side, Right or Left	60130343	60130344
	See Note(s) A	See Note(s) A
NOTE: A - Side Rail Shroud can accommodate the right or left side.		

Shrouds, Rear Top and Lower Rear, Nylon



SHEP_323A




Shrouds, Rear Top and Lower Rear, Nylon

Item Number	Note	Part Number	Description	Quantity Per	
				Package	Assembly
1	A	See Chart	Item #1 Rear Top Shroud w/ Hardware Kit	1	1
6	B	60124360	Kit, Shroud, Rear Lower w/o Lights	1	1
7			Shroud, Rear Lower	1	1
8			Decal, Invacare Medallion, Reflective - Large	1	1
9			Reflector, Red	1	2
10			Reflector, Ring	1	2
11			Reflector, Red	1	2
12			Spacer, Reflector	1	2
13			Label, Warning Motor/SSD Gearbox Replacement	1	1
			Sheet, Insert Parts Service	1	1
14	C	60124349	Kit, Shroud, Rear Lower w/ Lights	1	1
7			Shroud, Rear Lower	1	1
8			Decal, Invacare Medallion, Reflective - Large	1	1
9			Reflector, Red	1	2
13			Label, Warning Motor/SSD Gearbox Replacement	1	1
15			Harness, Rear Lights, AVIVA FX	1	1
			Sheet, Insert Parts Service	1	1
NOTE: A - Includes items 2-5 B - Includes items 7-13 C - Includes items 7-9, 13 and 15.					


Item #1 Rear Top Shroud w/ Hardware Kit

Description	Standard Black	Antracite Black	Mercury Blue	Bengal Red	Silver Star	Midnight Purple
Rear Top Shroud w/ Hardware Kit	60124343-NLA	60124344-NLA	60124345-NLA	60124346-NLA	60124347-NLA	60124348-NLA
	See Note(s) A	See Note(s) A	See Note(s) A	See Note(s) A	See Note(s) A	See Note(s) A
NOTE: A - Includes items 2-5. Nylon Colors are No Longer Available. Select from Polycarbonate Colors.						

Shrouds, Rear Top and Lower Rear, Polycarbonate



SHEP_323A



Shrouds, Rear Top and Lower Rear, Polycarbonate

Item Number	Note	Part Number	Description	Quantity Per			
				Package	Assembly		
1	A	See Chart	Item #1 Rear Top Shroud w/ Hardware Kit	1	1		
1	A	See Chart	Item #1 Rear Top Shroud w/ Hardware Kit (continued)	1	1		
6	B	60130235	Kit, Shroud, Rear Lower w/o Lights	1	1		
7			Shroud, Rear Lower	1	1		
8			Decal, Invacare Medallion, Reflective - Large	1	1		
9			Reflector, Red	1	2		
10			Reflector, Ring	1	2		
11			Reflector, Red	1	2		
12			Spacer, Reflector	1	2		
13			Label, Warning Motor/SSD Gearbox Replacement	1	1		
			Sheet, Insert Parts Service	1	1		
14			C	60130236	Kit, Shroud, Rear Lower w/ Lights	1	1
7					Shroud, Rear Lower	1	1
8					Decal, Invacare Medallion, Reflective - Large	1	1
9	Reflector, Red	1			2		
13	Label, Warning Motor/SSD Gearbox Replacement	1			1		
15	Harness, Rear Lights, AVIVA FX	1			1		
			Sheet, Insert Parts Service	1	1		
NOTE: A - Includes items 2-5 B - Includes items 7-13 C - Includes items 7-9, 13 and 15.							


Item #1 Rear Top Shroud w/ Hardware Kit

Description	Obsidian Black	Ruby Red	Denim Navy	Peacock Blue	Purple Passion	Lime Twist
Rear Top Shroud w/ Hardware Kit	60130237	60130238	60130239	60130240	60130241	60130242
	See Note(s) A	See Note(s) A	See Note(s) A	See Note(s) A	See Note(s) A	See Note(s) A
NOTE: A - Includes items 2-5.						

Item #1 Rear Top Shroud w/ Hardware Kit (continued)

Description	Pink Flamingo	Arctic Blue
Rear Top Shroud w/ Hardware Kit	60130243	60130244
	See Note(s) A	See Note(s) A
NOTE: A - Includes items 2-5.		

Fender w/ Hardware, Drive Wheel, Nylon




Fender w/ Hardware, Drive Wheel, Nylon

Item Number	Note	Part Number	Description	Quantity Per	
				Package	Assembly
1	A	See Chart	Item #1 Fender w/ Hardware, Drive Wheel Kit (Pair)	2	1
6	A	60124351	Kit, Fender Brackets Right and Left (Pair)	2	1
NOTE: A - NOTE: The top two Motor Mounting Screws that attach the Fender Bracket are longer than the two bottom Motor Mounting Screws.					

Item #1 Fender w/ Hardware, Drive Wheel Kit (Pair)

Description	Standard Black	Anthracite Black	Mercury Blue	Bengal Red	Silver Star	Midnight Purple
Fender w/ Hardware, Drive Wheel Kit (Pair)	60124352-NLA See Note(s) A	60124353-NLA See Note(s) A	60124354-NLA See Note(s) A	60124355-NLA See Note(s) A	60124356-NLA See Note(s) A	60124357-NLA See Note(s) A
NOTE: A - Includes items 2-5 to accommodate both the Right and Left side. Nylon Colors are No Longer Available. Select from Polycarbonate Colors.						

Fender w/ Hardware, Drive Wheel, Polycarbonate



Fender w/ Hardware, Drive Wheel, Polycarbonate

Item Number	Note	Part Number	Description	Quantity Per	
				Package	Assembly
1	A	See Chart	Item #1 Fender w/ Hardware, Drive Wheel Kit	1	1
1	A	See Chart	Item #1 Fender w/ Hardware, Drive Wheel Kit (continued)	1	1
6	A	60124351	Kit, Fender Brackets Right and Left (Pair)	2	1
NOTE: A - NOTE: The top two Motor Mounting Screws that attach the Fender Bracket are longer than the two bottom Motor Mounting Screws.					


Item #1 Fender w/ Hardware, Drive Wheel Kit

Description	Obsidian Black	Ruby Red	Peacock Blue	Pink Flamingo	Purple Passion	Denim Navy
Fender w/ Hardware, Drive Wheel Kit	60129165	60129166	60129167	60129168	60129169	60129170
	See Note(s) A	See Note(s) A	See Note(s) A	See Note(s) A	See Note(s) A	See Note(s) A
NOTE: A - Includes items 2-5 to accommodate one side only.						

Item #1 Fender w/ Hardware, Drive Wheel Kit (continued)

Description	Lime Twist	Artic Blue
Fender w/ Hardware, Drive Wheel Kit	60129171	60129172
	See Note(s) A	See Note(s) A
NOTE: A - Includes items 2-5 to accommodate one side only.		

Light System




Light System

Item Number	Note	Part Number	Description	Quantity Per	
				Package	Assembly
1	A	60130236	Kit, Shroud, Rear Lower w/ Lights	1	1
			Shroud, Rear Lower	1	1
			Decal, Invacare Medallion, Reflective - Large	1	1
			Reflector, Red	1	2
			Label, Warning Motor/SSD Gearbox Replacement	1	1
			Harness, Rear Lights, AVIVA FX	1	1
			Sheet, Insert Parts Service	1	1
2		60124848	Kit, Cable, Power Light	1	1
3		60124382	Kit, Cable, Front Lights	1	1
4	B,C	60124884	Kit, Head Light Assembly w/ Hardware - Right	1	1
5			Light Assembly - Right	1	1
6			Locknut (M6)	1	2
7			Screw, Socket Button Head Cap (M6 x 1 x 18mm)	1	2
4	B,C	60124885	Kit, Head Light Assembly w/ Hardware - Left	1	1
5			Light Assembly - Left	1	1
6			Locknut (M6)	1	2
7			Screw, Socket Button Head Cap (M6 x 1 x 18mm)	1	2
NOTE: A - For exploded detail view of the Rear Lower Shroud w/ Lights see page titled "Shrouds, Top Rear and Lower Rear" B - Includes items 5-7 C - NOTE: The top two Motor Mounting Screws that attach the Fender Bracket and Light Assembly are longer than the two bottom Motor Mounting Screws.					


Motors, 4-Pole, SSD Gearbox and Mounting Hardware

When installing a new Motor or Motor/Gearbox Assembly, the Motor will need to be calibrated with Invacare Tech. Services. An LAK is required to complete the calibration.




4-Pole Motor/SSD Gearbox Assembly


4-Pole Motor Assembly



Motor/SSD Gearbox Band Clamp and Hardware Kit



SSD Gearbox Assembly



4-Pole Motor/SSD Gearbox Assembly, Mounting Bracket and Fender Assembly shown for reference only

SHEP_325A

Motors, 4-Pole, SSD Gearbox and Mounting Hardware

When installing a new Motor or Motor/Gearbox Assembly, the Motor will need to be calibrated with Invacare Tech. Services. An LAK is required to complete the calibration.

Item Number	Note	Part Number	Description	Quantity Per	
				Package	Assembly
1	A,E,F,H	60124863	Kit, 4-Pole 5K Motor/23:1 SSD Gearbox Assembly - Right	1	1
			4-Pole 5K Motor/23:1 SSD Gearbox Assembly - Right	1	1
			Decal, Motor Warning	1	1
			Label, Push/Drive	1	1
			Grease, Super Lube (1cc Tube)	1	1
			Sheet, Insert Parts Service	1	1
1	A,E,F,H	60124859	Kit, 4-Pole 5K Motor/23:1 SSD Gearbox Assembly - Left	1	1
			4-Pole 5K Motor/23:1 SSD Gearbox Assembly - Left	1	1
			Decal, Motor Warning	1	1
			Label, Push/Drive	1	1
			Grease, Super Lube (1cc Tube)	1	1
			Sheet, Insert Parts Service	1	1
1	A,E,H,I	60124887	Kit, 4Pole Low Speed Motor/18:1 Gearbox Assembly - Right	1	1
			4-Pole Low Speed Motor/18:1 SSD Gearbox Assembly - Right	1	1
			Decal, Motor Warning	1	1
			Label, Push/Drive	1	1
			Grease, Super Lube (1cc Tube)	1	1
			Sheet, Insert Parts Service	1	1
1	A,E,H,I	60124886	Kit, 4Pole Low Speed Motor/18:1 Gearbox Assembly - Left	1	1
			4-Pole Low Speed Motor/18:1 SSD Gearbox Assembly - Left	1	1
			Decal, Motor Warning	1	1
			Label, Push/Drive	1	1
			Grease, Super Lube (1cc Tube)	1	1
			Sheet, Insert Parts Service	1	1
6		1168883	Package, Coupling/Cup Assembly	10	1
2	B,E,F,H	60124867	Kit, Motor Assembly, 4-Pole 5K	1	1
3			Motor Assembly, 4-Pole 5K	1	1
4			Decal, Motor Warning	1	1
5			Label, Push/Drive	1	1
7			Clamp, Band	1	1
8			Locknut, DIN 985-M6- (A2K)-8	1	1
9			Washer (M6 6.6 x 12 x 1.8mm)	1	1
10			Spacer, Clamp, Molded	1	1
11			Bolt, Round Head Square Neck (M6 x 1 x 35mm)	1	1

Motors, 4-Pole, SSD Gearbox and Mounting Hardware

When installing a new Motor or Motor/Gearbox Assembly, the Motor will need to be calibrated with Invacare Tech. Services. An LAK is required to complete the calibration.

Item Number	Note	Part Number	Description	Quantity Per	
				Package	Assembly
2	B,E,H,I	60124869	Sheet, Insert Parts Service	1	1
3			Kit, Motor Assembly, 4-Pole Low Speed	1	1
4			Motor Assembly, 4-Pole Low Speed	1	1
5			Decal, Motor Warning	1	1
7			Label, Push/Drive	1	1
8			Clamp, Band	1	1
9			Locknut, DIN 985-M6- (A2K)-8	1	1
10			Washer (M6 6.6 x 12 x 1.8mm)	1	1
11			Spacer, Clamp, Molded	1	1
13			Bolt, Round Head Square Neck (M6 x 1 x 35mm)	1	1
14			Sheet, Insert Parts Service	1	1
14	D	60129163	Kit, Brake Lever Handle and Cap w/ Hardware (Pair)	1	1
15			Lever, Brake, Black - Right	1	1
16			Lever, Brake, Black - Left	1	1
17			Cap, Brake Lever, Motor	1	2
18			Screw, LG Phillips Head (M3.5 x 0.6 x 8mm)	1	2
19			NLA - Screw, Set Socket Cup Point w/ Patch (M5 x .8 x 6mm)	1	2
20			Screw, LG Torx Socket Head Cap SC (M4.5 x 1.59 x 10mm)	1	2
7	E	60124874	Sheet, Insert Parts Service	1	1
8			Kit, Motor/SSD Gearbox Band Clamp and Hardware	1	1
9			Clamp, Band	1	1
10			Locknut, DIN 985-M6- (A2K)-8	1	1
11			Washer (M6 6.6 x 12 x 1.8mm)	1	1
20			Spacer, Clamp, Molded	1	1
20			Bolt, Round Head Square Neck (M6 x 1 x 35mm)	1	1
20	C,E,F	60124872	Sheet, Insert Parts Service	1	1
20			Kit, SSD Gearbox Assembly, 23:1 - Right	1	1
20			Clamp, Band	1	1
20			Locknut, DIN 985-M6- (A2K)-8	1	1
20			Washer (M6 6.6 x 12 x 1.8mm)	1	1
20			Spacer, Clamp, Molded	1	1
20			Bolt, Round Head Square Neck (M6 x 1 x 35mm)	1	1
20			Gearbox Assembly, SSD, 23:1 - RWD/FWD-Left & CWD-Right	1	1
20			Sheet, Insert Parts Service	1	1
20			Kit, SSD Gearbox Assembly, 23:1 -Left	1	1
20	C,E,F	60124870	Clamp, Band	1	1
20			Locknut, DIN 985-M6- (A2K)-8	1	1
20			Washer (M6 6.6 x 12 x 1.8mm)	1	1
20			Washer (M6 6.6 x 12 x 1.8mm)	1	1


Motors, 4-Pole, SSD Gearbox and Mounting Hardware

When installing a new Motor or Motor/Gearbox Assembly, the Motor will need to be calibrated with Invacare Tech. Services. An LAK is required to complete the calibration.

Item Number	Note	Part Number	Description	Quantity Per	
				Package	Assembly
20	C,E,I	60124873	Spacer, Clamp, Molded	1	1
			Bolt, Round Head Square Neck (M6 x 1 x 35mm)	1	1
			Gearbox Assembly, SSD, 23:1 - RWD/FWD-Right & CWD-Left	1	1
			Sheet, Insert Parts Service	1	1
			Kit, SSD Gearbox Assembly, 18:1 - Right	1	1
			Gearbox Assembly, SSD, 18:1 - RWD/FWD-Left & CWD-Right	1	1
			Clamp, Band	1	1
			Locknut, DIN 985-M6- (A2K)-8	1	1
			Washer (M6 6.6 x 12 x 1.8mm)	1	1
			Spacer, Clamp, Molded	1	1
20	C,E,I	60124871	Bolt, Round Head Square Neck (M6 x 1 x 35mm)	1	1
			Sheet, Insert Parts Service	1	1
			Kit, SSD Gearbox Assembly, 18:1 - Left	1	1
			Gearbox Assembly, SSD, 18:1 - RWD/FWD-Right & CWD-Left	1	1
			Clamp, Band	1	1
			Locknut, DIN 985-M6- (A2K)-8	1	1
			Washer (M6 6.6 x 12 x 1.8mm)	1	1
			Spacer, Clamp, Molded	1	1
			Bolt, Round Head Square Neck (M6 x 1 x 35mm)	1	1
			Sheet, Insert Parts Service	1	1
21	G	60124358	Kit, Motor/SSD Gearbox Mounting Hardware	1	1
			Screw, Socket Head Cap (M8 x 1.25 x 65mm)	1	2
			Screw, Socket Head Cap (M8 x 1.25 x 60mm)	1	2
			Sheet, Insert Parts Service	1	1
		1180609	Package, Grease, Super Lube (1cc Tube) (Pkg. of 10)	10	1

NOTE: A - Includes items 3-12, Super Lube Grease and Instruction Sheet. NOTE: Brake Lever Item #13 is "not" included.
 B - Includes items 3-5, 7-11, Super Lube Grease and Instruction Sheet. NOTE: Brake Lever Item #13 is "not" included.
 C - Includes items 7-12, Super Lube Grease and Instruction Sheet.
 D - Includes items 14-18 and sold as a pair.
 E - Before placing an order, it is recommended to call Invacare's Technical Service Department at (1-800-832-4707) to verify if the Motor, SSD Gearbox or the Motor/SSD Gearbox Assembly is needed.
 F - Used for HSMTR High Speed Motor Option.
 G - Includes items 22 and 23
 H - When installing a new Motor or Motor/Gearbox Assembly, the Motor will need to be calibrated with Invacare Tech. Services. A LAK is required to complete the calibration.
 I - Used for Models without HSMTR High Speed Motors

Caster - 8" Rear - Single Sided Fork, Nylon




Caster - 8" Rear - Single Sided Fork, Nylon

Item Number	Note	Part Number	Description	Quantity Per	
				Package	Assembly
1	A	60124368	Kit, Single Sided Fork Assembly - Right	1	1
1		60124369	Kit, Single Sided Fork Assembly - Left	1	1
2		60124367	Kit, Caster Mounting Hardware	1	1
3			Screw, DIN7991SVZ TX PB (M8 x 20 x 12.9mm)	1	1
4	B		Spacer, Black (9 x 27 x 5.5mm)	1	1
5		60124366	Kit, Wheel Assembly, Flat Free, Black (8x2.50-4) Right and Left (Pair)	2	1
6		C	See Chart	Item #6 Rim Cap Kit For Caster and Drive Wheel, Nylon	1
10	D	1110709	Cap, Head Tube, Black	1	1
11		60124370	Kit, Fork Mounting Hardware	1	1
12			Locknut (M14 x 2mm)	1	1
13			Washer, Nylon (5/8 x 1-5/16 x 1/8")	1	1
14			Spacer, Stem (5/8 x 1-1/8 x 3/16")	1	1
15		40117X006	Bearing (5/8")	1	2
16		1110717	Ring, Tolerance	1	2
NOTE: A - Includes items 3 and 4. B - Right and Left side Wheel Assembly sold as pair only. C - Includes items 7-9. D - Includes items 12-16.					

Item #6 Rim Cap Kit For Caster and Drive Wheel, Nylon

Description	Standard Black	Anthracite Black	Mercury Blue	Bengal Red	Silver Star	Midnight Purple
Rim Cap, Caster and Drive Wheel Kit	60124853-NLA See Note(s) A	60124854-NLA See Note(s) A	60124855-NLA See Note(s) A	60124856-NLA See Note(s) A	60124857-NLA See Note(s) A	60124858-NLA See Note(s) A
NOTE: A - Front Drive Wheel Center Rim Cap remains Black in color. Individual Rims and Caps are not sold separately. Nylon Colors are No Longer Available. Select from Polycarbonate Colors.						

Caster - 8" Rear - Single Sided Fork, Polycarbonate



Caster - 8" Rear - Single Sided Fork, Polycarbonate

Item Number	Note	Part Number	Description	Quantity Per	
				Package	Assembly
1	A	60124368	Kit, Single Sided Fork Assembly - Right	1	1
1		60124369	Kit, Single Sided Fork Assembly - Left	1	1
2		60124367	Kit, Caster Mounting Hardware	1	1
3			Screw, DIN7991SVZ TX PB (M8 x 20 x 12.9mm)	1	1
4	C		Spacer, Black (9 x 27 x 5.5mm)	1	1
5		60124366	Kit, Wheel Assembly, Flat Free, Black (8x2.50-4) Right and Left (Pair)	2	1
6		See Chart	Item #6 Caster Rim Cap Kit, Polycarbonate	1	1
6		See Chart	Item #6 Caster Rim Cap Kit (continued), Polycarbonate	1	1
7	B	1110709	Cap, Head Tube, Black	1	1
8		60124370	Kit, Fork Mounting Hardware	1	1
9			Locknut (M14 x 2mm)	1	1
10			Washer, Nylon (5/8 x 1-5/16 x 1/8")	1	1
11			Spacer, Stem (5/8 x 1-1/8 x 3/16")	1	1
12		40117X006	Bearing (5/8")	1	2
13		1110717	Ring, Tolerance	1	2
NOTE: A - Includes items 3 and 4. B - Includes items 9-11. C - Right and Left side Wheel Assembly sold as pair only					


Item #6 Caster Rim Cap Kit, Polycarbonate

Description	Obsidian Black	Ruby Red	Peacock Blue	Pink Flamingo	Purple Passion	Denim Navy
Caster Rim Cap Kit	60129184 See Note(s) A	60129185 See Note(s) A	60129186 See Note(s) A	60129187 See Note(s) A	60129188 See Note(s) A	60129189 See Note(s) A
NOTE: A - Includes 1 ea. Caster Rim Cap						

Item #6 Caster Rim Cap Kit (continued), Polycarbonate

Description	Lime Twist	Artic Blue	Banana Split	Orange Marmalade	Morning Mist
Caster Rim Cap Kit	60129190 See Note(s) A	60129191 See Note(s) A	60129192 See Note(s) A	60129193 See Note(s) A	60129194 See Note(s) A
NOTE: A - Includes 1 ea. Caster Rim Cap					

Drive Wheels and Mounting Hardware (14" x 3"), Nylon



SHEP_328A


Drive Wheels and Mounting Hardware (14" x 3"), Nylon

Item Number	Note	Part Number	Description	Quantity Per	
				Package	Assembly
1	A	60124365	Kit, Drive Wheel, Flat Free w/ Hub Assembly, Black (14" x 3")	1	1
2	B	See Chart	Item #2 Rim Cap Kit For Caster and Drive Wheel, Nylon	1	1
6	C	60129154	Kit, Drive Wheel Mounting Hardware	1	1
7			Locknut (M12 x 1.5)	1	1
8			Washer, Belleville (12 x 23 x 1.5mm)	1	1
9			Key (6mm x 6mm x 40mm)	1	1
			Sheet, Insert Parts Service	1	1
9		1023264	Key (6mm x 6mm x 40mm)	1	1
NOTE: A - Sold Complete. B - Includes items 3-5. C - Includes items 7-9.					

Item #2 Rim Cap Kit For Caster and Drive Wheel, Nylon

Description	Standard Black	Anthracite Black	Mercury Blue	Bengal Red	Silver Star	Midnight Purple
Rim Cap, Caster and Drive Wheel Kit	60124853-NLA	60124854-NLA	60124855-NLA	60124856-NLA	60124857-NLA	60124858-NLA
	See Note(s) A	See Note(s) A	See Note(s) A	See Note(s) A	See Note(s) A	See Note(s) A
NOTE: A - Front Drive Wheel Center Rim Cap remains Black in color. Individual Rims and Caps are not sold separately. Nylon Colors are No Longer Available. Select from Polycarbonate Colors.						

Drive Wheels and Mounting Hardware (14" x 3"), Polycarbonate



Drive Wheels and Mounting Hardware (14" x 3"), Polycarbonate

Item Number	Note	Part Number	Description	Quantity Per	
				Package	Assembly
1	A	60124365	Kit, Drive Wheel, Flat Free w/ Hub Assembly, Black (14" x 3")	1	1
2		60129164	Kit, Drive Wheel Hub Cap, Black	1	1
3		See Chart	Item #3 Drive Wheel Rim Cap Kit, Polycarbonate	1	1
3		See Chart	Item #3 Drive Wheel Rim Cap Kit (continued), Polycarbonate	1	1
4	B	60129154	Kit, Drive Wheel Mounting Hardware	1	1
5			Locknut (M12 x 1.5)	1	1
6			Washer, Belleville (12 x 23 x 1.5mm)	1	1
7			Key (6mm x 6mm x 40mm)	1	1
			Sheet, Insert Parts Service	1	1
7		1023264	Key (6mm x 6mm x 40mm)	1	1
NOTE: A - Sold Complete. B - Includes items 5-7.					


Item #3 Drive Wheel Rim Cap Kit, Polycarbonate

Description	Obsidian Black	Ruby Red	Peacock Blue	Pink Flamingo	Purple Passion	Denim Navy
Drive Wheel Rim Cap Kit	60129173 See Note(s) A	60129174 See Note(s) A	60129175 See Note(s) A	60129176 See Note(s) A	60129177 See Note(s) A	60129178 See Note(s) A
NOTE: A - Includes 1 ea. Drive Wheel Rim Cap.						

Item #3 Drive Wheel Rim Cap Kit (continued), Polycarbonate

Description	Lime Twist	Artic Blue	Banana Split	Orange Marmalade	Morning Mist
Drive Wheel Rim Cap Kit	60129179 See Note(s) A	60129180 See Note(s) A	60129181 See Note(s) A	60129182 See Note(s) A	60129183 See Note(s) A
NOTE: A - Includes 1 ea. Drive Wheel Rim Cap.					

Ultra Low Maxx-"STANDARD" - Seat Tilt/Seat Tilt w/ Elevate, Main System Mounting and Interface Hardware (ULMT/ ULMTR/AELM)




Ultra Low Maxx-"STANDARD" - Seat Tilt/Seat Tilt w/ Elevate, Main System Mounting and Interface Hardware (ULMT/ ULMTR/AELM)

Item Number	Note	Part Number	Description	Quantity Per	
				Package	Assembly
1	A	60102606	Kit, Tilt Module Assembly w/ Actuator	1	1
2	A	60103390	Kit, Tilt/Elevate Module Assembly w/ Actuators	1	1
3		60102635	Kit, Bumper and Rear Wiring Cover w/ Mounting Hardware, Elevate	1	1
4	B	60124363	Kit, Seat Mounting Hardware, Standard	1	1
5			Screw, Socket Head Cap (M6 x 1 x 20mm)	1	12
6			Washer, Split Lock, Black Zinc (M6)	1	12
7			Washer, Plain, Black Zinc (1/4 x 9/16 x 1/16")	1	12
			Sheet, Insert Parts Service	1	1
8	C	60124361	Kit, Seat Mounting Bracket Hardware, Standard	1	1
9			Seat Rail Mounting Bracket - Right	1	1
10			Seat Rail Mounting Bracket - Left	1	1
11			Screw, Socket Head Cap w/ Patch (M8 x 1.25 x 16mm)	1	8
12			Bracket, Vertical Seat Mounting	1	4
			Sheet, Insert Parts Service	1	1

NOTE: A - Kit is sold as a complete assembly only.
B - Includes items 5-7
C - Includes items 9-12


Ultra Low Maxx-"STANDARD" - Main System Mounting and Interface Hardware (ULMR/ASBA-2)



Ultra Low Maxx-"STANDARD" - Main System Mounting and Interface Hardware (ULMR/ASBA-2)

Item Number	Note	Part Number	Description	Quantity Per	
				Package	Assembly
1	A	60124944	Kit, Non Power Seat Mounting Bracket	1	1
2		60124363	Kit, Seat Mounting Hardware, Standard	1	1
3			Screw, Socket Head Cap (M6 x 1 x 20mm)	1	12
4			Washer, Split Lock, Black Zinc (M6)	1	12
5			Washer, Plain, Black Zinc (1/4 x 9/16 x 1/16")	1	12
	B		Sheet, Insert Parts Service	1	1
6		60124361	Kit, Seat Mounting Bracket Hardware, Standard	1	1
7			Seat Rail Mounting Bracket - Right	1	1
8			Seat Rail Mounting Bracket - Left	1	1
9			Screw, Socket Head Cap w/ Patch (M8 x 1.25 x 16mm)	1	8
10			Bracket, Vertical Seat Mounting	1	4
			Sheet, Insert Parts Service	1	1
NOTE: A - Includes items 3-5 B - Includes items 7-10					


Ultra Low Maxx-"ELOW" - Seat Tilt/Seat Tilt w/ Elevate, Main System Mounting and Interface Hardware (ULMT/ ULMTR/AELM)



Ultra Low Maxx-"ELOW" - Seat Tilt/Seat Tilt w/ Elevate, Main System Mounting and Interface Hardware (ULMT/ ULMTR/AELM)

Item Number	Note	Part Number	Description	Quantity Per	
				Package	Assembly
1	A	60102606	Kit, Tilt Module Assembly w/ Actuator	1	1
2	A	60103390	Kit, Tilt/Elevate Module Assembly w/ Actuators	1	1
3		60102635	Kit, Bumper and Rear Wiring Cover w/ Mounting Hardware, Elevate	1	1
4	B	60124362	Kit, Seat Mounting Hardware, Low	1	1
5			Screw, Socket Head Cap (M6 x 1 x 20mm)	1	12
6			Washer, Split Lock, Black Zinc (M6)	1	12
			Sheet, Insert Parts Service	1	1
7	C	60124359	Kit, Seat Mounting Bracket Hardware, Low	1	1
8			Bracket, Seat Mounting, Low - Right	1	1
9			Bracket, Seat Mounting, Low - Left	1	1
10			Screw, Socket Head Cap w/ Patch (M8 x 1.25 x 16mm)	1	4
			Sheet, Insert Parts Service	1	1
NOTE: A - Kit is sold as a complete assembly only. B - Includes items 5 and 6 C - Includes items 8-10					


Ultra Low Maxx-"ELOW" - Main System Mounting and Interface Hardware (ULMR/ASBA-2)




Ultra Low Maxx-"ELOW" - Main System Mounting and Interface Hardware (ULMR/ASBA-2)

Item Number	Note	Part Number	Description	Quantity Per	
				Package	Assembly
1	A	60124944	Kit, Non Power Seat Mounting Bracket	1	1
2		60124362	Kit, Seat Mounting Hardware, Low	1	1
3			Screw, Socket Head Cap (M6 x 1 x 20mm)	1	12
4			Washer, Split Lock, Black Zinc (M6)	1	12
	B		Sheet, Insert Parts Service	1	1
5		60124359	Kit, Seat Mounting Bracket Hardware, Low	1	1
6			Bracket, Seat Mounting, Low - Right	1	1
7			Bracket, Seat Mounting, Low - Left	1	1
8			Screw, Socket Head Cap w/ Patch (M8 x 1.25 x 16mm)	1	4
			Sheet, Insert Parts Service	1	1
NOTE: A - Includes items 3 and 4 B - Includes items 6-8					

Ultra Low Maxx - Anterior/Posterior Tilt (5FAA/ 5SA)



Ultra Low Maxx Seating
shown for reference only




1588265.iso
SHEP_176A

Ultra Low Maxx - Anterior/Posterior Tilt (5FAA/ 5SA)

Item Number	Note	Part Number	Description	Quantity Per	
				Package	Assembly
1		60102637	Kit, Pre-Tilt Mounting 5 Deg., Anterior/Posterior-Tilt	1	1


Ultra Low Maxx - Seat Plates Positioning Hardware (ULMT/ AELM) (also mounts to Tilt Assembly)



Ultra Low Maxx - Seat Plates Positioning Hardware (ULMT/ AELM) (also mounts to Tilt Assembly)

Item Number	Note	Part Number	Description	Quantity Per	
				Package	Assembly
1		60102638	Kit, Seat Plates Positioning Hardware	1	1

Canes, Back and Back Mount Rails and Hardware (ASBA-2/ ULMT)




1641966.iso
SHEP_207A

Canes, Back and Back Mount Rails and Hardware (ASBA-2/ ULMT)

Item Number	Note	Part Number	Description	Quantity Per	
				Package	Assembly
1		60102595	Kit, Rehab Back Mount Rails and Hardware (Pair)	1	1
2		60102640	Kit, Back Cane Mounting Hardware (Pair)	1	1
3	A	60103415	Kit, Back Cane Assembly, 5 Degree Angle Bend (20")	1	1
3	A	60103416	Kit, Back Cane Assembly, 5 Degree Angle Bend (22")	1	1
3	A	60103417	Kit, Back Cane Assembly, 5 Degree Angle Bend (24")	1	1
4		60103411	Kit, Back Cane Assembly, Straight (18")	1	1
4		60103412	Kit, Back Cane Assembly, Straight (20")	1	1
4		60103413	Kit, Back Cane Assembly, Straight (22")	1	1
4		60103414	Kit, Back Cane Assembly, Straight (24")	1	1
NOTE: A - Not used for model pick SLINGBK Sling Back Upholstery.					


Mini-Set Clamp Only Hardware Kit (ULMT/ ASBA-2/ SLINGBK/ 16-22" Wide)



SHEP_224A

Item Number	Note	Part Number	Description	Quantity Per	
				Package	Assembly
1		60102675	Kit, Mini-Set Clamp Only Hardware (Pair)	1	1

Sling Back Upholstery - Nylon (ASBA-2 w/ SLNGBK)



SHEP_195A

Item Number	Note	Part Number	Description	Quantity Per	
				Package	Assembly
		See Chart	Upholstery Kit - Sling Back - Nylon w/ Mounting Hardware	1	1
	A	1026298	Package, Screw with Washer (#10 x 7/8") (Pair)	2	1

NOTE: A - Not Shown


Sling Back Upholstery - Nylon

(ASBA-2 w/ SLNGBK)

Upholstery Kit - Sling Back - Nylon w/ Mounting Hardware

Width	18" - 19" High	20" High	21" - 22" High	23" - 24" High
16" W	60105204	60105205	60105206	60105207
17" W	60105208	60105209	60105210	60105211
18" W	60105212	60105213	60105214	60105215
19" W	60105216	60105217	60105218	60105219
20" W	60105220	60105221	60105222	60105223
21" W	60105224	60105225	60105226	60105227
22" W	60105228	60105229	60105230	60105231

Ultra Low Maxx - High Back Backrest Assembly (HB16/ 18/ 20" H)



1588261.iso
SHEP_192A

Ultra Low Maxx - High Back Backrest Assembly (HB16/ 18/ 20" H)

Item Number	Note	Part Number	Description	Quantity Per	
				Package	Assembly
1		See Chart	Item #1 High Back Backrest Assembly 16"-18" Wide	1	1
1		See Chart	Item #1 High Back Backrest Assembly 19"-22" Wide	1	1
2		See Chart	Item #2 High Back Back Cover and foam 16"-18" Wide	1	1
2		See Chart	Item #2 High Back Back Cover and foam 19"-22" Wide	1	1
3		See Chart	Item #3 Privacy Flap	1	1

Item #1 High Back Backrest Assembly 16"-18" Wide

Overall Back Height	Material Back Height	16" Wide	17" Wide	18" Wide
TH20, TH21, TH22 TRH22, TRH23	16"	60102571	60102572	60102573
TH23 TRH24, TRH25	18"	60102782	60102783	60102784
TH24, TH25 TRH26, TRH27	20"	60102775	60102776	60102777

Item #1 High Back Backrest Assembly 19"-22" Wide

Overall Back Height	Material Back Height	19" Wide	20" Wide	21" Wide	22" Wide
TH20, TH21, TH22 TRH22, TRH23	16"	60102574	60102786	60102787	60102788
TH23 TRH24, TRH25	18"	60102785	60102779	60102780	60102781
TH24, TH25 TRH26, TRH27	20"	60102778	60102772	60102773	60102774

Ultra Low Maxx - High Back Backrest Assembly (HB16/ 18/ 20" H)

Item #2 High Back Back Cover and foam 16"-18" Wide

Overall Back Height	Material Back Height	16" Wide	17" Wide	18" Wide
TH20, TH21, TH22 TRH22, TRH23	16"	60102641	60102644	60102647
TH23 TRH24, TRH25	18"	60102642	60102645	60102648
TH24, TH25 TRH26, TRH27	20"	60102643	60102646	60102649

NOTE: Back cover and foam are sold complete.

Item #2 High Back Back Cover and foam 19"-22" Wide


Overall Back Height	Material Back Height	19" Wide	20" Wide	21" Wide	22" Wide
TH20, TH21, TH22 TRH22, TRH23	16"	60102650	60102653	60102656	60102659
TH23 TRH24, TRH25	18"	60102651	60102654	60102657	60102660
TH24, TH25 TRH26, TRH27	20"	60102652	60102655	60102658	60102661

NOTE: Back cover and foam are sold complete.

Item #3 Privacy Flap

Chair Width	Material Width	Description	Part Number
16-17" Wide	13" Wide	Privacy Flap (13" W x 16" H)	PFLAP-S
18-19" Wide	15" Wide	Privacy Flap (15" W x 16" H)	PFLAP-M
20-21" Wide	17" Wide	Privacy Flap (17" W x 16" H)	PFLAP-L
22" Wide	19" Wide	Privacy Flap (19" W x 16" H)	Currently Not Available

Ultra Low Maxx - Rehab Back Assembly (RHBK)



Ultra Low Maxx - Rehab Back Assembly (RHBK)

Item Number	Note	Part Number	Description	Quantity Per	
				Package	Assembly
1		See Chart	Item #1 Rehab Back Assembly 16"-18" Wide	1	1
1		See Chart	Item #1 Rehab Back Assembly 19"-22" Wide	1	1
2		See Chart	Item #2 Rehab Back Cover and Foam 16"-18" Wide	1	1
2		See Chart	Item #2 Rehab Back Cover and Foam 19"-22" Wide	1	1
3		See Chart	Item #3 Privacy Flap	1	1

Ultra Low Maxx - Rehab Back Assembly (RHBK)

Item #1 Rehab Back Assembly 16"-18" Wide

Overall Back Height	Material Back Height	16" Wide	17" Wide	18" Wide
TH18, TH19 TRH20, TRH21	14"	60102607	60102608	60102609
TH20, TH21 TRH22, TRH23	16"	60102612	60102613	60102614
TH22, TH23 TRH24, TRH25	18"	60102617	60102618	60102619
TH24, TH25 TRH26, TRH27	20"	60102622	60102623	60102624

Item #1 Rehab Back Assembly 19"-22" Wide

Overall Back Height	Material Back Height	19" Wide	20" Wide	21" Wide	22" Wide
TH18, TH19 TRH20, TRH21	14"	60102610	60102611	60102789	60102790
TH20, TH21 TRH22, TRH23	16"	60102615	60102616	60102795	60102796
TH22, TH23 TRH24, TRH25	18"	60102620	60102621	60102791	60102792
TH24, TH25 TRH26, TRH27	20"	60102625	60102626	60102793	60102794

Ultra Low Maxx - Rehab Back Assembly (RHBK)

Item #2 Rehab Back Cover and Foam 16"-18" Wide

Overall Back Height	Material Back Height	16" Wide	17" Wide	18" Wide
TH18, TH19 TRH20, TRH21	14"	60102851	60102852	60102853
TH20, TH21 TRH22, TRH23	16"	60102858	60102859	60102860
TH22, TH23 TRH24, TRH25	18"	60102865	60102866	60102867
TH24, TH25 TRH26, TRH27	20"	60102872	60102873	60102874

NOTE: Back cover and Foam are sold complete.

Item #2 Rehab Back Cover and Foam 19"-22" Wide

Overall Back Height	Material Back Height	19" Wide	20" Wide	21" Wide	22" Wide
TH18, TH19 TRH20, TRH21	14"	60102854	60102855	60102856	60102857
TH20, TH21 TRH22, TRH23	16"	60102861	60102862	60102863	60102864
TH22, TH23 TRH24, TRH25	18"	60102868	60102869	60102870	60102871
TH24, TH25 TRH26, TRH27	20"	60102875	60102876	60102877	60102878


NOTE: Back cover and foam are sold complete.

Ultra Low Maxx - Rehab Back Assembly (RHBK)

Item #3 Privacy Flap

Chair Width	Material Width	Description	Part Number
16-17" Wide	13" Wide	Privacy Flap (13" W x 16" H)	PFLAP-S
18-19" Wide	15" Wide	Privacy Flap (15" W x 16" H)	PFLAP-M
20-21" Wide	17" Wide	Privacy Flap (17" W x 16" H)	PFLAP-L
22" Wide	19" Wide	Privacy Flap (19" W x 16" H)	Currently Not Available

Recessed Planar Back Assembly




Recessed Planar Back Assembly

Item Number	Note	Part Number	Description	Quantity Per	
				Package	Assembly
1		60105189	Recessed Planar Back Assembly (16" Wide)	1	1
1		60105190	Recessed Planar Back Assembly (17" Wide)	1	1
1		60105191	Recessed Planar Back Assembly (18" Wide)	1	1
1		60105192	Recessed Planar Back Assembly (19" Wide)	1	1
1		60105193	Recessed Planar Back Assembly (20" Wide)	1	1
1		60105194	Recessed Planar Back Assembly (21" Wide)	1	1
1		60105195	Recessed Planar Back Assembly (22" Wide)	1	1
2		See Chart	Item #2 Privacy Flap	1	1

#2 Privacy Flap

Chair Width	Material Width	Description	Part Number
16-17" Wide	13" Wide	Privacy Flap (13" W x 16" H)	PFLAP-S
18-19" Wide	15" Wide	Privacy Flap (15" W x 16" H)	PFLAP-M
20-21" Wide	17" Wide	Privacy Flap (17" W x 16" H)	PFLAP-L
22" Wide	19" Wide	Privacy Flap (19" W x 16"H)	Currently Not Available

Ultra Low Maxx - Backrest Recliner Actuator, Shroud (L-Cover) and Mounting Hardware (ULMR/ ULMTR) (When ASL Options were ordered On-Chair)




Ultra Low Maxx - Backrest Recliner Actuator, Shroud (L-Cover) and Mounting Hardware (ULMR/ ULMTR) (When ASL Options were ordered On-Chair)

Item Number	Note	Part Number	Description	Quantity Per	
				Package	Assembly
1	A	60123530	Kit, Recline Actuator w/ Slotted Spring Pin	1	1
2		60102672	Kit, Shroud, Backrest Recline	1	1
3		60102721	Kit, Knob, Hand Wheel (M6 x 15)	1	1
4		60102546	Kit, Adapter Cable Assembly, DX2-ACT4 To Molex 2 Pin	1	1
			Adapter Cable Assembly, DX2-ACT4 To Molex 2 Pin	1	1
			Sheet, Insert Parts Service	1	1
NOTE: A - Use existing Recline Actuator Mounting Hardware when attaching the Recline Actuator.					

Ultra Low Maxx - Backrest Recliner Actuator, Shroud and Mounting Hardware

(ULMR/ ULMTR) (When Non-ASL Options were Ordered On-Chair)




Ultra Low Maxx - Backrest Recliner Actuator, Shroud and Mounting Hardware


(ULMR/ ULMTR) (When Non-ASL Options were Ordered On-Chair)

Item Number	Note	Part Number	Description	Quantity Per	
				Package	Assembly
1	A	60123530	Kit, Recline Actuator w/ Slotted Spring Pin	1	1
2		60125831	Kit, Recline Shroud	1	1
3	B	60125830	Kit, Backrest Bracket Mounting Hardware, Recline for DLX-ACT400/ GLM-CONX4	1	1
4		60102546	Kit, Adapter Cable Assembly, DX2-ACT4 To Molex 2 Pin	1	1
			Adapter Cable Assembly, DX2-ACT4 To Molex 2 Pin	1	1
			Sheet, Insert Parts Service	1	1
NOTE: A - Use existing Recline Actuator Mounting Hardware when attaching the Recline Actuator. B - Includes all hardware shown including, LiNX Electronics Buss Bracket, Recline Shroud/Buss Bracket and Tie Wraps (7 ea.)					

Recline Sharps Covers Kit (ULMR/ ULMTR)



Recline Seating System
shown for reference only




SHEP_225A

Recline Sharps Covers Kit (ULMR/ ULMTR)

Item Number	Note	Part Number	Description	Quantity Per	
				Package	Assembly
1		60102639	Kit, Recline Sharps Covers	1	1

Headrest Assembly and Headrest Pads (MEHSP-6,10,14), (MEHOP-9, 12), (MEH-4P), (CAH)




SHEP_229A

Headrest Assembly and Headrest Pads (MEHSP-6,10,14), (MEHOP-9, 12), (MEH-4P), (CAH)

Item Number	Note	Part Number	Description	Quantity Per	
				Package	Assembly
1	A	60102835	Kit, Motion Concepts Multi-Axis Headrest Mounting Hardware	1	1
2	B	60102826	Kit, Elan Multi-Axis Articulating Headrest Mounting Hardware	1	1
3	D	60103328	Kit, Concepts Auto Style Headrest Pad and Mounting Hardware	1	1
4	E	60102825	Kit, Matrx Elan Standard Headrest Pad (6" W)	1	1
4	F	60102830	Kit, Matrx Elan Standard Headrest Pad (10" W)	1	1
4	G	60102833	Kit, Matrx Elan Standard Headrest Pad (14" W)	1	1
5	H	60102834	Kit, Matrx Elan Occipital All Positioning Headrest Pad (9" W)	1	1
5	I	60102836	Kit, Matrx Elan Occipital All Positioning Headrest Pad (12" W)	1	1
6	J	60102829	Kit, Matrx Elan 4-Point Headrest Pad (11" W x 10" H)	1	1
NOTE: A - For model CMAH B - For model MEMAH D - For model CAH E - For model MEHSP-6 F - For model MEHSP-10 G - For model MEHSP-14 H - For model MEHOP-9 I - For model MEHOP-12 J - For model MEH-4P					


Headrest Extension Kit, Elite, Low Profile Recline



Headrest Extension Kit, Elite, Low Profile Recline

Item Number	Note	Part Number	Description	Quantity Per	
				Package	Assembly
1		60125833	Headrest Extension, Elite, Low Profile Recline	1	1
2			Plate, Headrest Extension	1	1
3			Locknut (M6)	1	2
4			Screw, Button Socket Cap (M6-1.0 x 16mm)	1	2


Headrest Extension Kit, Reclining Rehab Back



Headrest Extension Kit, Reclining Rehab Back

Item Number	Note	Part Number	Description	Quantity Per	
				Package	Assembly
1	A	60125832	Kit, Headrest Extension, Reclining Rehab Back	1	1
2			Shim	1	1
3			Plate, Headrest Extension	1	1
4			Screw, Button Socket Cap (M6-1.0 x 20mm)	1	2
5			Locknut (M6)	1	2
NOTE: A - Includes items 2-5					

Arm Assembly, Dual Post, Quick Height Adjust, Reduced Extension (ADPAH-S)




**Dual Post Arm Assembly,
Quick Height Adjust,
Reduced Extension (Pair)
(9-13" H)
ADPAH-S**

SHEP_240A

Arm Assembly, Dual Post, Quick Height Adjust, Reduced Extension (ADPAH-S)

Item Number	Note	Part Number	Description	Quantity Per	
				Package	Assembly
1	A	60102665	Kit, Arm Assembly, Dual Post, Quick Height Adjust, Reduced Extension (Pair) (9-13" H)	2	1
2	B	60102711	Kit, Arm Wrap w/ Pouch, Dual Post Adjustable Hight Armrest	1	1
NOTE: A - Arms include the Right and Left side Arm Assembly and Mounting Hardware B - Kit contains only one Arm Wrap.					

Arm Assembly, Dual Post, Quick Height Adjust, Reduced Extension (ADPAH-M)




SHEP_257A

Arm Assembly, Dual Post, Quick Height Adjust, Reduced Extension (ADPAH-M)

Item Number	Note	Part Number	Description	Quantity Per	
				Package	Assembly
1	A	60102663	Kit, Arm Assembly, Dual Post, Quick Height Adjust, Reduced Extension (Pair) (12-3/4-15-3/4" H)	2	1
2	B	60102711	Kit, Arm Wrap w/ Pouch, Dual Post Adjustable Hight Armrest	1	1
NOTE: A - Arms include the Right and Left side Arm Assembly and Mounting Hardware B - Kit contains only one Arm Wrap.					

Arm Assembly, Dual Post, Quick Height Adjust, Reduced Extension (ADPAH-T)




Dual Post Arm Assembly, Quick Height Adjust, Reduced Extension (Pair) (15-3/4-18-3/4" H) ADPAH-T

SHEP_258A

Arm Assembly, Dual Post, Quick Height Adjust, Reduced Extension (ADPAH-T)

Item Number	Note	Part Number	Description	Quantity Per	
				Package	Assembly
1	A	60102662	Kit, Arm Assembly, Dual Post, Quick Height Adjust, Reduced Extension (Pair) (15-3/4-18-3/4" H)	2	1
2	B	60102711	Kit, Arm Wrap w/ Pouch, Dual Post Adjustable Hight Armrest	1	1
NOTE: A - Arms include the Right and Left side Arm Assembly and Mounting Hardware B - Kit contains only one Arm Wrap.					


Ultra Low Maxx - Arm Assembly, Recline, Flip Back, Adjustable Height, Reduced Extension (REFB-L/ T)



Ultra Low Maxx - Arm Assembly, Recline, Flip Back, Adjustable Height, Reduced Extension (REFB-L/ T)

Item Number	Note	Part Number	Description	Quantity Per	
				Package	Assembly
1	A	60102666	Kit, Arm Assembly, Recline, Flip Back, Adjustable Height, Reduced Extension (Pair) (9-13" H)	2	1
2	A	60102715	Kit, Arm Assembly, Recline, Flip Back, Adjustable Height, Reduced Extension (Pair) (13-16" H)	2	1
NOTE: A - Arms include the Right and Left side Arm Assembly and Mounting Hardware					


Arm Assembly, Seat Mount Flip Back Cantilever, Reduced Extension (ULMFB)



Arm Assembly, Seat Mount Flip Back Cantilever, Reduced Extension (ULMFB)

Item Number	Note	Part Number	Description	Quantity Per	
				Package	Assembly
1	A	60102664	Kit, Arm Assembly, Seat Mount Flip Back Cantilever, Reduced Extension (Pair) (9-12" H)	2	1
NOTE: A - Arms include the Right and Left side Arm Assembly and Mounting Hardware					

Arm Assembly, Recline Adjustable Height, Angle Cantilever, Reduced Extension (MXFB-S)




Reclining Arm Assembly, Flip Back, Adjustable Height, Straight, Reduced Extension (9-13" H) MXFB-S

Arm Assembly, Recline Adjustable Height, Angle Cantilever, Reduced Extension (MXFB-S)

Item Number	Note	Part Number	Description	Quantity Per	
				Package	Assembly
1	A	60102667	Arm Assembly, Recline Adjustable Height, Straight, Angle Cantilever, Reduced Extension (Pair) (9-13" H)	2	1
NOTE: A - Arms include the Right and Left side Arm Assembly and Mounting Hardware					

Arm Assembly, Recline Adjustable Height, Angle Cantilever, Reduced Extension (MXFB-C)




Reclining Arm Assembly, Flip Back, Adjustable Height, Curved, Reduced Extension (12-1/2-16" H) MXFB-C

Arm Assembly, Recline Adjustable Height, Angle Cantilever, Reduced Extension (MXFB-C)

Item Number	Note	Part Number	Description	Quantity Per	
				Package	Assembly
1	A	60102714	Kit, Arm Assembly, Recline Adjustable Height, Curved, Angle Cantilever, Reduced Extension (Pair) (12-1/2-16" H)	2	1
NOTE: A - Arms include the Right and Left side Arm Assembly and Mounting Hardware					

Arm Pad Assembly, Modular

(MISK-LD/ RD/ LF/ RF; MSTEEX-LD/ RD/ LF/ RF)



1588505.iso
SHEP_246A

Arm Pad Assembly, Modular


(MISK-LD/ RD/ LF/ RF; MSTEEX-LD/ RD/ LF/ RF)

Item Number	Note	Part Number	Description	Quantity Per	
				Package	Assembly
1-2		See Chart	Item #1-2 Arm Pad Assembly, Modular	1	1
3		60102710	Kit, Desk to Full Armrest Extension	1	1

Item #1-2 Arm Pad Assembly, Modular

Item Number	Model #	Part Number	Description
1	MISK-LD	60102752	Modular Arm Pad Assembly, Integral Soft Skin Foam, Desk Length - Left
1	MISK-RD	60102754	Modular Arm Pad Assembly, Integral Soft Skin Foam, Desk Length - Right
1	MSTEEX-LD	60102753	Modular Arm Pad Assembly, Startex Covered Visco Foam, Desk Length - Left
1	MSTEEX-RD	60102755	Modular Arm Pad Assembly, Startex Covered Visco Foam, Desk Length - Right
2	MISK-LF	60102756	Modular Arm Pad Assembly, Integral Soft Skin Foam, Full Length - Left
2	MISK-RF	60102758	Modular Arm Pad Assembly, Integral Soft Skin Foam, Full Length - Right
2	MSTEEX-LF	60102757	Modular Arm Pad Assembly, Startex Covered Visco Foam, Full Length - Left
2	MSTEEX-RF	60102759	Modular Arm Pad Assembly, Startex Covered Visco Foam, Full Length - Right


Arm Trough Assembly, Ergonomic (ERGOL/ R; FHJL-LT/ RT)



Arm Trough Assembly, Ergonomic (ERGOL/ R; FHJL-LT/ RT)

Item Number	Note	Part Number	Description	Quantity Per	
				Package	Assembly
1		60102766	Kit, Arm Trough, Ergonomic - Left	1	1
1		60102767	Kit, Arm Trough, Ergonomic - Right	1	1
2		60102710	Kit, Desk to Full Armrest Extension	1	1
3		60102768	Kit, Flat Hand Pad	1	1


Arm Pad Assembly, Standard (STAP-LF/ LD/RF/RD)



Arm Pad Assembly, Standard (STAP-LF/ LD/RF/RD)

Item Number	Note	Part Number	Description	Quantity Per	
				Package	Assembly
1		60102761	Kit, Arm Pad w/ Screws, Full, Standard	1	1
2		60102710	Kit, Desk to Full Armrest Extension	1	1
3		60102760	Kit, Arm Pad w/ Screws, Desk, Standard	1	1



Arm Pad Assembly, Waterfall , Desk (WFALL-RD/LD), Full Length (WFALL-RF/LF)



Arm Pad Assembly, Waterfall , Desk (WFALL-RD/LD), Full Length (WFALL-RF/LF)

Item Number	Note	Part Number	Description	Quantity Per	
				Package	Assembly
1	A	60106010	Kit, Arm Pad w/ Screws, Full, Waterfall	1	1
2		60102710	Kit, Desk to Full Armrest Extension	1	1
3	A	60106009	Kit, Arm Pad w/ Screws, Desk, Waterfall	1	1
NOTE: A - Arm Pad can assemble to the right or left side of the wheelchair.					

Arm Spacer Kits (1/2" and 1") For Tilt or Recline Armrest (OAS1, OAS12)




Arm Spacer Kits (1/2" and 1") For Tilt or Recline Armrest (OAS1, OAS12)

Item Number	Note	Part Number	Description	Quantity Per	
				Package	Assembly
1		60105188	Kit, Offset Spacer, Dual Post (Pair) (1/2")	2	1
2		60105187	Kit, Offset Spacer, Dual Post (Pair) (1")	2	1
3		60105186	Kit, Offset Spacer, Recline (Pair) (1/2")	2	1
4		60105185	Kit, Offset Spacer, Recline (Pair) (1")	2	1

Elbow Block Bar and Clamp Assembly (MASL/ MASR)

Elbow Block Clamp used for Elbow Block Bar Hardware models
(VEXSHEXS-L, VEXSHEXL-L, VEXLHEXS-L, VEXLHEXL-L and
VEXSHEXS-R, VEXSHEXL-R, VEXLHEXS-R, VEXLHEXL-R)



SHEP_249A


Elbow Block Bar and Clamp Assembly (MASL/ MASR)

Elbow Block Clamp used for Elbow Block Bar Hardware models
(VEXSHEXS-L, VEXSHEXL-L, VEXLHEXS-L, VEXLHEXL-L and
VEXSHEXS-R, VEXSHEXL-R, VEXLHEXS-R, VEXLHEXL-R)


Item Number	Note	Part Number	Description	Quantity Per	
				Package	Assembly
1		60102688	Kit, Elbow Block Clamp Assembly	1	1
2		60105133	Kit, Elbow Block Bar Mounting (3" Vertical x 5" Horizontal)	1	1
2		60105134	Kit, Elbow Block Bar Mounting (3" Vertical x 7" Horizontal)	1	1
2		60105135	Kit, Elbow Block Bar Mounting (5" Vertical x 5" Horizontal)	1	1
2		60105136	Kit, Elbow Block Bar Mounting (5" Vertical x 7" Horizontal)	1	1

Elbow Block/ Trunk/ Lateral Pad Assembly


Used for the following models EBLT-XS/ EBRT-XS/ L-XS/ R-XS, EBLT-S/ EBRT-S/ L-S/ R-S, EBLT-M/ EBRT-M/ L-M/ R-M, EBLT-L/ EBRT-L/ L-L/ R-L




SHEP_250A




EBLT-XS/ EBRT-XS/ L-XS/ R-XS



EBLT-S/ EBRT-S/ L-S/ R-S



EBLT-M/ EBRT-M/ L-M/ R-M




EBLT-L/ EBRT-L/ L-L/ R-L

Elbow Block/ Trunk/ Lateral Pad Assembly

Used for the following models EBLT-XS/ EBRT-XS/ L-XS/ R-XS, EBLT-S/ EBRT-S/ L-S/ R-S, EBLT-M/ EBRT-M/ L-M/ R-M, EBLT-L/ EBRT-L/ L-L/ R-L

Item Number	Note	Part Number	Description	Quantity Per	
				Package	Assembly
1		60102832	Kit, Elbow Block, Trunk, Lateral Pad Assembly (4-1/4" W x 3-1/4" H)	1	1
1		60102828	Kit, Elbow Block, Trunk, Lateral Pad Assembly (5-1/4" W x 3-1/4" H)	1	1
1		60102831	Kit, Elbow Block, Trunk, Lateral Pad Assembly (5-3/4" W x 3-3/4" H)	1	1
1		60102827	Kit, Elbow Block, Trunk, Lateral Pad Assembly (6-1/2" W x 4-1/4" H)	1	1


Elbow Block/ Trunk/ Lateral Pad Assembly Fixed/Swing Away Hardware (FLATL/ FLATR/ LMSA/ RMSA)



Elbow Block/ Trunk/ Lateral Pad Assembly Fixed/Swing Away Hardware (FLATL/ FLATR/ LMSA/ RMSA)

Item Number	Note	Part Number	Description	Quantity Per	
				Package	Assembly
1		60105132	Kit, Swingaway/Adjustable Lateral Support Hardware - Left	1	1
1		60105131	Kit, Swingaway/Adjustable Lateral Support Hardware - Right	1	1
2		60105130	Sleeve, Lateral	1	1
3		60105129	Kit, Fixed Lateral Support Hardware - Right/Left	1	1
4		60105128	Kit, Lateral Link Extension, Telescoping (Pair) (up to 5-1/2" medial offset)	2	1

Matrix Swing Away Lateral Hardware, Standard, Elite/Elite TR Back (SAEHBL/ SAEHBR)






Matrix Swing Away Lateral Hardware, Standard, Elite/Elite TR Back - Right/Left (SAEHBL/ SAHEBR)

Matrx Swing Away Lateral Hardware, Standard, Elite/Elite TR Back (SAEHBL/ SAEHBR)

Item Number	Note	Part Number	Description	Quantity Per	
				Package	Assembly
1		60105127	Kit, Matrx Swing Away Lateral Hardware, Standard, Elite/Elite TR Back - Right/Left	1	1




Matrix Fixed/ Swing Away Lateral Hardware, Standard and Offset, Elite TR Back (PBLSL/R) (PBOL-L/R) (SAEHBL/R) (OS AHL/R)



Matrx Fixed/ Swing Away Lateral Hardware, Standard and Offset, Elite TR Back (PBLSL/R) (PBOL-L/R) (SAEHBL/R) (OSAHL/R)

Item Number	Note	Part Number	Description	Quantity Per	
				Package	Assembly
1	A	60105139	Kit, Matrx Swing Away Lateral Hardware, Offset, Elite/Elite TR Back - Left	1	1
1	B	60105140	Kit, Matrx Swing Away Lateral Hardware, Offset, Elite/Elite TR Back - Right	1	1
2	C	60105127	Kit, Matrx Swing Away Lateral Hardware, Standard, Elite/Elite TR Back - Right/Left	1	1
3	D	60104858	Kit, Matrx Fixed Lateral Hardware, Standard, Elite/Elite TR Back - Right/Left	1	1
4	E	60104859	Kit, Matrx Fixed Lateral Hardware, Offset, Elite/Elite TR Back - Left	1	1
4	F	60104860	Kit, Matrx Fixed Lateral Hardware, Offset, Elite/Elite TR Back - Right	1	1
NOTE: A - Used for model OSAHL B - Used for model OSAHR C - Used for model SAEHBL/SAEHBR D - Used for model PBLSL/PBLSR E - Used for model PBOL-L F - Used for model PBOL-R					


Lift-Off (LHSFL/ LHSFR) and Fixed Lateral Hip Support Hardware - Right/Left



Lift-Off (LHSFL/ LHSFR) and Fixed Lateral Hip Support Hardware - Right/Left

Item Number	Note	Part Number	Description	Quantity Per	
				Package	Assembly
1	A	60105197	Kit, Lift-Off Lateral Hip Support Hardware - Right/Left	1	1
2	A	60105196	Kit, Fixed Lateral Hip Support Hardware - Right/Left	1	1
NOTE: A - Used for either the Right or Left Side of the Seating System.					



Lateral Hip Support Pads (RHS-XS/ S/ M/ L) (LHS-XS/ S/ M/ L)



Lateral Hip Support Pads (RHS-XS/ S/ M/ L) (LHS-XS/ S/ M/ L)

Item Number	Note	Part Number	Description	Quantity Per	
				Package	Assembly
1	A	60102748	Kit, Lateral Hip Support Pad (5-1/2" W x 4" H)	1	1
2	A	60102749	Kit, Lateral Hip Support Pad (7" W x 6" H)	1	1
3	A	60102750	Kit, Lateral Hip Support Pad (9" W x 4" H)	1	1
4	A	60102751	Kit, Lateral Hip Support Pad (13" W x 4" H)	1	1
NOTE: A - Includes the Lateral Hip Pad and Hip Pad Cover. Used for Right or Left side application.					

MAXX Style Quick Release Removable Hardware For Hip Supports (MQRHL/ MQRHR)




SHEP_255A

MAXX Style Quick Release Removable Hardware For Hip Supports (MQRHL/ MQRHR)

Item Number	Note	Part Number	Description	Quantity Per	
				Package	Assembly
1		60102687	Kit, Hip Support Quick Release Hardware - Left	1	1
2		60102671	Kit, Hip Support Quick Release Hardware - Right	1	1

70 Degree HD Swing Away Front Rigging (M70HD)





70 Degree HD Swing Away Front Rigging (M70HD)


Item Number	Note	Part Number	Description	Quantity Per	
				Package	Assembly
1	A	60105198	Front Rigging Support Assembly, Down Tube HD, Cam-Lock, M70HD	1	1
2	B	60105200	Kit, Extension Tube Assembly w/ Hardware, M70HD (9" L)	1	1
NOTE: A - Item #1 can be used for both Right or Left side. B - Item #2 Extension Tube has a 13-17" Range					

Basic Fixed Center Mount Platform 70/ 90/ 97 Degrees (FCMT70/ 90/ 97)

Support Hanger
Bracket Assembly,
FCMT 70/90/97 Degree
shown for reference
only



Assembled View




SHEP_233A

Basic Fixed Center Mount Platform 70/ 90/ 97 Degrees (FCMT70/ 90/ 97)

Item Number	Note	Part Number	Description	Quantity Per	
				Package	Assembly
1		60105142	Basic Fixed Center Mount Footplate Assembly, 70/90/97 Degree (12" W x 11" D)	1	1
2		60105143	Flat Calf Pad Assembly, FCMT70, FCMT90, FCMT97	1	1

Center Mount Rigging (CTMT-FJ/CTMT-F)



Center Mount Rigging (CTMT-FJ/CTMT-F)

Item Number	Note	Part Number	Description	Quantity Per	
				Package	Assembly
1	A	60105180	Kit, Center Mount Front Rigging, Junior (CTMTJ/CTMT-FJ) ASBA/ADJASBA (9-1/2" - 11" Range)	1	1
2			Rack, Center Mount Rigging Angle Adjustable	1	1
3			Support, Footrest, Center Mount, Black - Short	1	1
4			Bearing, Flange (1/2 x 13/16 x 1/2")	1	2
5			Cleat, Button Head	1	2
6			Screw, Socket, Button Head w/ Patch (5/16-18 x 3/4")	1	6
7			Tube, Extension, Footrest, Short, Black - Right	1	1
7			Tube, Extension, Footrest, Short, Black - Left	1	1
8			Screw, Shoulder Hex Socket Head (5/16" x 2-1/2", 1/4-20)	1	1
9			Washer (1/4 x 1/2 x 1/16")	1	1
10			Locknut (1/4-20)	1	1
11			NLA - Washer, Wave (17/32 x 3/4 x 1/64")	1	2
12			Screw, Socket Button Head (1/4-20 x 1-1/2")	1	2
13			Washer, Nylon (5/16 x 11/16 x 1/16")	1	2
14			Washer (1/2 x 7/8 x 3/64")	1	2
15			Spacer, Nylon (1/2 x 3/4 x 1/4")	1	8
1	A	60105182	Kit, Center Mount Front Rigging (CTMT) (CTMT-F), ASBA/ADJASBA (11" - 15" Range)	1	1
2			Rack, Center Mount Rigging Angle Adjustable	1	1
3			Support, Footrest, Center Mount, Black - Regular	1	1
4			Bearing, Flange (1/2 x 13/16 x 1/2")	1	2
5			Cleat, Button Head	1	2
6			Screw, Socket, Button Head w/ Patch (5/16-18 x 3/4")	1	6
7			Tube, Extension, Footrest, Medium, Black - Right	1	1
7			Tube, Extension, Footrest, Medium, Black - Left	1	1
8			Screw, Shoulder Hex Socket Head (5/16" x 2-1/2", 1/4-20)	1	1
9			Washer (1/4 x 1/2 x 1/16")	1	1
10			Locknut (1/4-20)	1	1
11			NLA - Washer, Wave (17/32 x 3/4 x 1/64")	1	2
12			Screw, Socket Button Head (1/4-20 x 1-1/2")	1	2
13			Washer, Nylon (5/16 x 11/16 x 1/16")	1	2
14			Washer (1/2 x 7/8 x 3/64")	1	2
15			Spacer, Nylon (1/2 x 3/4 x 1/4")	1	8
16	B	60105181	Kit, Center Mount Front Rigging w/Calf Pad - Short (CTMTJ/CTMT-FJ) ASBA/ADJASBA (9-1/2" - 11" Range)	1	1
2			Rack, Center Mount Rigging Angle Adjustable	1	1
3			Support, Footrest, Center Mount, Black - Short	1	1
4			Bearing, Flange (1/2 x 13/16 x 1/2")	1	2
5			Cleat, Button Head	1	2

Center Mount Rigging (CTMT-FJ/CTMT-F)

Item Number	Note	Part Number	Description	Quantity Per	
				Package	Assembly
6			Screw, Socket, Button Head w/ Patch (5/16-18 x 3/4")	1	6
7			Tube, Extension, Footrest, Short, Black - Right	1	1
7			Tube, Extension, Footrest, Short, Black - Left	1	1
8			Screw, Shoulder Hex Socket Head (5/16" x 3", 1/4-20)	1	1
9			Washer (1/4 x 1/2 x 1/16")	1	1
10			Locknut (1/4-20)	1	1
11			NLA - Washer, Wave (17/32 x 3/4 x 1/64")	1	2
12			Screw, Socket Button Head (1/4-20 x 1-1/2")	1	2
13			Washer, Nylon (5/16 x 11/16 x 1/16")	1	2
14			Washer (1/2 x 7/8 x 3/64")	1	2
15			Spacer, Nylon (1/2 x 3/4 x 1/4")	1	8
17			Screw, Socket Button Head w/ Patch (1/4-20 x 3/4")	1	6
18			Washer, Nylon (5/16 x 7/8 x 1/8")	1	2
19			Bracket, Calf Pad, Black	1	1
20			Nut, T (5/16-18)	1	2
21			Adjuster, Calf Pad	1	2
22			Bracket, Mounting, Calf Pad	1	2
23			Locknut (1/4-20)	1	2
24			Screw, Socket, Button Head w/ Patch (5/16-18 x 1/2")	1	2
25			Pad, Calf (4-1/2" W x 7-1/4" L)	1	2
16	B	60105183	Kit, Center Mount Front Rigging w/Calf Pad - Medium (CTMT) (CTMT-F), ASBA/ADJASBA (11" - 15" Range)	1	1
2			Rack, Center Mount Rigging Angle Adjustable	1	1
3			Support, Footrest, Center Mount, Black - Regular	1	1
4			Bearing, Flange (1/2 x 13/16 x 1/2")	1	2
5			Cleat, Button Head	1	2
6			Screw, Socket, Button Head w/ Patch (5/16-18 x 3/4")	1	6
7			Tube, Extension, Footrest, Medium, Black - Right	1	1
7			Tube, Extension, Footrest, Medium, Black - Left	1	1
8			Screw, Shoulder Hex Socket Head (5/16" x 3", 1/4-20)	1	1
9			Washer (1/4 x 1/2 x 1/16")	1	1
10			Locknut (1/4-20)	1	1
11			NLA - Washer, Wave (17/32 x 3/4 x 1/64")	1	2
12			Screw, Socket Button Head (1/4-20 x 1-1/2")	1	2
13			Washer, Nylon (5/16 x 11/16 x 1/16")	1	2
14			Washer (1/2 x 7/8 x 3/64")	1	2
15			Spacer, Nylon (1/2 x 3/4 x 1/4")	1	8
17			Screw, Socket Button Head w/ Patch (1/4-20 x 3/4")	1	6
18			Washer, Nylon (5/16 x 7/8 x 1/8")	1	2

Center Mount Rigging (CTMT-FJ/CTMT-F)

Item Number	Note	Part Number	Description	Quantity Per	
				Package	Assembly
19			Bracket, Calf Pad, Black	1	1
20			Nut, T (5/16-18)	1	2
21			Adjuster, Calf Pad	1	2
22			Bracket, Mounting, Calf Pad	1	2
23			Locknut (1/4-20)	1	2
24			Screw, Socket, Button Head w/ Patch (5/16-18 x 1/2")	1	2
25			Pad, Calf (4-1/2" W x 7-1/4" L)	1	2
26		See Chart	#26 Extension Tube Kits	1	1
27		See Chart	#27 Footplate Kits	1	1
30	C	60102816	Kit, Calf Pad and Adjusting Hardware (Pair)	1	1
17			Screw, Socket Button Head w/ Patch (1/4-20 x 3/4")	1	2
20			Nut, T (5/16-18)	1	2
21			Adjuster, Calf Pad	1	2
22			Bracket, Mounting, Calf Pad	1	2
23			Locknut (1/4-20)	1	2
24			Screw, Socket, Button Head w/ Patch (5/16-18 x 1/2")	1	2
25			Pad, Calf (4-1/2" W x 7-1/4" L)	1	2
31	D	1163374	Kit, Footrest Pad - Small (Pair)	1	1
29			Pad, Footplate, Small - Right	1	1
29			Pad, Footplate, Small - Left	1	1
31	D	1163375	Kit, Footrest Pad - Medium (Pair)	1	1
29			Pad, Footplate, Medium - Right	1	1
29			Pad, Footplate, Medium - Left	1	1
31	D	1163376	Kit, Footrest Pad - Large (Pair)	1	1
29			Pad, Footplate, Large - Right	1	1
29			Pad, Footplate, Large - Left	1	1
	E	1165505	Kit, Hardware, Calf Pad Adjustment	1	1
17			Screw, Socket Button Head w/ Patch (1/4-20 x 3/4")	1	2
20			Nut, T (5/16-18)	1	2
21			Adjuster, Calf Pad	1	2
22			Bracket, Mounting, Calf Pad	1	2
23			Locknut (1/4-20)	1	2
24			Screw, Socket, Button Head w/ Patch (5/16-18 x 1/2")	1	2
	F	60102817	Kit, Bracket, Calf Pad Mounting w/ Hardware	1	1
8			Screw, Shoulder Hex Socket Head (5/16" x 3", 1/4-20)	1	1
9			Washer (1/4 x 1/2 x 1/16")	1	1
10			Locknut (1/4-20)	1	1
18			Washer, Nylon (5/16 x 7/8 x 1/8")	1	2
19			Bracket, Calf Pad, Black	1	1

Center Mount Rigging (CTMT-FJ/CTMT-F)

NOTE: A - Includes items 2-15.
 B - Includes items 2-15 and 17-25.
 C - Includes items 17, 20-25.
 D - Includes a pair of item 29.
 E - Not shown. Includes items 17 and 20-24.
 F - Not shown. Includes items 8-10, 18 and 19.

CTMT Rigging 12-1/4"-15" STF range does not work for 16" Deep on models ULMT, ULMTR and ULMR.

#26 Extension Tube Kits

Front Rigging Side	Short 9-1/4" - 12" Range	Medium 12-1/4" - 15" Range
Right	1170677 See Note(s) A	1170679 See Note(s) A
Left	1170678 See Note(s) A	1170680 See Note(s) A


NOTE: A - Includes items 5-7. The Extension Tube Kits are used for adjustable width footplate hardware.

#27 Footplate Kits

Footplate Side	Small 5-1/4" W x 7-1/2" L	Medium 5-1/2" W x 9-1/2" L	Large 5-7/8" W x 11-1/2" L
Right	1170683 See Note(s) A	1170685 See Note(s) A	1170687 See Note(s) A
Left	1170684 See Note(s) A	1170686 See Note(s) A	1170688 See Note(s) A

NOTE: A - Includes items 10-18 and the Adjustable Width Footplate Instruction Sheet (Not Shown). The footplate mounting hardware included is for adjustable width.


Center Mount Rigging Mounting Hardware (CTMT-FJ/CTMT-F)



Center Mount Rigging Mounting Hardware (CTMT-FJ/CTMT-F)

Item Number	Note	Part Number	Description	Quantity Per	
				Package	Assembly
1	A	60102591	CTMT Front Rigging Mounting Hardware, Ultra Low Maxx	1	1
NOTE: A - CTMT Rigging 12-1/4"-15" STF range does not work for 16" Deep on models ULMT, ULMTR and ULMR.					


Center Mount Rigging Support Hanger Brackets and Mounting Hardware (CTMT-FJ/F, FCMT70/ 90/ 97, LNX)



Center Mount Rigging Support Hanger Brackets and Mounting Hardware (CTMT-FJ/F, FCMT70/ 90/ 97, LNX)

Item Number	Note	Part Number	Description	Quantity Per	
				Package	Assembly
1	A	60102716	Kit, Center Mount Rigging Support Hanger Mounting Hardware (16-22" W)	1	1
2-3		60103320	Kit, Support Hanger Bracket Assembly w/ Mounting Hardware, Cable Release, CTMT-F/FJ	1	1
4		60105160	Support Hanger Bracket Assembly, FCMT 70/90/97 Degree	1	1
NOTE: A - Item #1 used for all Center Mount Riggings shown. LNX, CTMT-Fj/F and FCMT70/90/97.					

Center Mount Legrest Power (LNX)



Center Mount Legrest Power (LNX)

Item Number	Note	Part Number	Description	Quantity Per	
				Package	Assembly
1	D	60102726	Kit, LNX Power Legrest w/ Independant Footplates - Small	1	1
1	D	60103356	Kit, LNX Power Legrest w/ Independant Footplates - Medium	1	1
1	D	60102727	Kit, LNX Power Legrest w/ Independant Footplates - Large	1	1
2		60123531	Kit, LNX Actuator w/ Elev./Extending Assembly	1	1
3		60102709	Kit, Flat Calf Pad Assembly, LNX	1	1
4a	A,E	60102545	Kit, Actuator Cable, LNX, (MPELR - Left)	1	1
			Actuator Cable, LNX, (MPELR - Left)	1	1
			Sheet, Insert Parts Service	1	1
4b	F	60132729	Kit, Harness Wire Extension, 650mm 25-19/32"	1	1
			Harness Wire Extension, 650mm 25-19/32"	1	1
			Sheet, Insert Parts Service	1	1
5	B	60102713	Kit, Footboard Assembly, Center Mount, Fixed Heel Cup (12" W x 9" D)	1	1
6	B	60105162	Center Mount Flip-Up Foot Platform Assembly, LNX (12" W x 8-1/2" D)	1	1
6	B	60105161	Center Mount Flip-Up Foot Platform Assembly, LNX (12" W x 10-3/4" D)	1	1
6	B	60105163	Center Mount Flip-Up Foot Platform Assembly, LNX (17" W x 13" D)	1	1
7	C	60105164	Center Mount Foot Platform Mount Assembly, Flex Heel Cup	1	1
8		60105165	Independent Footplate Mounting Hardware, Center Mount (Pair)	1	1
9		1163374	Kit, Footrest Pad - Small (Pair)	1	1
			Pad, Footplate, Small - Right	1	1
			Pad, Footplate, Small - Left	1	1
9		1163375	Kit, Footrest Pad - Medium (Pair)	1	1
			Pad, Footplate, Medium - Right	1	1
			Pad, Footplate, Medium - Left	1	1
9		1163376	Kit, Footrest Pad - Large (Pair)	1	1
			Pad, Footplate, Large - Right	1	1
			Pad, Footplate, Large - Left	1	1
10		See Chart	#10 Footplate Kits	1	1
NOTE: A - The Left Side Actuator Cable is used for both LNX and MPELR Power Legrest. B - Item's 5 and 6 only include the lower Footboard Assembly portion of the Power LNX Legrest as shown. C - Item 6 and 7 assemble together. D - Kit includes the Actuator Motor. E - Used before 01/24/22 F - Used after 01/23/22 for both LNX and MPELR Power Legrest on both right and left side.					



#10 Footplate Kits

Center Mount Legrest Power (LNX)

#10 Footplate Kits

Footplate Side	Small 5-1/4" W x 7-1/2" L	Medium 5-1/2" W x 9-1/2" L	Large 5-7/8" W x 11-1/2" L
Right	1170683 See Note(s) A	1170685 See Note(s) A	1170687 See Note(s) A
Left	1170684 See Note(s) A	1170686 See Note(s) A	1170688 See Note(s) A
NOTE: A - Includes the Footplate, attaching hardware and the Adjustable Width Footplate Instruction Sheet (Not Shown).			

Independent Front Rigging Cross Brace Assembly (M70HD/ MELR/ MPELR), (16" W - 22" W)




SHEP_234B

Independent Front Rigging Cross Brace Assembly (M70HD/ MELR/ MPELR), (16" W - 22" W)

Item Number	Note	Part Number	Description	Quantity Per	
				Package	Assembly
1		60102718	Kit, Cross Brace Assembly, Independent Front Rigging	1	1


Independent Front Rigging Cam Lock Receiver, Straight (M70HD/ MELR/ MPELR)




Independent Front Rigging Cam Lock Receiver, Straight (M70HD/ MELR/ MPELR)

Item Number	Note	Part Number	Description	Quantity Per	
				Package	Assembly
1		60105141	Independent Front Rigging Cam Lock Receiver, Straight (Pair)	2	1

Maxx Style Pivot Plus Elevating Legrest



SHEP_238C




Maxx Style Pivot Plus Elevating Legrest

Item Number	Note	Part Number	Description	Quantity Per	
				Package	Assembly
1		60102846	Kit, Maxx Style Pivot Plus Elevating Legrest, Manual MELR - Left	1	1
1		60102845	Kit, Maxx Style Pivot Plus Elevating Legrest, Manual MELR - Right	1	1
2		60102720	Kit, Maxx Style Pivot Plus Elevating Legrest, Power MPELR - Left	1	1
2		60102719	Kit, Maxx Style Pivot Plus Elevating Legrest, Power MPELR - Right	1	1
3		60105167	Kit, Extension Tube Assembly w/ Hardware, MELR, MPELR, MSKTHL-S (11" L)	1	1
3		60105168	Kit, Extension Tube Assembly w/ Hardware, MELR, MPELR, MSKTHL-M (14" L)	1	1
3		60105169	Kit, Extension Tube Assembly w/ Hardware, MELR, MPELR, MSKTHL-L (17" L)	1	1
4	A,C	60102545	Kit, Actuator Cable, LNX, (MPELR - Left)	1	1
			Actuator Cable, LNX, (MPELR - Left)	1	1
			Sheet, Insert Parts Service	1	1
4	B,C	60105177	Kit, Actuator Cable (MPELR - Right)	1	1
			Actuator Cable (MPELR - Right)	1	1
			Sheet, Insert Parts Service	1	1
5	D	60132729	Kit, Harness Wire Extension, 650mm 25-19/32"	1	1
			Harness Wire Extension, 650mm 25-19/32"	1	1
			Sheet, Insert Parts Service	1	1
NOTE: A - The Left Side Actuator Cable is used for both LNX and MPELR Power Legrest. B - The Right Side Actuator Cable is used for MPELR Power Legrest only. C - Used before 01/24/22 D - Used after 01/23/22 for MPELR Power Legrest on both right and left side.					

Phono Jack Adapter Box Kit

For (MPELR) Maxx Style Pivot Plus Power Elevating Legrest




Phono Jack Adapter Box Kit

For (MPELR) Maxx Style Pivot Plus Power Elevating Legrest

Item Number	Note	Part Number	Description	Quantity Per	
				Package	Assembly
1	A,B	NPN	Kit, Phono Jack Adapter Box, MPELR	1	1
NOTE: A - Includes items 2-6 and accommodates both the right and left side of the Power Seating system. B - NOTE: "NPN" denotes that Phono Jack Adapter Box Kit is currently unavailable. For first time installation the wheelchair must be returned to Invacare to be installed.					


Adjustable Angle Footplate, Pivot and Multi-Axis, Independent Front Riggings



Adjustable Angle Footplate, Pivot and Multi-Axis, Independent Front Riggings

Item Number	Note	Part Number	Description	Quantity Per	
				Package	Assembly
1		60105170	Kit, Footplate Assembly, Heavy Duty w/ Heel Loop and Caster Wheel - Left (5-1/4" W x 8" L)	1	1
1		60105171	Kit, Footplate Assembly, Heavy Duty w/ Heel Loop and Caster Wheel - Right (5-1/4" W x 8" L)	1	1
1		60105172	Kit, Footplate Assembly, Heavy Duty w/ Heel Loop and Caster Wheel - Left (6-1/4" W x 8" L)	1	1
1		60105173	Kit, Footplate Assembly, Heavy Duty w/ Heel Loop and Caster Wheel - Right (6-1/4" W x 8" L)	1	1
2		60105174	Kit, Independent Front Rigging Pivot Mount Assembly - Right/Left (Pair)	2	1
3	A	60105175	Kit, Independent Front Rigging Multi-Axis Adjustable Angle Footplate Mount Assembly - Right/Left	1	1
NOTE: A - Item #3 can be used on the Right or Left side.					

Front Riggings - Accessories, Air Mesh Calf Panel CPMD/ CPLG









SHEP_237A


Front Riggings - Accessories, Air Mesh Calf Panel CPMD/ CPLG

Item Number	Note	Part Number	Description	Quantity Per	
				Package	Assembly
1	A	60105176	Kit, Calf Panel, Mesh, Medium CPMD (35-1/2" L, Overall)	1	1
1	A	60102594	Kit, Calf Panel, Mesh, Large CPLG (39-1/2" L, Overall)	1	1
NOTE: A - Used for M70HD Front Riggings.					

Remotes and Displays





SHEP_345A



Remotes and Displays

Item Number	Note	Part Number	Description	Quantity Per	
				Package	Assembly
1	A	60102503	Kit, LiNX Remote DLX-REM110	1	1
2	B	60102501	Kit, LiNX Remote DLX-REM210	1	1
2	C	60102502	Kit, LiNX Remote DLX-REM211	1	1
3	D	60124979	Kit, LiNX Remote DLX-REM215	1	1
3	E	60102504	Kit, LiNX Remote DLX-REM216	1	1
4	E	60102505	Kit, LiNX Remote DLX-REM400 Touch Screen	1	1
5	E	60108460	Kit, LiNX Remote DLX-REM400 Touch Screen w/ Toggles	1	1
6	E	60102506	Kit, LiNX Remote DLX-REM500 Display Only	1	1
7	F	GPL55470	Knob, Joystick (JS4REM050)	1	1
8		GPA55475	Knob, Joystick (DLX-REM400-A)	1	1
NOTE: A - Drive Only. B - Drive and Actuators 1 or 2. C - Drive, Actuators 2 or more D - Drive, 1 or 2 Actuators, Lights E - Drive, Actuators 2 or more, Lights F - Used for all Remotes except DLX-REM400					


Toggle Switch Kit (SW400TGL) For LiNX DLX-REM400 Remote



Toggle Switch Kit (SW400TGL) For LiNX DLX-REM400 Remote

Item Number	Note	Part Number	Description	Quantity Per	
				Package	Assembly
1	A	60108451	Kit, Toggle, DLX-REM400 Touch Screen (SW400TGL)	1	1
NOTE: A - Kit includes all items shown for Item #1.					

Joystick Tubes, Offset and Short




SHEP_208A

Joystick Tubes, Offset and Short

Item Number	Note	Part Number	Description	Quantity Per	
				Package	Assembly
1		60102482	Kit, Joystick Tube Assembly, Offset	1	1
2	A	60102483	Kit, Joystick Tube Assembly, Short	1	1
NOTE: A - Short Joystick Tube used for Quad Link Assembly and REM500 Display.					


Quad Link Assembly, Remote



Quad Link Assembly, Remote

Item Number	Note	Part Number	Description	Quantity Per	
				Package	Assembly
1	A	60104402	Kit, Quad Link Assembly and Mounting Hardware - Right	1	1
1	A	60104403	Kit, Quad Link Assembly and Mounting Hardware - Left	1	1
NOTE: A - Includes the Quad Link, Screws, Short Joystick Tube and Height Adjustable Clamp.					


LiNX DLX-CR400 Compact Remote (CREM) and Mounting Hardware



LiNX DLX-CR400 Compact Remote (CREM) and Mounting Hardware

Item Number	Note	Part Number	Description	Quantity Per	
				Package	Assembly
1		60102510	Kit, LiNX DLX-CR400 Compact Remote	1	1
2		60102583	Kit, LiNX CR400 (CREM)/ CR400-LF (CREM-LF) Compact Remote Mounting Hardware	1	1
3		60102522	Kit, Adapter, Nucleus Swing Away Midline Mount for ASL660	1	1
4		60102530	Kit, LiNX Bus Cable (PWH) 1500mm 59-1/16" (GLM-BUS150)	1	1
			LiNX Bus Cable (PWH) 1500mm 59-1/16" (GLM-BUS150)	1	1
			Sheet, Insert Parts Service	1	1
5		GPL55470	Knob, Joystick (JS4REM050)	1	1


LiNX DLX-CR400 Compact Remote Low Force (CREM-LF) and Mounting Hardware



LiNX DLX-CR400 Compact Remote Low Force (CREM-LF) and Mounting Hardware

Item Number	Note	Part Number	Description	Quantity Per	
				Package	Assembly
1		60102509	Kit, LiNX DLX-CR400LF Compact Remote Low Force	1	1
2		60102583	Kit, LiNX CR400 (CREM)/ CR400-LF (CREM-LF) Compact Remote Mounting Hardware	1	1
3		60102522	Kit, Adapter, Nucleus Swing Away Midline Mount for ASL660	1	1
4		60102530	Kit, LiNX Bus Cable (PWH) 1500mm 59-1/16" (GLM-BUS150)	1	1
			LiNX Bus Cable (PWH) 1500mm 59-1/16" (GLM-BUS150)	1	1
			Sheet, Insert Parts Service	1	1
5		GPL55470	Knob, Joystick (JS4REM050)	1	1

Swingaway Mounting Hardware, DLX-REM 500




Swingaway Mounting Hardware, DLX-REM 500

Item Number	Note	Part Number	Description	Quantity Per	
				Package	Assembly
1		60102485	Kit, Swingaway Mounting Hardware, DLX-REM 500 - Right	1	1
1		60102484	Kit, Swingaway Mounting Hardware, DLX-REM 500 - Left	1	1

LiNX Power Module 120AL and Mounting Hardware

NOTE: When installing a new Power Module, the Motors will need to be calibrated with Invacare Tech. Services. An LAK is required to complete the calibration.




LiNX Power Module 120AL and Mounting Hardware

NOTE: When installing a new Power Module, the Motors will need to be calibrated with Invacare Tech. Services. An LAK is required to complete the calibration.

Item Number	Note	Part Number	Description	Quantity Per	
				Package	Assembly
1	A	60124961	Kit, LiNX Power Module 120AL and Mounting Hardware	1	1
2			LiNX Power Module 120AL	1	1
3			Locknut (M5 x .8mm)	1	2
4		60124373	Kit, Battery Connector Harness, GP24	1	1
NOTE: A - Includes items 2 and 3					


G-Trac-LX - LiNX Gyro Module with Mounting Hardware



G-Trac-LX - LiNX Gyro Module with Mounting Hardware

Item Number	Note	Part Number	Description	Quantity Per	
				Package	Assembly
1	A	60124376	Kit, Module, Gyro w/ Hardware	1	1
2			Module, GYRO LiNX (310mm)	1	1
3			Nut (M6 x 1mm)	1	1
NOTE: A - Includes items 2 and 3					

Backrest Bracket Mounting Hardware, Recline - DLX-ACT400, DLX-IN500, DLX-OUT500, GLM-CONX4 (When ASL Options were ordered On-Chair)




SHEP_346B

Backrest Bracket Mounting Hardware, Recline - DLX-ACT400, DLX-IN500, DLX-OUT500, GLM-CONX4 (When ASL Options were ordered On-Chair)

Item Number	Note	Part Number	Description	Quantity Per	
				Package	Assembly
1	A	60102541	Kit, Backrest Bracket Mounting Hardware, Recline for DLX-IN500, CONX4, DLX ACT400	1	1
2		60102721	Kit, Knob, Hand Wheel (M6 x 15)	1	1
3	D,E	60102474	Kit, LiNX DLX-IN500 Input Module (INPUT)	1	1
4	E	60102475	Kit, LiNX 4-Channel Actuator Module (DLX-ACT400)	1	1
5	C	60108452	Kit, LiNX DLX-OUT500 Output Module (OUT500)	1	1
6		60102473	Kit, LiNX CONX4 4-way Expansion Block	1	1
7		60102672	Kit, Shroud, Backrest Recline	1	1
	B	60102523	Kit, Sip-N-Puff, 16" Fixed (PKG32666) Includes: (3) Air Tube Straws, (2) Saliva Traps, Tubing, Pivot Shaft/Flexible Metal Gooseneck Assembly, and Extruded Figure "8" Clamp to fit 7/8" Wheelchair Frame Tubing. Extruded Back Cane Adapter Clamp.	1	1
NOTE: A - Includes all hardware shown including item #2 and Tie Wraps (7 ea.). B - Not Shown C - Item #5 includes LiNX DLX-OUT500 Output Module, 2 ea. Washer (M4 4.3 x 9 x .8), 2 ea. Locknut (M4), 1 ea. LiNX Bus Cable (PWH) 300mm 11-13/16" (GLM-BUS030) D - Item #4 LiNX DLX-IN500 Input Module does not include the LiNX Bus Cable (PWH) 300mm 11-13/16" (GLM-BUS030). E - Item #3 and #4 switched mounting locations after 01/23/22					


Backrest Bracket Mounting Hardware, Recline for (DLX-ACT400, GLM-CONX4) (When Non-ASL Options were ordered On-Chair)



Backrest Bracket Mounting Hardware, Recline for (DLX-ACT400, GLM-CONX4) (When Non-ASL Options were ordered On-Chair)

Item Number	Note	Part Number	Description	Quantity Per	
				Package	Assembly
1	A	60125830	Kit, Backrest Bracket Mounting Hardware, Recline for DLX-ACT400/ GLM-CONX4	1	1
2		60102473	Kit, LiNX CONX4 4-way Expansion Block	1	1
3		60102475	Kit, LiNX 4-Channel Actuator Module (DLX-ACT400)	1	1
4		60125831	Kit, Recline Shroud	1	1
NOTE: A - Includes all hardware shown including, LiNX Electronics Buss Bracket, Recline Shroud/Buss Bracket, Strap Clamp and Tie Wraps (7 ea.)					


Backrest Bracket Mounting Hardware, Tilt or ASBA-2-DLX-IN200, DLX-ACT400, DLX-IN500, DLX-OUT500, GLM-CONX4



Backrest Bracket Mounting Hardware, Tilt or ASBA-2-DLX-IN200, DLX-ACT400, DLX-IN500, DLX-OUT500, GLM-CONX4

Item Number	Note	Part Number	Description	Quantity Per	
				Package	Assembly
1	C	60102496	Kit, Backrest Bracket Mounting Hardware, Tilt and Sling Back- DLX-ACT400, DLX-IN500, GLM-CONX4	1	1
2	A	60102475	Kit, LiNX 4-Channel Actuator Module (DLX-ACT400)	1	1
2	A,F	60124980	Kit, LiNX DLX-IN200 Input Module (INPUT)	1	1
3	E	60102474	Kit, LiNX DLX-IN500 Input Module (INPUT)	1	1
4	D	60108452	Kit, LiNX DLX-OUT500 Output Module (OUT500)	1	1
5		60102473	Kit, LiNX CONX4 4-way Expansion Block	1	1
6	G	60132728	Kit, Bracket Assembly, LiNX, Black	1	1
7		60102472	Kit, Shroud Cover, Module	1	1
	B	60102523	Kit, Sip-N-Puff, 16" Fixed (PKG32666) Includes: (3) Air Tube Straws, (2) Saliva Traps, Tubing, Pivot Shaft/Flexible Metal Gooseneck Assembly, and Extruded Figure "8" Clamp to fit 7/8" Wheelchair Frame Tubing. Extruded Back Cane Adapter Clamp.	1	1
NOTE: A - IN200 orACT400 are attached in the same location. B - Not Shown. C - Includes all hardware shown. D - Item #4 includes LiNX DLX-OUT500 Output Module, 2 ea. Washer (M4 4.3 x 9 x .8), 2 ea. Locknut (M4), 1 ea. LiNX Bus Cable (PWH) 300mm 11-13/16" (GLM-BUS030) E - Item #3 LiNX DLX-IN500 Input Module does not include the LiNX Bus Cable (PWH) 300mm 11-13/16" (GLM-BUS030). F - No SNP. Switch Inputs only. G - Item 6 used after 01/23/22					


Backrest Bracket Mounting Hardware, Sling Back Only for - DLX-IN200, DLX-IN500, DLX-OUT500, GLM-CONX4



Backrest Bracket Mounting Hardware, Sling Back Only for - DLX-IN200, DLX-IN500, DLX-OUT500, GLM-CONX4

Item Number	Note	Part Number	Description	Quantity Per	
				Package	Assembly
1		60102477	Kit, Backrest Bracket Mounting Hardware, Sling Back for - DLX-ACT400, DLX-IN500, GLM-CONX4	1	1
2		60124980	Kit, LiNX DLX-IN200 Input Module (INPUT)	1	1
3		60102474	Kit, LiNX DLX-IN500 Input Module (INPUT)	1	1
4	B,C	60108452	Kit, LiNX DLX-OUT500 Output Module (OUT500)	1	1
5		60102473	Kit, LiNX CONX4 4-way Expansion Block	1	1
6	D	60132728	Kit, Bracket Assembly, LiNX, Black	1	1
7		60102472	Kit, Shroud Cover, Module	1	1
	A	60102523	Kit, Sip-N-Puff, 16" Fixed (PKG32666) Includes: (3) Air Tube Straws, (2) Saliva Traps, Tubing, Pivot Shaft/Flexible Metal Gooseneck Assembly, and Extruded Figure "8" Clamp to fit 7/8" Wheelchair Frame Tubing. Extruded Back Cane Adapter Clamp.	1	1
NOTE: A - Not Shown B - Item #4 includes LiNX DLX-OUT500 Output Module, 2 ea. Washer (M4 4.3 x 9 x .8), 2 ea. Locknut (M4), 1 ea. LiNX Bus Cable (PWH) 300mm 11-13/16" (GLM-BUS030) C - Item #3 LiNX DLX-IN500 Input Module does not include the LiNX Bus Cable (PWH) 300mm 11-13/16" (GLM-BUS030). D - Item 6 used after 01/23/22					


LiNX CONX4 4-way Expansion Block, Shroud Cover and Hardware



LiNX CONX4 4-way Expansion Block, Shroud Cover and Hardware

Item Number	Note	Part Number	Description	Quantity Per	
				Package	Assembly
1		60102489	Kit, LiNX CONX4 and Mounting Hardware	1	1
2		60102473	Kit, LiNX CONX4 4-way Expansion Block	1	1
3	A,B	60102487	Kit, Seat Frame Shroud Assembly	1	1
NOTE: A - Use existing hardware when attaching the Seat Frame Shroud. B - A Foam Piece was added to the right and left side of the Seat Frame Shroud Assembly. The addition of the two foam pieces attached to the Shroud is backward compatible.					


C629 Junction Box, Shroud Cover and Hardware



C629 Junction Box, Shroud Cover and Hardware

Item Number	Note	Part Number	Description	Quantity Per	
				Package	Assembly
1	A	60124982	Kit, Junction Box, C629	1	1
2	B	60102487	Kit, Seat Frame Shroud Assembly	1	1
3	C	60124849	Kit, C629 Junction Box Extension Cable, Long 1010mm 39-49/64"	1	1
			C629 Junction Box Extension Cable, Long 1010mm 39-49/64"	1	1
			Sheet, Insert Parts Service	1	1
4	D	60124879	Kit, C629 Junction Box Extension Cable 400mm 15-3/4"	1	1
			C629 Junction Box Extension Cable 400mm 15-3/4"	1	1
			Sheet, Insert Parts Service	1	1
NOTE: A - C629 Junction Box includes the Hook and Loop Adhesive Strips. B - Use existing hardware when attaching the Seat Frame Shroud. C - Long Extension Cable used with Elevate Option. D - Extension Cable used with Non-Elevate Option.					

Cables - Bus (PWH), Adapter, Sensor, Actuator



Cables - Bus (PWH), Adapter, Sensor, Actuator

Item Number	Note	Part Number	Description	Quantity Per	
				Package	Assembly
1		60124878	Kit, Cable, Accessory w/ Plug	1	1
			Cable, Accessory w/ Plug 420mm 16-17/32"	1	1
			Sheet, Insert Parts Service	1	1
2		60124378	Harness, Wire, TYCO-ACT 270mm 10-5/8"	1	1
			Harness, Wire, TYCO-ACT 270mm 10-5/8"	1	1
			Sheet, Insert Parts Service	1	1
2		60102526	Kit, Adapter Cable Assembly, TYCO-ACT 500mm 19-11/16"	1	1
			Adapter Cable Assembly, TYCO-ACT 500mm 19-11/16"	1	1
			Sheet, Insert Parts Service	1	1
3a		60102546	Kit, Adapter Cable Assembly, DX2-ACT4 To Molex 2 Pin	1	1
			Adapter Cable Assembly, DX2-ACT4 To Molex 2 Pin	1	1
			Sheet, Insert Parts Service	1	1
3a	C	60102545	Kit, Actuator Cable, LNX, (MPELR - Left)	1	1
			Actuator Cable, LNX, (MPELR - Left)	1	1
			Sheet, Insert Parts Service	1	1
3a	C	60105177	Kit, Actuator Cable (MPELR - Right)	1	1
			Actuator Cable (MPELR - Right)	1	1
			Sheet, Insert Parts Service	1	1
3b	D	60132729	Kit, Harness Wire Extension, 650mm 25-19/32"	1	1
			Harness Wire Extension, 650mm 25-19/32"	1	1
			Sheet, Insert Parts Service	1	1
4		60102528	Kit, LiNX Bus Cable (PWH) 300mm 11-13/16" (GLM-BUS030)	1	1
			LiNX Bus Cable (PWH) 300mm 11-13/16" (GLM-BUS030)	1	1
			Sheet, Insert Parts Service	1	1
4		60102529	Kit, LiNX Bus Cable (PWH) 700mm 27-9/16" (GLM-BUS070)	1	1
			LiNX Bus Cable (PWH) 700mm 27-9/16" (GLM-BUS070)	1	1
			Sheet, Insert Parts Service	1	1
4		60102536	Kit, LiNX Bus Cable (PWH) 1200mm 47-1/4" (GLM-BUS120)	1	1
			LiNX Bus Cable (PWH) 1200mm 47-1/4" (GLM-BUS120)	1	1
			Sheet, Insert Parts Service	1	1
4		60102530	Kit, LiNX Bus Cable (PWH) 1500mm 59-1/16" (GLM-BUS150)	1	1
			LiNX Bus Cable (PWH) 1500mm 59-1/16" (GLM-BUS150)	1	1
			Sheet, Insert Parts Service	1	1
4		60102531	Kit, LiNX Bus Cable (PWH) 1700mm 66-59/64" (GLM-BUS170)	1	1

Cables - Bus (PWH), Adapter, Sensor, Actuator


Item Number	Note	Part Number	Description	Quantity Per	
				Package	Assembly
4		60124987	LiNX Bus Cable (PWH) 1700mm 66-59/64" (GLM-BUS170)	1	1
			Sheet, Insert Parts Service	1	1
			Kit, LiNX Bus Cable (PWH) 1000mm 39-3/8" (GLM-BUS100)	1	1
			LiNX Bus Cable (PWH) 1000mm 39-3/8" (GLM-BUS100)	1	1
4		60126411	Sheet, Insert Parts Service	1	1
			Kit, LiNX Bus Cable (PWH) 500mm 19-11/16" (GLM-BUS 050)	1	1
			LiNX Bus Cable (PWH) 500mm 19-11/16" (GLM-BUS 050)	1	1
			Sheet, Insert Parts Service	1	1
5a		60124379	Kit, Harness, Wire Extension 1040mm 40-15/16"	1	1
			Harness, Wire Extension 1040mm 40-15/16"	1	1
5a	C	60124384	Sheet, Insert Parts Service	1	1
			Harness, Wire Extension 650mm 25-9/32"	1	1
5b	D	60132730	Harness, Wire Extension 650mm 25-9/32"	1	1
			Sheet, Insert Parts Service	1	1
6	A	60124849	Kit, Harness Wire Extension, 1473mm 57-63-64"	1	1
			Harness Wire Extension, 1473mm 57-63-64"	1	1
			Sheet, Insert Parts Service	1	1
			Kit, C629 Junction Box Extension Cable, Long 1010mm 39-49/64"	1	1
7	B	60124879	C629 Junction Box Extension Cable, Long 1010mm 39-49/64"	1	1
			Sheet, Insert Parts Service	1	1
			Kit, C629 Junction Box Extension Cable 400mm 15-3/4"	1	1
			C629 Junction Box Extension Cable 400mm 15-3/4"	1	1
8		60126412	Sheet, Insert Parts Service	1	1
			Kit, LiNX Bus Extension Cable 640mm 25-13/64" (GLM-BUSEXT 064)	1	1
			LiNX Bus Extension Cable 640mm 25-13/64" (GLM-BUSEXT 064)	1	1
			Sheet, Insert Parts Service	1	1

NOTE: A - Long Extension Cable used with Elevate Option.
 B - Extension Cable used with Non-Elevate Option.
 C - Used before 01/24/22
 D - Used after 01/23/22

NOTE: Refer to the Wiring Diagrams for detail application and configurations of Cables and Harnesses located at www.invacare.com/homecare under Product Literature Search/ Service Manual. 60108795 "LiNX Control System, AVIVA FX", 60107431 "AVIVA FX Base", 60108797 "Ultra Low Maxx by Motion Concepts Modular Positioning System".

Wire Guide Bracket and Mounting Hardware

ULMT/ULMTR with AELM




Wire Guide Bracket and Mounting Hardware

ULMT/ULMTR with AELM

Item Number	Note	Part Number	Description	Quantity Per	
				Package	Assembly
1	A	60125834	Kit, Wire Guide Bracket and Mounting Hardware	1	1
2			Bracket, Wire Guide	1	1
3			Screw, TORX CS-ME A2K DIN 7500 (M3 x 10mm)	1	2
4			Tie Wrap (11-1/2" L)	1	2
			Sheet, Insert Parts Service	1	1
NOTE: A - Includes items 2-4					

Tilt Angle Sensor Assembly w/ Hardware




Tilt Angle Sensor Assembly w/ Hardware

Item Number	Note	Part Number	Description	Quantity Per	
				Package	Assembly
1	A	60125835	Kit, Tilt Angle Sensor w/ Hardware	1	1
2			Tilt Angle Sensor w/ Connector Assembly	1	1
3			Screw, Socket Head Cap (M4 x 0.7 x 10mm)	1	1
			Sheet, Insert Parts Service	1	1
4	B	60124877	Kit, Tilt Angle Sensor Bracket w/ Hardware, ULMTR	1	1
5			Screw, Socket Head Cap w/ Patch (M8 x 1.25 x 16mm)	1	1
6			Nut, Tee (M8 x 1.25mm) 2 Hole	1	1
7			Plate, Tilt Sensor Mounting	1	1
	C	60124958	Sheet, Insert Parts Service	1	1
8			Kit, Tilt Angle Sensor Bracket, ULMT	1	1
			Bracket, Sensor Angle Non-Recline	1	1
			Sheet, Insert Parts Service	1	1
NOTE: A - Includes items 2 and 3 B - Used for ULMTR. Includes items 5-7 C - Used for Non-Recline Rehab Seating System.					

ASL 611 Mount For The Egg Switch

ASL 611 Mount is used to attach an Egg Switch when the Hook & Loop Adhesive method is not preferred





SHEP_217A

ASL 611 Mount For The Egg Switch

ASL 611 Mount is used to attach an Egg Switch when the Hook & Loop Adhesive method is not preferred

Item Number	Note	Part Number	Description	Quantity Per	
				Package	Assembly
1	A	60102517	Kit, Switch, Black Egg Reset w/ Hook/Loop Adhesive Back Fastener & Attaching Screws (EGSBLK)	1	1
2	A	60102525	Kit, Mount, ASL 611	1	1
NOTE: A - Individual components are not available					

Dual Toggle Switch (2SW) & Monoport "Y" Phono Cable (MPY)




Egg Switch (EGSBLK)

SHEP_343A

Dual Toggle Switch (2SW) & Monoport "Y" Phono Cable (MPY)

Item Number	Note	Part Number	Description	Quantity Per	
				Package	Assembly
1	A	1091496	Kit, Joystick Tube Toggle Switch Mounting Hardware	1	1
2			NLA - Half Clamp, Threaded Holes	1	1
3			Half Clamp, Non-Threaded Holes	1	1
4			Screw, Hex Head w/ Patch (1/4-20 x 1/2")	1	2
5			NLA - Bracket, Toggle Switch Mounting, Black	1	1
6			NLA - Screw, Hex Head (1/4-20 x 1-1/8")	1	2
7		60102516	Kit, Switch, Dual Function Toggle (2SW)	1	1
8		1140027	Cable, Monoport "Y" Phono (MPY)	1	1
9		60102517	Kit, Switch, Black Egg Reset w/ Hook/Loop Adhesive Back Fastener & Attaching Screws (EGSBLK)	1	1
NOTE: A - Includes items 2-6					

Dual (2PB)/ Quad (QPBM) Push Button Switch w/ Multi-Mount Hardware




SHEP_341A

Dual (2PB)/ Quad (QPBM) Push Button Switch w/ Multi-Mount Hardware

Item Number	Note	Part Number	Description	Quantity Per	
				Package	Assembly
1	A	60124955	Kit, Dual Push Button Switch (2PB) and Multi-Mount Hardware	1	1
2			Kit, Switch, Dual Push Button (2PB)	1	1
			Hardware, Multi-Mount, Dual Push Button Switch (2PB)	1	1
			Tie Wrap (11-1/2" L)	1	10
			Manual, Push Button Switches Service	1	1
1	B	60124953	Kit, Quad Push Button Switch (QPB) and Multi-Mount Hardware	1	1
3			Kit, Switch, Quad Push Button (QPB)	1	1
			Hardware, Multi-Mount, Quad Push Button Switch (QPB)	1	1
			Tie Wrap (11-1/2" L)	1	10
			Manual, Push Button Switches Service	1	1
NOTE: A - Kit includes items 2, ALL Multi-Mount Bracketry and Hardware shown, Tie Wraps (not shown) and Service Manual (not shown). B - Kit includes items 3, ALL Multi-Mount Bracketry and Hardware shown, Tie Wraps (not shown) and Service Manual (not shown). NOTE: The Dual or Quad Push Button Switch can assemble to any Multi-Mount Bracketry combination shown.					


4-Way Toggle Switch (FWTM) and Mounting Hardware




4-Way Toggle Switch (FWTM) and Mounting Hardware

Item Number	Note	Part Number	Description	Quantity Per	
				Package	Assembly
1		60124978	Kit, 4-Way Toggle W/ 10 Pin Molex	1	1
2		60102550	Kit, Toggle Switch Multi-Mount Hardware Assembly	1	1

LiNX DLX-ACU200 Compact Attendant Remote



**LiNX DLX-ACU200
ULMR/ULMTR
Spreader Bar Mounting Hardware**




**LiNX DLX-ACU200
ULMT/ASBA-2
Back Cane Mounting Hardware**




LiNX DLX-ACU200 Compact Attendant Remote

Item Number	Note	Part Number	Description	Quantity Per	
				Package	Assembly
1	A	60102512	Kit, LiNX DLX-ACU200 Compact Attendant Remote	1	1
2		60102521	Kit, LiNX DLX-ACU200 Compact Attendant Remote Mounting Hardware, Recline	1	1
3		60105178	LiNX DLX-ACU200 Compact Attendant Remote Mounting Hardware, Non-Recline	1	1
NOTE: A - A bus cable is needed when connecting the Compact Attendant Remote to the wheelchair.					

LiNX DLX-SW100 10-Way Switch and Mounting Hardware



Remote Bracket Assembly
shown for reference only




SHEP_221A

LiNX DLX-SW100 10-Way Switch and Mounting Hardware

Item Number	Note	Part Number	Description	Quantity Per	
				Package	Assembly
1	A	60102471	Kit, LiNX DLX-SW100 10-Way Switch	1	1
2		60102542	Kit, LiNX DLX-SW100 Mounting Hardware, Front, REM400	1	1
3		60102538	Kit, LiNX DLX-SW100 Mounting Hardware, Front, REM211/REM216	1	1
4		60102540	Kit, LiNX DLX-SW100 Mounting Hardware, Rear Shroud	1	1
5	A	60102539	Kit, Accessory Cable,LiNX, w/ 10-Way Switch and Lift AELM	1	1
			Accessory Cable,LiNX, w/ 10-Way Switch and Lift AELM	1	1
			Sheet, Insert Parts Service	1	1
NOTE: A - Item #5 Accessory Cable is always used with a LiNX DLX-SW100 10-Way Switch.					



Accessories - USB Charger and LiNX Access Key (LAK)



Accessories - USB Charger and LiNX Access Key (LAK)

Item Number	Note	Part Number	Description	Quantity Per	
				Package	Assembly
1		60102519	Kit, LiNX USB Charger and Mounting Hardware	1	1
2		60102470	Kit, Key, LiNX Access	1	1

Accessories - Remote (Joystick) Tops





SHEP_223C

Accessories - Remote (Joystick) Tops

Item Number	Note	Part Number	Description	Quantity Per	
				Package	Assembly
1		60102547	Kit, Handle, Joystick, 3" U Shaped	1	1
1		60102549	Kit, Handle, Joystick, 4" U Shaped	1	1
2		60102588	Kit, Handle, Joystick, Dome Rubber	1	1
3		60106526	Kit, Handle, Joystick, Grooved Mushroom	1	1
4		60102524	Kit, Chin Cup Assembly - 1826	1	1
	A	60102586	Kit, Extension, Joystick, "T" Style (1560) (4" W x 5-1/2" L)	1	1
	A	60102585	Kit, Extension, Joystick, Straight Style (1561) (1/2" W x 5" L)	1	1

NOTE: A - Not Shown

Positioning Straps




SHEP_228A

Positioning Straps

Item Number	Note	Part Number	Description	Quantity Per	
				Package	Assembly
1		60102674	Kit, 2 Point Padded Positioning Strap Assembly (65" L)	1	1
2		60102673	Kit, 4 Point Padded Positioning Strap Assembly (65" L)	1	1
3		60102722	Kit, Positioning Strap Assembly w/ Hardware, Push Button Auto Style (48" L)	1	1
3		60102723	Kit, Positioning Strap Assembly w/ Hardware, Push Button Auto Style (60" L)	1	1
3		60102724	Kit, Positioning Strap Assembly w/ Hardware, Push Button Auto Style (69" L)	1	1

Accessories - TRRO Pelvic Belt and Hardware




Accessories - TRRO Pelvic Belt and Hardware

Item Number	Note	Part Number	Description	Quantity Per	
				Package	Assembly
1	A	60124876	Kit, TRRO Pelvic Belt Hardware (Pair)	1	1
2	B,C	1168513	Kit, Belt, Pelvic Restraint w/ Auto Style Push Button Buckle (Bushing Hole to Bushing Hole Length 37")	1	1
2	B,D	1168514	Kit, Belt, Pelvic Restraint w/ Auto Style Push Button Buckle (Bushing Hole to Bushing Hole Length 55")	1	1
NOTE: A - Sold complete. Individual components are non-saleable. Accommodates both right and left side. B - TRRO Pelvic Belt Hardware is not included with item #2. C - For 16" wide model D - For 17" and wider models					

Accessories - Vent Tray Assembly w/ Shroud

Vent Tray Assembly used on Non-Elevate Recline Only and Tilt/Recline Seating Systems.



Accessories - Vent Tray Assembly w/ Shroud

Vent Tray Assembly used on Non-Elevate Recline Only and Tilt/Recline Seating Systems.

Item Number	Note	Part Number	Description	Quantity Per	
				Package	Assembly
1		60105138	Kit, Vent Tray Assembly w/ Shroud	1	1

Invacare Corporation

www.invacare.com



One Invacare Way
Elyria, Oh USA
44036-2125
800-333-6900

Technical Support
800-832-4707
Customer Service
800-333-6900

Literature Orders
800-333-6900

© 2022 Invacare Corporation. All rights reserved. Republication, duplication or modification in whole or in part is prohibited without prior written permission from Invacare. Trademarks are identified by ™ and ®. All trademarks are owned by or licensed to Invacare Corporation or its subsidiaries unless otherwise noted. (Revised 10/11/22)

Form No. 19_118