


Rear Wheels & Anti-Tippers



Item Number	Note	Part Number	Description	Quantity Per	
				Package	Assembly
1	A	1144834-NLA	NLA - Wheel Assembly - Rear (L-4) (Single)	1	2
2		1144823	Wheel - Rear (L-3)	1	2
3		1144824	Kit, Anti-Tip Wheels Includes: (2) ea. Nylon Nut, (4) ea. Bearings, (2) ea. Bolts, and (2) ea. Anti-Tip Wheels	1	1
NOTE: A - No Longer Available. Order item #2 as a replacement					