



Steering Assembly

Panther LX4



Panther MX4



Steering Assembly

Item Number	Note	Part Number	Description	Quantity Per	
				Package	Assembly
1		1088726-NLA	NLA - Steering Assembly	1	1
2		1088727-NLA	NLA - Steering Assembly	1	1
3		1088787-NLA	NLA - Tie Rods, Panther LX4 (Pair - 1 Long, 1 Short)	2	1
4		1088766-NLA	NLA - Tie Rods, Panther MX4 (Pair - 1 Long, 1 Short)	2	1
	A	1088798-NLA	NLA - Steering Hardware - Panther LX4 Includes: 2 each Bolt (M6 x 1.0 x 35) 2 each Nut, Self-Lock (6 x 10 x 8) 2 each Nut, Self-Lock (8 x 13) 2 each Bolt (M8 x 1.25 x 20) 1 each U-Plate	1	1
	A	1091141-NLA	NLA - Steering Hardware - Panther MX4 Includes: 2 each Bolt (M6 x 1.0 x 35) 2 each Nut, Self-Lock (6 x 10 x 5) 1 each Screw, Hex (M8 x 1.25 x 30) 1 each Washer (8.1 x 13 x 3) 4 each Nut, Self-Lock (8 x 1.25 x 13) 2 each Screw, Hex (M8 x 1.25 x 35) 2 each Spacer (12 x 8) 1 each Screw, Hex (M8 x 1.25 x 60) 1 each Spacer (12 x 35)	1	1
NOTE: A - Not Shown					