


Motor & Transaxle



Item Number	Note	Part Number	Description	Quantity Per	
				Package	Assembly
1		1144788-NLA	NLA - Kit, Motor & Brake Includes: (1) ea. Motor, (1) ea. Guide, (4) ea. Bolts, and (2) ea. Keys	1	1
2		1144789	Kit, Brake Includes: (1) ea. Brake and (2) ea. Screws	1	1
3	A,A1	1144787	Kit, Transaxle Includes: (1) ea. Differential, 2 ea. Keys, 2 ea. E Rings, 2 ea. Plain Washers, 2 ea. Bolts (M8 x 20L), & 2 ea. Flat Washer	1	1
<p>NOTE: A - For all manufacture dates, an additional bracket is required when replacing the Transaxle. This bracket can only be ordered thru TAG. NOTE: TAG no longer sells the Bracket.</p> <p>A1 - The Bracket from TAG must be on the Scooter for the Transaxle to fit properly. Order a Transaxle replacement for only Scooters equipped with a previously purchased Bracket from TAG.</p>					