


Trans-Axle Assembly

Lynx LX3/LX3Plus


Panther LX4/MX4

Lynx LX3/LX3Plus

Panther LX4



Panther MX4



Trans-Axle Assembly

Lynx LX3/LX3Plus

Panther LX4/MX4

Item Number	Note	Part Number	Description	Quantity Per	
				Package	Assembly
1		1088731-NLA	NLA - Trans-Axle w/ Motor and Brake - Lynx LX3/ LX3Plus and Panther LX4	1	1
1		1088732-NLA	NLA - Trans-Axle w/ Motor and Brake - Panther MX4	1	1
2		1088734-NLA	NLA - Motor and Brake - Lynx LX3/ LX3 Plus and Panther LX4	1	1
2		1088736-NLA	NLA - Motor and Brake - Panther MX4	1	1
3		1088739-NLA	NLA - Lever, Brake Release - Lynx LX3/ LX3 Plus and Panther LX4	1	1
3		1088740-NLA	NLA - Lever, Brake Release - Panther MX4	1	1
4		1091136-NLA	NLA - Hardware, Trans-Axle - Lynx LX3/ LX3 Plus and Panther LX4 Includes: 2 each Rubber 2 each Bracket 4 each Bolt (M8 x 1-1/4 x 45) 4 each Washer (8 x 22 x 2) 4 each Washer, Spring (8 x 13) 4 each Nut, Self-Lock (8 x 13)	1	1
4		1091142-NLA	NLA - Hardware, Trans-Axle - Panther MX4 Includes: 1 each Rubber 1 each Bracket 1 each Bracket 1 each Rubber 2 each Bolt (M8 x 1.25 x 45) 2 each Screw (M8 x 1.25 x 50) 4 each Washer (8 x 22 x 2) 2 each Washer, Spring (8 x 13) 4 each Nut, Self-Lock (8 x 1.25 x 13)	1	1
5		1099710-NLA	NLA - Trans-Axle - Lynx LX3/ LX3 Plus and Panther LX4	1	1
5		1099711-NLA	NLA - Trans-Axle - Panther MX4	1	1
	A,B	1099736-NLA	NLA - Brush, Motor (Pair)	2	1
NOTE: A - Not Shown B - Brushes are for use on Motors manufactured before December 2002; any Motor Manufactured after December 2002, a new Motor will need to be ordered.					