





Seat Assembly

Lynx SX3P/SX3/LX3/LX3 Plus

Panther LX4/MX4

1  Arms not shown	2  Arms not shown	3  Arms not shown
Contour Cloth Seat (Includes: Arms and Seat Plate)	Deluxe Vinyl Seat (Includes: Arms and Seat Plate)	Deluxe Vinyl Extra Wide Seat (Includes: Arms and Seat Plate)



Seat Assembly

Lynx SX3P/SX3/LX3/LX3 Plus

Panther LX4/MX4

Item Number	Note	Part Number	Description	Quantity Per	
				Package	Assembly
1	A,C	1088753-NLA	NLA - Seat Assembly, Contour Cloth - Lynx SX3P/ SX3/ LX3	1	1
2	A	1088754-NLA	NLA - Seat Assembly, Deluxe Vinyl - Lynx LX3/ LX3 Plus	1	1
2	A	1091861-NLA	NLA - Seat Assembly, Deluxe Vinyl with Sliders - Panther LX4	1	1
3	A	1088755-NLA	NLA - Seat Assembly, Deluxe Extra Wide Vinyl w/ Sliders - Panther MX4	1	1
4		1088759-NLA	NLA - Headrest, Vinyl (Black and White) - Lynx LX3/ LX3 Plus and Panther LX4/ MX4	1	1
4		1106211-NLA	NLA - Headrest, Vinyl (Gray) - Lynx LX3/ LX3 Plus and Panther LX4/ MX4	1	1
5		1088757-NLA	NLA - Upholstery w/ Foam, Contour Seat and Back, Lynx SX3P/ SX3/ LX3 and Zoom-3 - Cloth	1	1
5		1088758-NLA	NLA - Upholstery with Foam, Contour Seat, and Back	1	1
6		1088756-NLA	NLA - Handle, Seat Swivel with Pin	1	1
7		1102926-NLA	NLA - Post, Adjustable Height Seat	1	1
	B	P715	Strap, Seat Positioning - All Models	1	1
NOTE: A - Seat Assembly includes Arms and Seat Plate B - Not Shown C - No Longer Available.					