



Headlight Assembly

Lynx LX3/LX3^{Plus}
Panther LX4



Panther MX4



Headlight Assembly

Item Number	Note	Part Number	Description	Quantity Per	
				Package	Assembly
1		1088771-NLA	NLA - Head Light Assembly - Lynx LX3/ LX3 Plus/ Panther LX4	1	1
1		1088772-NLA	NLA - Head Light Assembly - Panther MX4	1	1
2		1088773-NLA	NLA - Bulbs, Head Light - Lynx LX3/ LX3 Plus and Panther LX4	2	1
2		1091868-NLA	NLA - Bulbs, Head Light - Panther MX4	2	1
	A	1088735-NLA	NLA - Head Light Hardware - Lynx LX3/ LX3 Plus Includes: 1 each Plate 2 each Bolt, (M6 x 1.0 x 35) 2 each Nut, Self-Lock (6 x 10 x 8)	1	1
	A	1091144-NLA	NLA - Head Light Hardware - Panther MX4 Incudes: 1 each Screw (M4 x 0.7 x 15) 1 each Washer, Spring (4 x 7 x 1.25) 3 each Washer (4 x 12 x 0.8) 1 each Screw (M4 x 0.7 x 80) 2 each Washer (6 x 18 x 1.6) 2 each Washer, Spring (6 x 10) 2 each Screw, Nylock (M6 x 1.0 x 16)	1	1
NOTE: A - Not Shown					