


# Back Cane Attaching Hardware

For ADJASBA Seats (16"-24" Wide) with TRRO Option



## Back Cane Attaching Hardware

### For ADJASBA Seats (16"-24" Wide) with TRRO Option

| Item Number | Note | Part Number | Description   | Quantity Per |          |
|-------------|------|-------------|---|--------------|----------|
|             |      |             |   | Package      | Assembly |
| 1           | A    | 1137234-NLA | NLA - Kit, Back Cane Mounting Hardware For TRRO Option (Pair) | 2            | 1        |
| 2           |      |             | Washer, Flat (21/64 x 53/64 x 7/64")                          | 1            | 4        |
| 3           |      |             | Bolt, Round Head Square Neck (5/16-18 x 1-1/4")               | 1            | 2        |
| 4           |      |             | Bolt, Round Head Square Neck (5/16-18 x 1")                   | 1            | 2        |
| 5           |      |             | Locknut (5/16-18)   | 1            | 4        |
| 6           |      |             | Screw w/ Lug (5/16-24 x 1-3/8")                               | 1            | 2        |
| 7           |      |             | Locknut (5/16-24)   | 1            | 2        |

NOTE: A - Includes items 2-7, Right and Left sides