



Arm Mounting Swingaway Hardware

ADJASBA/ADJRECL Seats



ADJASBA Seat Frame 12"-16"
wide shown for reference only

Hex Cap Screw
shown for reference
only. See ADJASBA
Seat Frame, Seat Pan
and Hardware page



ADJASBA Seat Frame
16"-24" wide shown
for reference only

Hex Cap Screw
shown for reference only.
See ADJASBA Seat Frame,
Seat Pan and Hardware page

ADJRECL Manual Recline
Rear Arm Socket shown

Arm Mounting Swingaway Hardware ADJASBA/ADJRECL Seats

Item Number	Note	Part Number	Description	Quantity Per	
				Package	Assembly
1	A	1068488	Kit, Swingaway Armrest Mounting Hardware	1	1
2		1019590	Screw, Socket Button Head (1/4-20 x 1-5/8")	1	1
3		1068025	Clevis, Pivot	1	1
4		1025195	Locknut (1/4-20)	1	1
5		1068039	Pin, Pull	1	1

NOTE: A - Includes items 2-5