


Packaging and Literature

For 3G Storm Series



Packaging and Literature For 3G Storm Series

Item Number	Note	Part Number	Description	Quantity Per	
				Package	Assembly
1	A	1030364	Carton, Cover	1	1
1	B	1031682	Carton, Cover	1	1
1	C	1036028	Carton, Cover	1	1
2		1043198-NLA	NLA - Tube, Paper	1	1
3		1015609-NLA	NLA - Bag, Poly (28" W x 44" D x 80" H)	1	1
4	D	1030362	Carton, Tray	1	2
4	E	1036027	Carton, Tray	1	2
5	F	See Note	See Footnote	1	1
6		1021149-NLA	NLA - Bag, Poly Zip-Lip (13" L x 9" W)	1	1
7		1143155-NLA	NLA - Owner's Manual, Formula CG Seating	1	1
8		00077X031-NLA	NLA - Warranty Card	1	1
9		50000X121	Strap, Battery Lift	1	1
10	H	1012001-NLA	NLA - Carton, Joystick	1	1
10		40131X011-NLA	NLA - Carton, Front Riggings	1	1
11		40131X161	Insert, Plain Pad	1	1
12	A	1115174	Carton, Corner Reinforcement	1	2
12	B	1129589-NLA	Carton, Corner Reinforcement	1	2
12	C	1129590-NLA	Carton, Corner Reinforcement	1	2
	G	1164348	Mailer, Return, Mercury Free	1	1
NOTE: A - 12"- 20" wide systems with base B - 16"- 20" wide systems medium/tall seat-to-floor heights with base C - 21"- 24" wide systems with base D - 12"- 20" wide systems with base E - 21"- 24" wide systems with base F - Skid not supplied G - For Massachusetts providers only. H - Use part #40131X011 for a Joystick also.					
NOTE: When adjusting or replacing parts, please refer to your Owner's Manual for assistance.					