


Aerosol Compressor - Envoy

Model - IRC1192, IRC1192D, IRC1194A, IRC1194B


Plug Styles are Shown for Reference Only



IRC1194A equipped with Australian Plug



IRC1194B equipped with British Plug



Aerosol Compressor - Envoy

Model - IRC1192, IRC1192D, IRC1194A, IRC1194B

Item Number	Note	Part Number	Description	Quantity Per	
				Package	Assembly
1	J	2002837-NLA	NLA - Envoy Shroud Assembly	1	1
2		2002376-NLA	NLA - Door, Storage Compartment	1	1
3		2002594-NLA	NLA - Cover, Filter	1	1
4		650391-NLA	NLA - Filter, Foam Intake	6	1
		1156679-NLA	NLA - Package, Intake Filter	12	1
5		2002380-NLA	NLA - Switch, Power, IRC1192, IRC1192D	1	1
5		2003025-NLA	NLA - Switch, Power, IRC1194A, IRC1194B	1	1
6	B,G	MS2200-NLA	NLA - Nebulizer, Disposable Sidestream	50	1
	A,B,I	MS2400-NLA	NLA - Nebulizer, Durable Sidestream	10	1
7		2002379-NLA	NLA - Fitting, Air Outlet	1	1
	A	2002593-NLA	NLA - Decal, Instruction - IRC1192/ IRC1192D	1	1
	A	2003031-NLA	NLA - Decal, Instruction - IRC1194A/ IRC1194B	1	1
	A	2002382-NLA	NLA - Owner's Manual - English/Spanish (IRC1192, IRC1192D)	1	1
	A	2003027-NLA	NLA - Owner's Manual - English/Spanish (IRC1194A, IRC1194B)	1	1
	A,H	1079289-NLA	NLA - Kit, Compressor Replacement	1	1
	A	2002930-NLA	NLA - Bumper Feet, Rubber	4	1
	A	2002387-NLA	NLA - Carton	1	1
	A	2002863-NLA	NLA - Insert, Carton	1	1

NOTE: A - Not Shown.
 B - MS2200, MS2400 can be used with either IRC1192, IRC1192A, IRC1192B, or IRC1192D. Standard for IRC1192, IRC1194A and IRC1194B is MS2200, Standard for IRC1192D is MS2400.
 G - Includes Compressor and Hardware.
 H - Includes Mouth Piece and 7 ft. Tubing.
 I - Individual item no longer available, after 10/12/04 replace entire unit.
 J - Item no longer available, refer to part #1156679 for replacement.