


# Seat Frame and Hardware

## 2nd Generation Storm, Ranger II MWD and Pronto R2

NOTE: Ranger II MWD and Pronto R2 share all the components on this page with the 2nd Generation Storm except items 7-9.



# Seat Frame and Hardware

## 2nd Generation Storm, Ranger II MWD and Pronto R2

**NOTE: Ranger II MWD and Pronto R2 share all the components on this page with the 2nd Generation Storm except items 7-9.**

Item Number	Note	Part Number	Description	Quantity Per	
				Package	Assembly
1	A	See Chart	#1 Seat Side Frame and Hardware	1	1
2		1068142-NLA	NLA - Frame, Seat, Black (14" - 20" W)	1	1
2		1068163-NLA	NLA - Seat Frame, 21" - 22" Wide	1	1
3		1025736	Locknut (1/4-20)	1	2
4		0050406	Washer, Flat (1/4 x 5/8 x 1/16")	1	4
5		0050506-NLA	NLA - Lockwasher, High Collar, Split (1/4")	1	4
6		0051006-NLA	NLA - Screw, Hex Head Full Thread (1/4-20 x 1-1/2")	1	2
7		0225100-NLA	NLA - Bumper, Tilt	1	2
8		1077358	Lockwasher, Split (3/8-16)	1	2
9		1077359	Nut, Hex (3/8-16)	1	2

NOTE: A - Includes seat side frame and items 3-6

### #1 Seat Side Frame and Hardware

Chair Depth	14"-15" Width Right	14"-15" Width Left	16"-22" Width Right	16"-22" Width Left
14"-15" Depth	1072933-NLA	1072934-NLA	N/A	N/A
16"-20" Depth	1072931-NLA See Note(s) A	1072932-NLA See Note(s) A	1072744-NLA	1072745-NLA
21"-22" Depth	N/A	N/A	1072807-NLA	1072808-NLA

NOTE: A - When part number 1072931 or 1072932 are required, call Invacare Customer Service and order as a "Special Order".