




Armrest Assembly

Non-Flip Back Arms Before 10/15/98



Armrest Assembly

Non-Flip Back Arms Before 10/15/98

Item Number	Note	Part Number	Description	Quantity Per	
				Package	Assembly
1		1041522	Arm Tube Assembly w/ Plug Buttons, Full Length, Black - Right/Left	1	1
1		1070859	Arm Tube Assembly w/ Plug Buttons, Desk Length, Black - Right/Left	1	1
2	A	1117328	Kit, Arm Pad w/ Screws, Full Length, Urethane Waterfall, Black (13" L x 2-5/8" W)	1	1
2		1110818	Kit, Arm Pad w/ Screws, Full Length, Urethane Waterfall, Black (13" L x 2-5/8" W) (Pair)	2	1
2	A	1117327	Kit, Arm Pad w/ Screws, Desk Length, Urethane Waterfall, Black (9" L x 2-5/8" W)	1	1
2		1110817	Kit, Arm Pad w/ Screws, Desk Length, Urethane Waterfall, Black (9" L x 2-5/8" W) (Pair)	2	1
3		1027095	Plug Button, Black (7/8")	1	2
4		1055949	Package, Screw, Phillips Pan Head (#10-32 x 1-3/8")	2	1
4	B	1094270-NLA	NLA - Package, Screw, Phillips Pan Head (#10-32 x 1-1/4")	2	1
5		1090768	Kit, Lock Handle w/Rivets (Pair)	2	1
6		47471X074	Spring, Arm Lock	1	1
7		1028594	Pin, Arm Lock, Black	1	1
8		1032245	Grip, Foam Hand (6-7/8" L)	1	1
9		1018225	Spring, Torsion with Button Lock	1	1
10		1028645	Rivet	1	6
11		1064130-NLA	NLA - Guard, Clothing	1	1
12		1035065	Cover, Clothing Guard w/ Pocket, FXDESK & ADJFULL/ADJDESK	1	1
13		0050406	Washer, Flat (1/4 x 5/8 x 1/16")	1	8
14		7930333	Locknut (1/4-20)	25	4
15		1024793-NLA	NLA - Screw, Hex Head (1/4-20 x 1-3/4")	1	4
16		0050506-NLA	NLA - Lockwasher, High Collar, Split (1/4")	1	4
17-21		See Chart	Item #'s 17-21 Armrest Assembly - Non-Flip Back Arms Before 10/15/98	1	1

NOTE: A - As of 1/3/02 upholstered waterfall armpads were replaced with black urethane waterfall armpads. Order in pairs with screws for chairs serviced before 1/3/02.
 B - For Black Urethane Waterfall Arm Pads Only

Item #'s 17-21 Arms

Side	17. Arm Tube w/ Arm Pad - Full Length	17. Arm Tube w/ Arm Pad - Desk Length	18. Arm Base Assembly	19. Armrest Socket - Front	20. Armrest Socket - Rear	21. Armrest Mtg. Hardware
Right 14" Deep	1110821-NLA See Note(s) A	1110819-NLA See Note(s) A	1072943-NLA	1072951-NLA	0121602-NLA	1072949-NLA
Left 14" Deep	1110822-NLA See Note(s) A	1110820-NLA See Note(s) A	1072944-NLA	1072952-NLA	0121502-NLA	1072950-NLA

Armrest Assembly

Non-Flip Back Arms Before 10/15/98

Item #'s 17-21 Arms

Side	17. Arm Tube w/ Arm Pad - Full Length	17. Arm Tube w/ Arm Pad - Desk Length	18. Arm Base Assembly	19. Armrest Socket - Front	20. Armrest Socket - Rear	21. Armrest Mtg. Hardware
Right 15" and Deeper	1110821-NLA See Note(s) A	1110819-NLA See Note(s) A	1057411-NLA	0121502-NLA	0121302-NLA	1072789-NLA
Left 15" and Deeper	1110822-NLA See Note(s) A	1110820-NLA See Note(s) A	1057412	0121602-NLA	0121402-NLA	1072930-NLA
NOTE: A - As of 1/3/02 upholstered waterfall armpads were replaced with black urethane waterfall armpads. Order right and left arm tube with arm pad for chairs serviced before 1/3/02.						