



Pivot/Slide Tube and Composite (1651) Footplates

For complete pivot/slide tube and footplates see Front Riggings Reference Chart page.



Pivot/Slide Tube and Composite (1651) Footplates

Item Number	Note	Part Number	Description	Quantity Per	
				Package	Assembly
1	A	1057852	Tube, Pivot/Slide, Curved End, Steel, Black (9-7/8")	1	1
1	A	1057854	Tube, Pivot/Slide, Curved End, 2" Longer, Steel, Black	1	1
1	A	1058418	Tube, Pivot/Slide, Curved End, 4" Longer, Steel, Black	1	1
2		1025457	Spacer, Coved (1/4 x 3/4 x 11/32")	1	2
3		1027097	Plug Button, Black (3/4")	1	1
4		7930323-NLA	NLA - Package, Locknut (1/4-20)	25	1
5		7930366	Package, Socket Head Screw (1/4-20 x 2-1/4")	25	1
6	C	1026879	Footplate, Composite, Small, Black (6" W x 6" L)	1	1
6	D	1026880	Footplate, Composite, Medium, Black (7" W x 6" L)	1	1
6	E	1026881	Footplate, Composite, Large, Black (8" W x 6" L)	1	1
6	B,F	1026879	Footplate, Composite, Small, Black (6" W x 6" L)	1	1
6	B,G	1026880	Footplate, Composite, Medium, Black (7" W x 6" L)	1	1
6	B,H	1026881	Footplate, Composite, Large, Black (8" W x 6" L)	1	1
NOTE: A - Curved B - 70STAPER and PW93ST Only C - 14" Wide D - 15-16" Wide E - 17-22" Wide F - 15-17" Wide G - 18-19" Wide H - 20-22" Wide					