

Packaging & Touch-Up Paint

Part Number	Kit Contains:
1148170-NLA	NLA - Packaging Kit Includes: (1 ea.) Carton, (1 ea.) Insert, (1 ea.) Foam Pad, and (2 ea.) Foam Block.
1149226-NLA See Note(s) A	Paint, Touch-Up, Blue
1169179-NLA See Note(s) B	Paint, Touch-Up, Onyx Blue
1149227-NLA	Paint, Touch-Up, Red
NOTE: A - To be used on scooters built before serial number 10BCT600001(L-3) and 10BCT600249 (L-4). B - To be used on scooters built after serial number 10BCT600000 (L-3) and 10BCT600248 (L-4).	