



Front Rigging Hangers and Hardware Hemi Style



Front Rigging Hangers and Hardware Hemi Style

Item Number	Note	Part Number	Description	Quantity Per		
				Package	Assembly	
1	A	1125726	Kit, Front Rigging Hanger Mounting Hardware (Pair)	1	1	
2			Screw, Socket Head Cap w/ Patch (1/4-20 x 1-1/2")	1	4	
3			NLA - Spacer (1/4 x 3/8 x 3/16")	1	4	
4			Cap, End, Side Rail, Black	1	2	
5			Block, Threaded	1	2	
6	B	1124723	Kit, Side Rail Cap Hardware (Pair)	1	1	
7			Screw, Socket Head Cap w/ Patch (1/4-20 x 1-1/2")	1	4	
8			Spacer (1/4 x 3/8 x 29/32")	1	4	
9			Cap, End, Side Rail, Black	1	2	
10			Plug Button, Black (1")	1	2	
11			Block, Threaded	1	2	
12			1143783	Front Rigging Hanger Assembly w/ Plug Button, Hemi, Black - Right	1	1
12			1143784	Front Rigging Hanger Assembly w/ Plug Button, Hemi, Black - Left	1	1
13		1134426	Package, Plug Button, Black (3/4")	10	1	

NOTE: A - Includes items 2-5
B - Includes items 7-11