


# MKIV Sip-N-Puff/ Digital Interface



# MKIV Sip-N-Puff/ Digital Interface

SNP4M On Chair  
1038635 Off Chair

Item Number	Note	Part Number	Description	Quantity Per	
				Package	Assembly
6		1038635-NLA	NLA - Sip-N-Puff w/ Digital Interface	1	1
7	A,B	1049593	Kit, SNP/DI, MK5 SAC, 4WSB, S4WSB, GTAM, AUX12M6, AUX34M6 Mtg. Hdw.	1	1
1			NLA - Half Clamp, Threaded Holes	1	1
2			Half Clamp, Non-Threaded Holes	1	1
3			NLA - Screw, Hex Head (1/4-20 x 1-1/8")	1	2
4			Screw, Phillips Pan Head (#4-40 x 3/8")	1	2
5			Plate, Mounting	1	1
NOTE: A - Includes items 1-5 B - As of 5/22/07 the Mounting Plate holes and the orientation of the Hardware was revised to attach the electronics closer to the chair. The Mounting Hardware Kit is backward compatible.					