

Invacare Bariatric Transport Chair TRHD22FR (Before 9/22/09)



Invacare Bariatric Transport Chair TRHD22FR (Before 9/22/09)

Item Number	Note	Part Number	Description	Quantity Per	
				Package	Assembly
1		1132515	Kit, Hand Grip (Pair)	1	1
2		1132507-NLA	NLA - Kit, Arm Pad w/ Screws (Pair)	1	1
3	B	1132509-NLA	NLA - Kit, Arm Assembly, Chrome - Right (Includes Arm, Clothing Guard, Arm Pad and Hardware)	1	1
3	B	1132512-NLA	NLA - Kit, Arm Assembly, Chrome - Left (Includes Arm, Clothing Guard, Arm Pad and Hardware)	1	1
4	C	1132510-NLA	NLA - Front Rigging with Hardware, Chrome - Right	1	1
4	C	1132511-NLA	NLA - Front Rigging with Hardware, Chrome - Left	1	1
5		1132506	Kit, Wheel, Front - Includes Fork and Front Wheel Assembly, (2) Bearings, Locknut and Head Tube Cap	1	2
6		1132505	Wheel Assembly with Hub Cap - Rear	1	2
7		1132513	Wheel Lock Assembly - Right - Includes Complete Wheel Lock Assembly and Mounting Hardware	1	1
7		1132514	Wheel Lock Assembly - Left - Includes Complete Wheel Lock Assembly and Mounting Hardware	1	1
	A	1141761	Package, Cap, End, Black, Head Tube	10	1
	A	1125095-NLA	NLA - Manual, Invacare Transport, Lightweight Transport, Bariatric Transport	1	1
8		1149223-NLA	NLA - Hand Brake w/Cable Assembly (Pair)	1	1
9		1149923	NLA - Upholstery Kit, Back w/Hardware	1	1
10		1149924-NLA	NLA - Upholstery Kit, Seat w/Hardware	1	1
11		1149928-NLA	NLA - Cap, Step Tube (Pair)	1	2
	D,E	1149926-NLA	NLA - Package, Spacers	25	1
	D,E	1149927-NLA	NLA - Package, Flanged Bearings	25	1

NOTE: A - Not Shown.
 B - Arm Assembly Kits in Silver Vein are currently not available. Order the right and left side Chrome Arm Kits to match.
 C - Front Rigging in Silver Vein is currently not available. Order the right and left side Chrome Front Riggings to match.
 D - Not shown.
 E - Used for Rear Wheel. Bearings and Spacers are No Longer Available. A new Rear Wheel Assembly item #6 must be ordered. Includes the Hub Cap, Flanged Bearings and Spacer also.