







Battery Chargers, MK6i Battery Cables & Sensor Cables



Battery Chargers, MK6i Battery Cables & Sensor Cables

Item Number	Note	Part Number	Description	Quantity Per	
				Package	Assembly
1	H	1053161-NLA	NLA - Charger, Battery, 8 Amp, 115vac, Off Board Dual Mode Round Plug, Lester	1	1
1		1123249	Charger, Battery, 8 Amp. 115vac (Off Board) Single Mode Round Plug w/ AC Power Line Cord	1	1
1		1123249E-NLA	NLA - WARR RPLMT - Charger, Battery, 8 Amp, 115VAC (Off Board) Single Mode Round Plug w/ AC Power Line Cord	1	1
2	G	1089627-NLA	NLA - Cord, AC Power, 115VAC, 2 Prong Off-Board (11' 6" L)	1	1
2		1121086	Cord, AC Power, 115VAC, 3 Prong On-Board (6' 6-1/2" L)	1	1
3	A	1136055	Cable and Bracket w/ Black Connector (GP24 Batteries & 4-Pole Motors)	1	1
3	B	1147019	Cable and Bracket w/ Red Connector (GP24 Batteries & GB TT Motors)	1	1
4	C	1140112	Cable, Battery Extension	1	1
5	D	1146573-NLA	NLA - Cable, Drive Lockout Sensor, MK6i	1	1
6	E	1144926	Cable, Drive Lockout Sensor, MK6i	1	1
7	F	1142231	Cable, GB TT Interface	1	1

NOTE: A - Used for GP24 Batteries only with 4-Pole Motors
 B - Used for GP24 Batteries only with GB TT Motors
 C - Used for 22NF Batteries with 4-Pole Motors only
 D - Used for chairs manufactured with ADJRECL Manual Recline or Full Recline Van Seats with 4-Pole Motors only.
 E - Sensor Cable used with item 7 Interface Cable for chairs manufactured with ADJRECL Manual Recline Seat Frame or Full Recline Van Seats with MK6i Electronics and GB TT Motors
 F - Interface Cable used with item 6 Sensor Cable for chairs manufactured with ADJRECL Manual Recline Seat Frame or Full Recline Van Seats with MK6i Electronics and GB TT Motors
 G - No longer available as a service part. Please use 1121086 as its replacement.
 H - Use Battery Charger part number 1123249 as a replacement.