


# Elevate Integrated Actuator Module (EIAM)

## 3G Storm w/ GB TT Motors & Formula CG Elevate



# Elevate Integrated Actuator Module (EIAM)

## 3G Storm w/ GB TT Motors & Formula CG Elevate

| Item Number | Note | Part Number | Description   | Quantity Per |          |
|-------------|------|-------------|---|--------------|----------|
|             |      |             |   | Package      | Assembly |
| 1           |      | 1140036     | Module, Elevate Integrated Actuator (EIAM) - MK6i (54" L) | 1            | 1        |
| 2           |      | 1052041     | Screw, Phillips Pan Head (#4-40 x 3/8")                   | 1            | 4        |