




Battery Chargers and Cables



Battery Chargers and Cables

Item Number	Note	Part Number	Description	Quantity Per	
				Package	Assembly
1	D	1053161-NLA	NLA - Charger, Battery, 8 Amp, 115vac, Off Board Dual Mode Round Plug, Lester	1	1
1	A	1123249	Charger, Battery, 8 Amp. 115vac (Off Board) Single Mode Round Plug w/ AC Power Line Cord	1	1
1	A	1123249E-NLA	NLA - WARR RPLMT - Charger, Battery, 8 Amp, 115VAC (Off Board) Single Mode Round Plug w/ AC Power Line Cord	1	1
2		1089627-NLA	NLA - Cord, AC Power, 115VAC, 2 Prong Off-Board (11' 6" L)	1	1
2	C	1121086	Cord, AC Power, 115VAC, 3 Prong On-Board (6' 6-1/2" L)	1	1
3		60102824	Kit, Cable, Battery Extension	1	1
4	B	1146573-NLA	NLA - Cable, Drive Lockout Sensor, MK6i	1	1
NOTE: A - CTE charger is standard on-chair. Used on Gel Cell or AGM (Absorbed Glass Mat) batteries. For wet cell batteries refer to the Lester charger part number 1053161. B - Drive Lockout Sensor Cable used for TDX SC model chairs manufactured with a ADJRECL Manual Recline Seat Frame C - No longer available as a service part. Please use 1121086 as its replacement. D - Use Battery Charger part number 1123249 as a replacement.					