




# Shrouds and Hardware (Elevate Only CGE) TDX SC and (TDX SPREE with or without Tilt)



**Top Shroud with TRBKTS Option**

WC19/Engage  
Docking Pins Warning  
Label for TDX SPREE  
is non-serviceable.



TRRO Adjustment Label  
is non-serviceable

**Top Shroud with TRRO Option**

## Shrouds and Hardware (Elevate Only CGE) TDX SC and (TDX SPREE with or without Tilt)

Item Number	Note	Part Number	Description	Quantity Per	
				Package	Assembly
1		1152728	Package, Screw, Thumb (10-32 x 7/8")	10	1
2	A	1152726	Shroud with Decals, Top Rear Elevate, Black	1	1
3		1119281-NLA	NLA - Plate, Name, Patent Numbers	1	1
4		1035900	Decal, Transportation	1	1
5	B	1154219-NLA	NLA - Shroud with Decals, Top Forward Elevate, Black	1	1
6		1079203	Decal, Warning Pinch Point	1	2
7	C	1149271	Shroud with Decals, Black - Rear	1	1
9	D,G	8881154217	Shroud with Decals, Top, Color (TRBKTS Option)	1	1
9	E,G	8881154218-NL A	NLA - Shroud with Decals, Top, Color (TRRO Option For TDX SPREE)	1	1
10		1083199-NLA	NLA - Label, Hook Tie-Down	1	4
11	F,G	8881151619	Shroud with Decal, Color - Front	1	1
<p>NOTE: A - Includes Shroud with items 3 and 4            B - Includes Shroud with item 6            C - Includes Shroud with item 8            D - Includes Shroud with item 10            E - Includes Shroud, item 10, WC 19/Engage Docking Pins Warning Label and TRRO Adjustment Label            F - Includes Shroud with item 8            G - Specify frame color when ordering. See page titled "Frame Finishes, Touch-Up Paints and Upholstery Options".</p> <p>Item #8 Reflective Decal is not available for sale individually.</p>					