




Arm Mounting Swingaway Hardware ADJASBA/ADJRECL Seats



Swingaway Armrest Mounting Hardware
ADJASBA Non-Manual Tilt Seat Frame &
ADJASBA Manual Tilt Seat Frame
(12"-15" Deep)
TDX SPREE



Swingaway Armrest Mounting Hardware
ADJASBA Non-Manual Tilt Seat Frame &
ADJASBA Manual Tilt Seat Frame
(16"-18" Deep)
TDX SPREE



Swingaway Armrest Mounting Hardware
ADJASBA Seat Frame (16"-20" Deep)
TDX SC

Arm Mounting Swingaway Hardware ADJASBA/ADJRECL Seats

Item Number	Note	Part Number	Description	Quantity Per	
				Package	Assembly
1	A	1068488	Kit, Swingaway Armrest Mounting Hardware	1	1
2		1019590	Screw, Socket Button Head (1/4-20 x 1-5/8")	1	1
3		1068025	Clevis, Pivot	1	1
4		1025195	Locknut (1/4-20)	1	1
5		1068039	Pin, Pull	1	1

NOTE: A - Includes items 2-5