



# Pivot/Slide Tube and Composite Footplates

NOTE: Heel loops are not included. See footrest and legrest accessories.

NOTE: Heel loops are not included. See footrest and legrest accessories.



# Pivot/Slide Tube and Composite Footplates

NOTE: Heel loops are not included. See footrest and legrest accessories.

Item Number	Note	Part Number	Description	Quantity Per	
				Package	Assembly
1	C	1120383	Tube, Pivot/Slide, Curved End, Black - Short (8-1/4")	1	1
1	D	1120382	Tube, Pivot/Slide, Curved End, Black - Long (9-7/8")	1	1
1	D	1120386	Tube, Pivot/Slide, Curved, 2" Longer, Black (11-7/8")	1	1
1	E	1120384-NLA	NLA - Tube, Pivot/Slide, Straight - Alum. (9-7/8" L)	1	1
1	E	1120385	Tube, Pivot/Slide, Straight - Black (11-7/8" L)	1	1
2		1025457	Spacer, Coved (1/4 x 3/4 x 11/32")	1	2
3		1027097	Plug Button, Black (3/4")	1	1
4		7930323-NLA	NLA - Package, Locknut (1/4-20)	25	1
5		7930366	Package, Socket Head Screw (1/4-20 x 2-1/4")	25	1
6	A	1026879	Footplate, Composite, Small, Black (6" W x 6" L)	1	1
6	B	1026880	Footplate, Composite, Medium, Black (7" W x 6" L)	1	1
7		See Chart	#7 Pivot/Slide Tube w/Composite Footplate Assemblies	1	1

NOTE: A - 13-14" Wide  
 B - 15-16" Wide  
 C - 60 Degree MFX  
 D - 60 Degree  
 E - 70 Degree NTaper

## #7 Pivot/Slide Tube w/Composite Footplate Assemblies

(without Heel Loops)

Front Rigging	Small 13"-14" Wide	Medium 15"-16" Wide
60 Degree (13"-19") -Right	1028391	1028387
60 Degree (13"-19") -Left	1028392	1028388
60 Degree (13"-19") 2" Longer -Right	1032627-NLA	1032625-NLA
60 Degree (13"-19") 2" Longer -Left	1032628-NLA	1032626-NLA
60 Degree MFX (9-3/4" - 14") -Right	1028389-NLA	1028385-NLA
60 Degree MFX (9-3/4" - 14") -Left	1028390-NLA	1028386-NLA
70 Degree (14"-19") -Right	1042784-NLA	1042786
70 Degree (14"-19") -Left	1042785-NLA	1042787
70 Degree (14"-19") 2" Longer -Right	1042790	1042792
70 Degree (14"-19") 2" Longer -Left	1042791	1042793