


Fortuna Model II

(Since 00-10-02)

Since serial number 662527



Fortuna Model II (Since 00-10-02)

Item Number	Note	Part Number	Description	Quantity Per	
				Package	Assembly
1		14535	Cover, Seat, White (5 mm)	--	--
1		14556	Cover, Seat, Blue (5 mm)	--	--
2		14593-NLA	NLA - Cover, Backrest, White (5 mm)	--	--
2		14595-NLA	NLA - Cover, Backrest, Blue (5 mm)	--	--
1+2		14980-NLA	NLA - Cover, Seat and Backrest, White (5 mm)	--	--
1+2		14981-NLA	NLA - Cover, Seat and Backrest, Blue (5 mm)	--	--
3		14592-NLA	NLA - Seat, Complete without Cover	--	--
4		14453-NLA	NLA - Flap, Side (525 x 160 mm)	--	--
5		14041	Rivet (3.2 x 9.5L)	--	--
6		11548	Hinge	--	--
7		10933	Cover, C-Rail	--	--
9		15072-NLA	Cups, Suction with Toggle (Pair)	--	--
10		15071	Cups, Suction, Screw-in (Pair)	--	--
11		15121	Cap for Suction Cup	--	--
12		14462-NLA	NLA - U-Clip	--	--
13a		14377-NLA	NLA - Backrest, Complete without Cover	--	--
13b		14452-NLA	NLA - Backrest (Frame Only)	--	--
14	B	15815-NLA	NLA - Actuator	--	--
15		15829	Clamp for Actuator	--	--
16		14468-NLA	NLA - Bolt, Motor	--	--
17		15183	Circuit Board	--	--
18		16093-NLA	NLA - Casing, Backrest	--	--
19		14457	Gasket, Casing	--	--
20		14451	Cup, Compression	--	--
21		14441	Flange for Compression Cup	--	--
22		14885-NLA	NLA - Clamp	--	--
23		15394-NLA	NLA - Hand Control with Spiral Cable and Screws	--	--
24		14470	Spindle Catch	--	--
25		14502	Ring, Rubber (60mm x 1.5mm)	--	--
26		14978	Locking Device	--	--
27		14929	Spring	--	--
28	C,D	15170	Kit, Conversion Height Adapter (2 cm)	--	--
29	C,D	14694	Adapters, Height (Set)	--	--
30		15239	Charger, Battery, Canada / USA	--	--
31		15635-NLA	NLA - Reduction	--	--
32		15224	Fuse, Glass (5 x 20) 5 AMP	--	--
33	E	12V-3	Battery, Complete, Sealed Lead Acid	--	--
	A	15298	Kit, Conversion without Height Adapters	--	--
	A	15877-NLA	NLA - Kit, Recirculation without Battery Charger, White - Fortuna	--	--
	A	15878-NLA	NLA - Kit, Recirculation without Battery Charger, Blue - Fortuna	--	--
	A	15243	Cover, Backrest, Hydraulically Recline, White (5 mm)	--	--
	A	15244	Cover, Backrest, Hydraulically Recline, Blue (5 mm)	--	--
	A	15466	Cover, Head Support, White (5 mm)	--	--

Fortuna Model II (Since 00-10-02)

Item Number	Note	Part Number	Description	Quantity Per	
				Package	Assembly
	A	15467	Cover, Head Support, Blue (5 mm)	--	--
NOTE: A - Not Shown. B - Includes Motor, Pulley Assembly and Spindle Foot. C - For height adaptation up to 4cm, you must order (1) Height Adapter Set Part # 14694 and (1) Conversion Kit Part # 15170. D - For height adaptation up to 6cm, you must order (2) Height Adapter Set Part # 14694 and (1) Conversion Kit Part # 15170. E - To use this as the replacement, the existing battery casing must be used. Contact customer service for more details.					