



# Wheel Locks and Hardware



## Wheel Locks and Hardware

Item Number	Note	Part Number	Description	Quantity Per	
				Package	Assembly
1	A,C	1061612-NLA	NLA - Wheel Lock Assembly - Right	1	1
1	A,C	1061613-NLA	NLA - Wheel Lock Assembly - Left	1	1
1	D	1141278-NLA	NLA - Wheel Lock Assembly - Right	1	1
1	D	1141279	Wheel Lock Assembly - Left	1	1
2		1136024-NLA	NLA - Package, Rubber Tip, Wheel Lock Handle, Black	10	1
3	B	1034136	Kit, Wheel Lock Attaching Hardware	2	1
4			NLA - Bolt, Hex Head Cap (1/4-20 x 3/4")	1	2
5			Washer (1/4 x 5/8 x 3/32")	1	2

NOTE: A - Wheel lock handle design changed after 2/17/99. Chairs manufactured before 2/18/99, order in pairs to match. These wheel locks were used 2/18/99-6/28/06  
 B - Item 3 is a pair number that includes 2 ea. of items 4 and 5  
 C - Before 6/29/06. No Longer Available, order a pair of Wheel Locks dated after 6/28/06  
 D - After 6/28/06