


# MKIV R.I.M. Head Control Mounting Hardware (1500M4) (After 8/15/95)



## MKIV R.I.M. Head Control Mounting Hardware (1500M4) (After 8/15/95)

Item Number	Note	Part Number	Description	Quantity Per	
				Package	Assembly
1	A	1055149	Cover, Nylon	1	1
2		1055148-NLA	NLA - Clamp with Knob, Two-Piece	1	1
3		1003645	Screw, Socket Button Head (#10-32 x 3/4")	1	4
4		7930326-NLA	NLA - Package, Screw, Hex Head (1/4-20 x 3/4")	25	8
5	A	1054931-NLA	NLA - Bracket, Attaching, Black (12 - 18" W)	1	1
6		1031095	Screws with Washers (#10-16 x 7/8")	4	1
6		1076284	Screw w/ Washer (#10 x 7/8")	1	4
7		1054929-NLA	NLA - Bracket, Wing, Black	1	2
8		1026907	Washer (1/4 x 1/2 x 1/16")	1	8
9		7930333	Locknut (1/4-20)	1	8
10		1026034-NLA	NLA - Washer, Nylon (3/16 x 7/16 x 1/8")	1	4
11		7930332-NLA	NLA - Locknut, High Profile (#10-32)	25	4
12	A	1056347	Hardware, Rim Control Mounting (RIMHW) (12-18" W)	1	1

NOTE: A - 12-16" Wide