


Fixed Back Assembly



Fixed Back Assembly

Item Number	Note	Part Number	Description	Quantity Per	
				Package	Assembly
1-5		440587	Kit, Fixed Angle Back	--	1
1-5		440107	Bracket, Back Post Mountain - Right	1	--
		440108	Bracket, Back Post Mountain - Left	1	--
2		318487	Screw, Phillips Pan Head (1/4-20 x 7/8")	4	--
3		1030949	Screw, Hex Head (1/4-20 x 1-1/2")	4	--
4		1025736	Locknut (1/4-20)	4	--
5		1025704	Washer, Flat (5/16 x 5/8 x 1/16")	4	--