




Pivot/Slide Tube w/ Footplate Assembly - Composite - 1651



Small Footplate



Medium Footplate



Large Footplate

Pivot/Slide Tube w/ Footplate Assembly - Composite - 1651

60 Degree	Small	Small	Medium	Medium	Large	Large	
13" - 19"	Long	2" Longer	Long	2" Longer	Long	2" Longer	
J Curved End	14"	14"	15" - 16"	15" - 16"	17" - 22"	17-22"	
	N/A	N/A	1028387	N/A	1028383	N/A	Right
	N/A	N/A	1028388	N/A	1028384	N/A	Left
60 Degree MFX	Small	Medium	Large				
9-3/4" - 14"	Short	Short	Short				
J Curved End	14"	15" - 16"	17" - 22"				
	N/A	1028385-NLA	1028381-NLA				
	N/A	1028386-NLA	1028382-NLA				
70 Degree	Small	Small	Medium	Medium	Large	Large	
14" - 19"	Short	2" Longer	Short	2" Longer	Short	2" Longer	
J Straight	14"	14"	15" - 16"	15" - 16"	17" - 22"	N/A	
	1042784-NLA	N/A	1042786	N/A	1042788	N/A	Right
	1042785-NLA	N/A	1042787	N/A	1042789		Left
70 Degree Tapered	Small	Small	Medium	Medium	Large	Large	
14-1/2" - 19-1/2"	Short	2" Longer	Short	2" Longer	Short	2" Longer	
J Straight	16" - 17"	16" - 17"	18" - 19"	18" - 19"	20" - 22"	20" - 22"	
	N/A	N/A	1042786	N/A	1042788	N/A	Right
	N/A	N/A	1042787	N/A	1042789	N/A	Left

Pivot/Slide Tube w/ Footplate Assembly - Composite - 1651

Item Number	Note	Part Number	Description	Quantity Per	
				Package	Assembly
1		1026879	Footplate, Composite, Small, Black (6" W x 6" L)	1	1
1		1026880	Footplate, Composite, Medium, Black (7" W x 6" L)	1	1
1		1026881	Footplate, Composite, Large, Black (8" W x 6" L)	1	1
2		60285X006	Screw, Socket Head (1/4-20 x 2-1/4")	1	1
		1025195	Locknut (1/4-20)	1	1
		1027097	Plug Button, Black (3/4")	1	1
		1025457	Spacer, Coved (1/4 x 3/4 x 11/32")	1	2
		1026703	Kit, Pivot/Slide Tube Mounting Hardware (Pair)	1	1
		See Chart	1651 Composite Footplates and Pivot Tubes Assemblies	1	1