



Hublock Replacement Kits - Before 11/01/07

Attendant and User Hand and Foot Operated

<p>12" Rear Wheel</p> <p>Foot - 1153658 Hand - 1153657</p> <p>12COM</p> <p>Mount Axle Bracket 4 holes over for non-vent tray & non- recliner chairs</p>	<p>18" Rear Wheel</p> <p>Foot - 1152207 Hand - 1153655</p> <p>18COM</p> <p>Mount Axle Bracket 3 holes over for non-vent tray & non- recliner chairs</p> <p>Low</p>	<p>Foot - 1153658 Hand - 1153657</p> <p>High</p>	<p>Vent Tray Wheelbase Position</p> <p>Foot - 1152207 Hand - 1153655</p>
<p>20" Rear Wheels</p> <p>Foot - 1152207 Hand - 1153655</p> <p>20COM</p> <p>Mount Axle Bracket 3 holes over for non-vent tray & non- recliner chairs</p>	<p>22" Rear Wheels</p> <p>Foot - 1152207 Hand - 1153655</p> <p>22COM</p> <p>Mount Axle Bracket 3 holes over for non-vent tray & non- recliner chairs</p>	<p>24" Rear Wheels</p> <p>Foot - 1152207 Hand - 1153655</p> <p>24COM7</p> <p>Mount Axle Bracket 3 holes over for non-vent tray & non- recliner chairs</p>	<p>Recliner Wheelbase Position</p> <p>Recliner is excluded with 18COMPA2XLE combo</p> <p>Foot - 1152207 Hand - 1153655</p>



Hublock Replacement Kits - Before 11/01/07

Attendant and User Hand and Foot Operated

Item Number	Note	Part Number	Description	Quantity Per				
				Package	Assembly			
1	A,B	1152207	Kit, Hublock Replacement, Foot Operated	1	1			
			Actuator Assembly, Foot Operated	1	1			
			Latch with Cable Assembly, Hublock - Right	1	1			
			Latch with Cable Assembly, Hublock - Left	1	1			
			Locknut (1/4-28)	1	6			
			NLA - Service Manual, Hublock	1	1			
			NLA - Kit, Hublock Replacement, Foot Operated	1	1			
1	A,B	1153658-NLA	Actuator Assembly, Foot Operated	1	1			
			Latch with Cable Assembly, Hublock - Right	1	1			
			Latch with Cable Assembly, Hublock - Left	1	1			
			Locknut (1/4-28)	1	6			
			NLA - Service Manual, Hublock	1	1			
			NLA - Kit, Hublock Replacement, Hand Operated	1	1			
			Actuator w/ Cable Assembly, Hand Operated	1	1			
2	A	1153657-NLA	Hublock Hand Operated Latch Assembly	1	1			
			Hublock Hand Operated Latch Assembly	1	1			
			Screw, Socket Button Head (#10-32 x 1")	1	2			
			Box, Junction Cable	1	2			
			Pin, Dowel (3/16 x 1/2")	1	1			
			Locknut (#10-32)	1	2			
			Bracket, Junction Cable, Black	1	1			
			Tie Wrap (5-5/8" L)	1	1			
			Locknut (1/4-28)	1	4			
			Box, Junction Cable	1	1			
			Instruction Sheet	1	1			
			2	A	1153655	Kit, Hublock Replacement, Hand Operated	1	1
						Actuator w/ Cable Assembly, Hand Operated	1	1
						Hublock Hand Operated Latch Assembly	1	1
						Hublock Hand Operated Latch Assembly	1	1
Screw, Socket Button Head (#10-32 x 1")	1	2						
Box, Junction Cable	1	2						
Pin, Dowel (3/16 x 1/2")	1	1						
Locknut (#10-32)	1	2						
Bracket, Junction Cable, Black	1	1						
Tie Wrap (5-5/8" L)	1	1						
Locknut (1/4-28)	1	4						
Box, Junction Cable	1	1						
Instruction Sheet	1	1						

NOTE: A - Refer to image in order to determine appropriate selection.

B - Part number 1130215 is no longer available from the Invacare Service Parts Department. Refer to the Invacare Website under "Product Literature Search". Enter the part number in the "Keywords" field and download a free PDF version copy.

Invacare does not recommend the use of service part number 1132357 and/or 1132358 for replacement hub lock cables for wheel chairs manufactured before November 1, 2007.