




Shroud Covers, Recline Actuator

Conv. & Contoura - TDX SP/SP-GT/SR (3G Storm After 12/9/07)



Formula CG Recline Actuator Shroud Cover
16"-20" Wide Conventional



Formula CG Recline Actuator Shroud Cover
16"-24" Wide Contoura

Formula CG Recline Actuator Shroud Cover
20"-24" Wide Conventional

Shroud Covers, Recline Actuator

Conv. & Contoura - TDX SP/SP-GT/SR (3G Storm After 12/9/07)

Item Number	Note	Part Number	Description	Quantity Per	
				Package	Assembly
1		1148698	Cover with Decal, Conventional (16 - 20" W)	1	1
1		1147139	Cover with Decal, Conventional (20 - 24" W)	1	1
1		1147959	Cover with Decal, Contoura (16 - 24" W)	1	1
3		1117240	Screw, Thumb (#10-32 x 1/2")	1	1

Item #2 Reflective Decal is not available for sale individually.