

Controller - Programmer - Shroud Cover - MK6i TRAM w/ Mtg. Hdwr. - 3G Storm w/ 4 Pole Motors and Formula CG

NOTE: After 12/9/07 the TRAM Module was relocated next to the Recline Actuator Conventional Motor for Tilt/Recline & Recline Only systems

Item Number	Note	Part Number	Description	Quantity Per	
				Package	Assembly
1	A	1144784-NLA	Kit, Controller Shroud with Hardware	1	1
2			Shroud, Controller	1	1
4			NLA - Retainer, Stud (1/4" Turn)	1	2
5			NLA - Washer, Nylon Retainer (1/4 x 1/2 x 1/32")	1	2
6			NLA - Washer, Nylon Cupped (9/32 x 19/32 x 3/32")	1	2
7			NLA - Spring, Torsion Ejector	1	2
8			NLA - Stud, Phillips Recess (1/4" Dia. x 1" Long, 1/4 Turn)	1	2
3	B	1104706-NLA	NLA - Kit, Controller Shroud Mounting Hardware (Pair)	1	1
4		1110509-NLA	NLA - Retainer, Stud (1/4" Turn)	1	2
5		1070580-NLA	NLA - Washer, Nylon Retainer (1/4 x 1/2 x 1/32")	1	2
6		1070581-NLA	NLA - Washer, Nylon Cupped (9/32 x 19/32 x 3/32")	1	2
7		1093312-NLA	NLA - Spring, Torsion Ejector	1	2
8		1110508-NLA	NLA - Stud, Phillips Recess (1/4" Dia. x 1" Long, 1/4 Turn)	1	2
9	C	1145381-NLA	NLA - Kit, Controller Bracket with Hardware	1	1
10			Bracket, Controller Mounting, Black	1	1
11			Locknut (#8-32)	1	2
12			Screw, Hex Head w/ Patch (5/16-18 x 1")	1	2
13			Screw, Phillips Pan Head (#4-40 x 3/8")	1	4
14			Locknut (#8-32)	1	2
15		1136898	Controller, MK690 w/ ACC	1	1
16		1140098	Module, Tilt Recline Actuator (TRAM) - MK6i (36" L)	1	1
16		1140098E-NLA	NLA - WARR RPLMT - Module, Tilt Recline Actuator (TRAM) - MK6i	1	1
17	G	1139985	Programmer, MK6i Remote	1	1
17	G	1142215	Kit, MK6i Remote Programmer w/ SD Professional Card & Card Holder (1813M6)	1	1
18	E	1140125-NLA	NLA - Jumper, DCI	1	1
19		1139997	Cover, SD Card Slot	1	1
	D	1095502-NLA	NLA - Package, Tie Wrap (11-1/2" L)	10	1
	F	1101400	Harness, ELRPW Power Legrest - Right	1	1

NOTE: A - Includes items 2 and 3.
 B - Includes items 4-8.
 C - Includes items 10-14.
 D - Not Shown
 E - DCI Jumper is used when an MPJ+ Joystick is used
 F - Not Shown. ELRPW Power Legrest Harness part #1101400 is used to connect between the TRAM Module and the Tilt Actuator as well
 G - The MK6i Programmer Manual (1141471) is not included and is strongly recommended the guide be ordered or downloaded from Literature Fulfillment.