

# Rollator (Four Wheeled) - Model 65350GR



## Rollator (Four Wheeled) - Model 65350GR

Item Number	Note	Part Number	Description	Quantity Per	
				Package	Assembly
1	A	1151023	Grip, Hand (Pair)	2	1
2	E	1153851-NLA	NLA - Brake Assembly - Left	1	1
3	E	1153850-NLA	NLA - Brake Assembly - Right	1	1
4		1151511-NLA	NLA - Backrest Assembly	1	1
5	B,G	1151021-NLA	NLA - Knob w/ Hardware (Pair)	1	1
5	H	1169007	Kit, Adjustable Knob (Pair)	2	1
6		1151015-NLA	NLA - Seat w/ Mounting Hardware	1	1
7		1151510	Basket	1	1
8		1151024-NLA	NLA - Package, Cap	6	1
9	C	1151022-NLA	NLA - Package, Springs, Brake (Pair)	2	1
10	D	1151019-NLA	NLA - Kit, Caster with Hardware - Rear (Pair)	2	1
11	F	1151018-NLA	NLA - Kit, Caster with Fork and Hardware - Front (Pair)	2	1

NOTE: A - Includes (2) Hand Grips.  
 B - Includes (2) Knobs with Attaching Hardware.  
 C - Includes (2) Springs.  
 D - Includes (2) Casters and Attaching Hardware.  
 E - Includes (1) Brake Handle, (1) Cable, (1) Brake Pad and Hardware.  
 F - Includes (2) Casters, (2) Forks and Attaching Hardware.  
 G - Part is no longer available as of 2/6/2012. Order Part number 1169007 as a replacement.  
 H - Washer and Locknut not required. Replacement Knob screws directly into adjustment sleeve.